Füllstandskontrolle
SIGNALBOX T1
SIGNALBOX T5

Level Control
SIGNALBOX T1
SIGNALBOX T5

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Thank you

EN

for purchasing this SCAT Europe product. Please read this operating manual carefully before use and especially observe the safety notices! For service and support, please contact our customer service or your local distributor. For a distributor in your country, please contact us on www.scat-europe.com.

SCAT Europe GmbH
Waldecker Straße 7
D-64546 Mörfelden
Germany

Phone: +49-(0)6105 - 305 586 - 0
Fax: +49-(0)6105 - 305 586 - 99
e-Mail: info@scat-europe.com
Web: www.scat-europe.com

Content

Device overview Signalbox T1	20
Device overview Signalbox T5	21
Technical data	22
Safety advice.....	22
First-Time-Operation.....	23
LED index.....	23
Level alert.....	24
Controlling external devices	26
Connection with Switch Amplifier.....	27
Switch Amplifier installation	28
Switch Amplifier installation	30
Troubleshooting	32
Warranty conditions.....	33


**Scan QR Code =
PDF Download**


Device overview Signalbox T1

EN


- 1** Power Green power-LED
- 2** Snooze & Reset Snooze & Reset button
- 3** Status LED Status channel 1


- 4** OUTPUT 1 Potential-free switch, channel 1
- 6** DC INPUT 24V electrical power supply
- 7** INPUT LEVELSENSOR Sensor input, channel 1

Device overview Signalbox T5

EN

- 1** Power Green power-LED
- 2** Snooze & Reset Snooze & Reset button
- 3** Status LEDs Status channel 1-5

- 4** OUTPUT 1-5 Potential-free switch, channel 1-5
- 5** OUTPUT ALL Potential-free switch, all channels
- 6** DC INPUT 24V electrical power supply
- 7** INPUT LEVELSENSOR Sensor input, channel 1-5

Technical data

Dimensions (W x H x D)	180 x 105 x 55 mm
Weight	approx. 620 g
Power supply	24V power supply
Switching power	24V/2A
Power consumption	max. 4W
Operating temperature	+5°C bis +40°C
Display	LEDs

Safety advice

- Always put the devices in a safe and dry place, and protect them from moisture.
- Do not drop the devices. Fix cables safely and do not pull them by force.
- **Important notice:** Before operation, please ensure correct function of the level sensor (Disc sensor or SafetyWasteCap). The floater on your SafetyWasteCap might be stuck or unmovable by sticky or viscous fluids. Check floater and level control mechanism constantly for dirt or contaminants. The floater must be freely movable. Clean and purge the floater if necessary with a suitable solvent or clean water.

First-Time-Operation

EN

1. Power supply

Connect the Signalbox to your power supply, using the delivered ac/dc adapter **6**. Now the power-LED lights up in green.

2. Connect signal cable

Connect the level sensor with the desired channel of the Signalbox **7** by using the signal cable. Now the Status-LED **3** lights up in green.

3. Optional: Control of external devices

Please follow the references on page 26.

LED index

Power LED

	Off	Device off
	Green	Device on

Status LEDs

	Off	No sensor connected
	Green	Sensor connected
	fast red flashing	Max. filling level reached
	slow red flashing	Max. filling level reached, Snooze-Mode is activated

Level alert

If a **container hits the max. filling level**, the Signalbox gives an **acoustic warning signal** and the corresponding **LED 3 flashes red in a fast interval**. Please proceed as follows:

1. Set the alarm in Snooze-Mode

Press the Snooze & Reset button 2. Now the acoustic warning signal becomes silent for 60 seconds and the LED of the corresponding channel **3** flashes red in a slow interval.

The container must be emptied now.

If the Signalbox will not be reseted within the next 60 seconds, the acoustic warning signal will start again with LED 3 flashing red in fast intervals regarding the filling level state.

If another container hits the max. filling level during this time (affects only Signalbox T5) the Signalbox will also alert acoustically and fastly flashing, even if a different channel is resting in Snooze-Mode. This container has to be emptied now also.

2. Empty the container

Dispose the container contents. Use the mandatory way of disposal!

Warning: Always wear adequate and suitable protective clothing when working with hazardous chemicals!

Attention: Do not dispose chemical waste through your canalization or domestic waste! Follow your local laws and directives when disposing fluid waste or other chemicals! Please ask your official safety inspector and your environmental authorities for more information and applicable law!

3. Reset

Hold the Snooze & Reset button 2 for 3 seconds to reset all channels. Now the Status LEDs 3 change to green light.

For those containers that still hit max. fill level, the alert will recur automatically.

Press and hold of the Snooze & Reset button is approved by a key tone.

Controlling external devices

EN

The OUTPUT **4** reports the status of the corresponding sensor through a potential-free contact, which can be used to control external devices. The status can only be reset by holding the Snooze & Reset button **2** for 3 seconds. The output status will also be maintained when the sensor is disconnected.

Sensor status	Status OUTPUT 4
Below max. filling level „Error-free condition“	Electrical contact between Pin 1 and 3 
Max. filling level exceeded „Error condition“	Electrical contact between Pin 3 and 4 

Pinout: 1=Brown, 3=Blue, 4=Black

Only Signalbox T5

The output „ALL“ **5** responds to all connected sensors, regardless of their number or channel.

Connection with Switch Amplifier

For operation in explosion-endangered areas (Ex zone) the SafetyWasteCap has to be connected to a Switch Amplifier. Signalboxes and Switch Amplifiers have to be installed outside the ex zone.

EN

The Disc Sensors 108 048 and 108 045 are not compatible with a Switch Amplifier. Please note that they only can be connected directly with the built-in signal box and are only to be used outside of explosion-endangered-areas.

S.C.A.T. Europe recommends the manufacturers Rechner Sensors Switch Amplifier. Tests with Signalbox and the use of the below listed device guarantees a trouble-free operation.

Switch Amplifier

Rechner Sensors - N-132/2-01 N00015

**The use of electronic components in explosion endangered areas is only allowed if they are protected with a Switch Amplifier. The installation has to be set up by qualified personnel only!
Installation instruction on page 28-31.**

Order the Switch Amplifier:

Part Number: 108 278 Rechner Sensors - N-132/2-01 N00015

+ Cable set for connection, Part Number: 108 219

SCAT Europe GmbH


Tel.: +49-(0)6105 - 305 586 - 0

Fax: +49-(0)6105 - 305 586 - 99

e-Mail: info@scat-europe.com

Switch Amplifier installation


EN

Switch Amplifier installation

EN

In use of intrinsically safe operation the cable should be marked in blue!

EN

Troubleshooting

EN

Alert rings out, although the max. filling level has not reached.

- Make sure that cables are connected and in good condition. If cables are damaged or broken replace them.
- Make sure that the level control is mounted correctly at the container. The cable connection must be tightened accurately.
- The floater on your SafetyWasteCap might be stuck or unmovable by sticky or viscous fluids. Check the floater mechanism for pollution or contaminants. The floater must be freely movable. Clean and purge the floater with a suitable solvent or clean water.

Warranty conditions

Our rigorous quality control ensures that you will receive high quality products without defects from us. If a product should nonetheless prove to be defective, you will receive a free of charge replacement of course. S.C.A.T. Europe only provides warranty according to current warranty conditions regulated by European law. We offer no warranty for products which are technically modified or damaged by the customer. In case of any questions, please contact our head office or your official distributor. Contact details can be found at:

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
User Manual

**PDF
Download**


The latest version of **User Manual** relating to your product is always available for you to download. Please visit www.scatt-europe.com or scan this **QR-Code** using the camera of your Smartphone/Tablet.

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A cartoon illustration of a hand with red polka dots on the back, wearing a blue sleeve, pointing towards a white rectangular sign. The sign contains the text 'Ihr Fachhändler vor Ort: Your local distributor:'.

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