## Product Instructions

SSL100 / SSL10BPW / SSL10NB / SSL10NB2G /  
SSL10NBDBL / SSL10NB6X9 / SSL10LET / SSL10LET2G / SSL10DE

-  (EN) **Sponge-Stick**
-  (FR) **Eponge Sponge-Stick**
-  (DE) **Sponge-Stick**
-  (IT) **Spugna con supporto in plastica**
-  (ES) **Sponge-Stick**
-  (NL) **Sponge-Stick**
-  (SV) **Sponge-Stick**
-  (DA) **Sponge-Stick**
-  (NO) **Sponge-Stick fuktet svaber**
-  (FI) **Sponge-Stick Sienitikku**
-  (PT) **Sponge-Stick**
-  (EL) **Σπόγγος με Χειρολαβή**
-  (PL) **Sztyft z gąbką**
-  (HU) **Sponge-Stick szivacsos mintavételi pálca**
-  (CS) **Tučinka s houbičkou**
-  (RO) **Burete de recoltare**
-  (RU) **Губка с держателем для отбора проб**
-  (TR) **Sponge-Stick**
-  (JA) **スポンジスティック**
-  (ZH) **海绵涂抹棒**
-  (TH) **สปอนจ์สติ๊ก**
-  (KO) **스폰지-스틱**

# Product Instructions

## Sponge-Stick

### PRODUCT DESCRIPTION

Dry or Pre-hydrated and biocide-free 3.8 x 7.6 cm (1.5 x 3 in.) sponge in sample bag. Pre-hydrated with specified agent. Various bag sizes. With or without gloves.

For Laboratory Use Only. Material Safety Data Sheet available.

### SAFETY

The user should read, understand, and follow all safety information in the instructions. Retain the safety instructions for future reference.

**⚠ WARNING** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

**NOTICE** Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

### ⚠ WARNING

#### To reduce the risks associated with environmental contamination:

- Follow current industry standards and local regulations for disposal of contaminated waste.

#### To reduce the risks associated with exposure to biohazards:

- Dispose of samples according to all applicable government regulations and applicable laboratory procedures.

#### To reduce the risk associated with false negatives resulting in the use of contaminated carcass or environmental surfaces for food or beverage products:

- Always reference package label for storage instruction and expiration date.
- Always reference product instruction for usage.

#### To reduce the risks associated with exposure to chemicals and biohazards:

- Always follow standard good laboratory safety practices (GLP<sup>1</sup> or ISO 17025<sup>2</sup>), including proper containment procedures, wearing appropriate protective apparel while handling testing materials and test samples.

### NOTICE

#### To reduce the risk of false-positive results due to cross contaminated carcass or environmental surfaces for food or beverage products that may result in re-testing or the rejection of food or beverage product:

- Do not touch the sponge or sponge stick to any unintended surface.
- Do not break the sponge from the sponge stick while sampling.
- Do not reach into the sponge bag.

#### To reduce the risk of cross contamination from reuse of sample handling device:

- Do not use the same sponge stick more than once.

### USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.3M.com/foodsafety](http://www.3M.com/foodsafety), or contact your local 3M representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any 3M Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

### LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, 3M DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any 3M Food Safety Product is defective, 3M or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify 3M within sixty days of discovery of any suspected defects in a product and return it to 3M. Please call Customer Service (1-800-328-1671 in the U.S.) or your official 3M Food Safety representative for a Returned Goods Authorization.

## LIMITATION OF 3M LIABILITY

3M WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall 3M's liability under any legal theory exceed the purchase price of the product alleged to be defective.

## STORAGE AND DISPOSAL

**Storage:** Refer to package label for storage information.

**Disposal:** Follow current industry standards and local regulations for disposal of contaminated waste. Consult the Material Safety Data Sheet for additional information.

Product Description	Volume	REF	Volume of Bag	Gloves
3M™ Sponge-Stick	0 mL	SSL100	0.89 L (30 oz.)	0
3M™ Sponge-Stick with Buffered Peptone Water Broth	10 mL	SSL10BPW	0.89 L (30 oz.)	0
3M™ Sponge-Stick with Lethen Broth	10 mL	SSL10LET	0.89 L (30 oz.)	0
3M™ Sponge-Stick with Lethen Broth, 2 Gloves	10 mL	SSL10LET2G	0.89 L (30 oz.)	2
3M™ Sponge-Stick with Neutralizing Buffer	10 mL	SSL10NB	0.89 L (30 oz.)	0
3M™ Sponge-Stick with Neutralizing Buffer, 2 Gloves	10 mL	SSL10NB2G	0.89 L (30 oz.)	2
3M™ Sponge-Stick with Neutralizing Buffer	10 mL	SSL10NB6X9	6" x 9" bag	0
3M™ Sponge-Stick with D/E Broth	10 mL	SSL10DE	0.89 L (30 oz.)	0

## INSTRUCTIONS FOR USE

1. Shake stick to end of bag.
2. Tear bag open.  
NOTE: Add sterile diluent or broth to dry sponge.
3. Squeeze bag to open. Aseptically grasp the stick above the thumb-stop line to remove the sponge.
4. Aseptically swab across the entire sampling surface<sup>3,4,5,6</sup>.
5. Turn sponge over. Change direction 90°. Aseptically swab the same sampled surface.
6. Aseptically place sponge into bag up to the thumb-stop.
7. Hold sponge in place inside bag. Bend stick to break. Allow sponge to drop in bag. Discard stick.
8. Fold bag to close.
9. Fold ends of blue wires inward.
10. Following user established procedures, remove any remaining enrichment broth or neutralizing solution residue from the sampled surface.

## REFERENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
4. US Food and Drug Administration - Bacteriological Analytical Method (available online at <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

Refer to the current version of the standard methods listed above.

## EXPLANATION OF SYMBOLS










Consult instructions for use

The lot in a box and the hourglass symbols are symbols that represent lot number and expiration date. The hourglass is followed by a year, month, and day which represents the expiration date (year, month, and day: YYYY- MM-DD). The entire line after the LOT represents the lot number (YYYY- MM AZ).

Storage conditions

## 3M Food Safety

### 3M United States

3M Center  
Bldg. 275-5W-05  
St. Paul, MN 55144-1000  
USA  
1-800-328-6553

### 3M Canada

Post Office Box 5757  
London, Ontario N6A 4T1  
Canada  
1-800-563-2921

### 3M Europe and MEA

3M Deutschland GmbH  
Carl-Schurz - Strasse 1  
D41453 Neuss/Germany  
+49-2131-14-3000

### 3M Latin America

3M Center  
Bldg. 275-5W-05  
St. Paul, MN 55144-1000  
USA  
1-954-340-8263

### 3M Asia Pacific

No 1, Yishun Avenue 7  
Singapore, 768923  
65-64508869

### 3M Japan

3M Health Care Limited  
33-1, Tamagawadai 2-chrome  
Setagaya-ku, Tokyo  
158-8583, Japan  
81-570-011-321

### 3M Australia

Bldg A, 1 Rivett Road  
North Ryde, NSW 2113  
Australia  
61 1300 363 878

### 3M Health Care

2510 Conway Ave  
St. Paul, MN 55144 USA  
[www.3M.com/foodsafety](http://www.3M.com/foodsafety)

© 2013, 3M. All rights reserved.  
3M is a trademark of 3M. Used under license in Canada.  
34-8711-5866-2

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## RÉFÉRENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. Exigences générales relatives aux compétences pour les analyses et les étalonnages en laboratoire.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3 : Microbiological Monitoring of the Food Processing Environment, 4th edition.
4. Food and Drug Administration (États-Unis) - Bacteriological Analytical Method (disponible en ligne en anglais à l'adresse suivante : <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04 et 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13 : Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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## RIFERIMENTI

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. Requisiti generali per la competenza dei laboratori di analisi e di calibrazione.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Capitolo 3: Microbiological Monitoring of the Food Processing Environment, 4a edizione.
4. US Food and Drug Administration - Bacteriological Analytical Method (disponibile online all'indirizzo <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
5. United States Department of Agriculture - Microbiological Lab Guidebook - Capitoli 4.04, 5.04 e 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Capitolo 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17a edizione.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
4. US Food and Drug Administration - Bacteriological Analytical Method (tilgængelig online på <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. Testaus- ja kalibrointilaboratorioiden pätevyys. Yleiset vaatimukset.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
4. US Food and Drug Administration - Bacteriological Analytical Method (käytettävissä osoitteessa <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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## REFERÊNCIAS

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. Requisitos gerais para a competência de laboratório de teste e calibração.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
4. US Food and Drug Administration - Bacteriological Analytical Method (disponível online em <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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- U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
- ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
- American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
- US Food and Drug Administration – Bacteriological Analytical Method (az interneten az alábbi címen elérhető: <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).
- United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
- American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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### REFERANSLAR

1. U.S. Food and Drug Administration. Code of Federal Regulations, Başlık 21, Bölüm 58. Good Laboratory Practice for Nonclinical Laboratory Studies.

2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.

4. US Food and Drug Administration - Bacteriological Analytical Method (available online at <http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>).

5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.

6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
4. US Food and Drug Administration - Bacteriological Analytical Method (オンラインでアクセス可能：  
<http://www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.html>) .
5. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
6. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

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