IKA

designed for scientists



Dispergierwerkzeuge/ Dispersion tools/ Outils dispersants

S 65 KD-G 65 G/M/F

S 65 KD-D 60/80

S 65 C-G 65 G/M/F

S 65 C-D 60/80

Technische Information	DE	3
Ursprungssprache		
Technical information	EN	8
Informations Technique	FR	13

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Content

	Page
Explication of warning symbols	8
Safety instructions	8
Product information	9
Application instruction	9
Spare parts list	9
Mounting the dispersing tool	10
Maintenance and cleaning	10
Warranty	11
Technical Data	12

Explication of warning symbols

-  **DANGER** Indicates an (extreme) hazardous situation, which, if not avoided, will result in death, serious injury.
-  **WARNING** Indicates a hazardous situation, which, if not avoided, can result in death, serious injury.
-  **CAUTION** Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
-  **NOTICE** Indicates practices which, if not avoided, can result in equipment damage.
-  **CAUTION** Burns hazard!

Safety instructions

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

 **WARNING** It is only allowed to use these tools on the drive **T65 basic/digital** with the safety switch **SI 400**, which make sure that the drive only starts when the tool is inside the vessel.

 **CAUTION** The dispersion tool may heats up during operation. In particular in the processing of hot medium, the tool must be cooled before it can be removed from the drive.

 **NOTICE** When in operation the dispersion tool can become hot. If the dispersion tool is not inserted into the drive flange correctly or the medium is hot, it can become extremely hot.

 **WARNING** Don't touch rotating parts during operation.

 **NOTICE** The vessel must be clamped, so that it cannot move during the operation. The **D 60** and **D 80** can destroy the vessel if the dissolver touches the vessel wall. The dissolver tools must have a distance of minimum 100 mm from the vessel wall.

 **CAUTION** Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There is a risk of splashing and evaporation of liquids; body parts, hair, clothing and jewelry getting caught.

- The work has to be interrupted immediately if you notice unusual noise and/or increasing emission of liquid from the upper, sidewise rinsing drilling. The bearings must be controlled, and replaced if necessary.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- Do not use any damaged dispersion tools.

 **NOTICE** Never run the dispersion tool dry. Without cooling by the medium the sealing will be destroyed.

- The dispersion tool is not suitable for continuous operation. Likewise the optimal dispersion duration and rotating frequency must be determined by attempts. Usually a few seconds are sufficient, maximum duration is one minute. Longer application times bring no improvements, increase however the sample temperature substantially.
- Dispersing causes the medium to heat up.
- Note the min./max. immersion depth in the "**Technical data**".
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator (integrated in the shaft tube). Abrasive media can round off these edges fast, whereby the disperse effect diminishes.

Product information

S65 KD: The shaft bearing of the dispersing tool consist of a ball bearing. The dispersing tool is sealed by means of a sealing ring.

S65 C: The shaft bearing of the dispersing tool consist of a ball bearing on the top and a ceramic slide bearing at the bottom.

All materials are FDA (Food and Drug Administration) conform.

Application instruction

Application: Wet crushing, dispersions, emulsifying (Batch-operation).

Spare parts list

S 65 KD:

Item	Designation
1	Plastic cap
2	Shaft
3	Shaft tube
4	Bearing
5	Sealing
6	Bushing
7	O-ring
8	Shaft sleeve
9	Coupling
11	O-ring
15	Nut
16	Dissolver
17	Rotor
18	Stator
19	Blade

S 65 C:


Item	Designation
1	Plastic cap
2	Shaft
3	Shaft tube
4	Bearing
6	Bushing
8	Shaft sleeve
9	Coupling
11	O-ring
13	Bushing
15	Nut
16	Dissolver
17	Rotor
18	Stator
19	Blade

Spare parts diagram see front page. For orders of spare parts, please specify the type of dispersion tool and the designation of the spare part.

Note to S 65 KD: During commissioning, after inserting a new sealing (5), can lead to increased wear.

Mounting the dispersing tool

WARNING The dispersion tool must only be mounted on the dispersion unit when the device is de-energized (e.g. by pulling the plug or pressing the mains switch) and the drive shaft is stationary.

- Open the handwheel bolt (**item 3**) so that the thread does not protrude into the hole.
- Insert the dispersing element (**item 2**) into the drive unit (**item 1**) as far as the stop. The dispersing element is fitted correctly if the recessed edge of the dispersing element is flush with the bottom edge of the drive unit.
- **Attention!** Now secure the dispersing element by screwing the handwheel bolt (**item 3**) tightly.

Maintenance and cleaning


Cleaning:

NOTICE The sterilization of the complete dispersion tool is **not** possible.


For precleaning the dispersion tool, operate it in a suitable cleaning fluid for dissolving the residuary of the substance. Due to the large flow rate, rotor and stator are cleaned to a large extent. The dispersion tool must be divided and cleaned immediately after working, so that the adhering residuary of the substance don't cultivate unwanted bacterial cultures.

For cleaning purposes disassemble the dispersion tool as described below (items see front page):


Tools for disassembling the dispersion tool:

Rotor wrench

Hook wrench


Shaft wrench

Note: The tools are included in the delivery of the drive.


Disassemble the rotor and dissolver:

G 65 G:


Insert the shaft wrench between the teeth of the stator.
Hold the shaft with the rotor wrench at the coupling section and unscrew the rotor in counterclockwise.

G 65 M/F:


Hold the stator with three pins of the shaft wrench.
Hold the shaft with the rotor wrench at the coupling section and unscrew the rotor in counterclockwise.

D 60/80:


Hold the dissolver with a wrench A/F of 17mm (**SW17**).
Hold the shaft with the rotor wrench at the coupling section and unscrew the rotor in counterclockwise.

Note: the user can also hold the shaft with the pin of the rotor wrench and unscrew the rotor with shaft wrench (this procedure will be necessary if the coupling (9) goes loose first.) or **SW 17** (for **D60/80**).

Disassemble the stator:

Hold the shaft tube with hand.
Unscrew the stator in counterclockwise (left hand thread) with the hook wrench.

The assembly of the dispersion tool takes place in reverse sequence.

NOTICE The rotor must not be tightened too strongly, otherwise the thread will be damaged.

Spare parts order:

Spare parts diagram see front page.

For orders of spare parts please specify the serial number on the shaft tube, the type of dispersion tool and the designation of the spare part.

Repair:

Please send dispersing tool for repair only after it has been cleaned and is free from any materials which may constitute a health hazard.

For repair, please request the "Decontamination Certificate" form **IKA®** or use the download printout of it from **IKA®** website: www.ika.com.

If you require servicing, return the instrument in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

Technical Data

		S65 KD/C-G 65G	S65 KD/C-G 65M	S65 KD/C-G 65F	S65 KD/C-D 60	S65 KD/C-D 80
Working range volumn	l	2–50	2–40	2–30	2–30	2–50
Stator / Rotor (Ø)	mm	65 / 53	65 / 58	65 / 58	- / 60	- / 80
Max. allowable speed	rpm	9500				
Max. circumferential speed	m/s	26.4	28.8	28.8	29.8	39.8
Min. / max. immersion depth	mm	90–450	80–450	80–450	90–450	90–450
Material in contact with medium		FFKM, PTFE, AISI 316L (for S 65 KD) Ceramic, AISI 316L (for S 65 C) (FFKM and PTFE conform to the FDA, Food and Drug Administration rules)				
pH range		2–13				
Suitable for solvents		yes				
Suitable for abrasive substance		no				
Sterilization methods		all methods				
Working range vacuum	mbar	100 (S 65 KD) / - (S 65 C)				
Temperature range upto	°C	180				
Ultimate fineness dispersions	µm	25–75	20–50	5–20	-	-
Ultimate fieness emulsions	µm	5–25	5–15	1–10	-	-

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IKA India Private Limited
IKA England LTD.
IKA Vietnam Company Limited
IKA Works (Thailand) Co. Ltd.