# IKA

designed for scientists

## IKA ETS-D6

Betriebsanleitung	DE	<b>6</b>
Ursprungssprache		
Operating instructions	EN	<b>14</b>
Mode d'emploi	FR	<b>22</b>
使用说明	ZH	<b>30</b>
Instrucciones de manejo	ES	<b>38</b>
Veiligheidsaanwijzingen	NL	<b>39</b>
Avvertenze per la sicurezza	IT	<b>41</b>
Säkerhetsanvisningar	SV	<b>42</b>
Sikkerhedshenvisninger	DA	<b>44</b>
Sikkerhetsveiledning	NO	<b>45</b>
Turvallisuusohjeet	FI	<b>47</b>
Instruções de serviço	PT	<b>48</b>
Wskazówki bezpieczeństwa	PL	<b>50</b>
Bezpečnostní pokyny a informace	CS	<b>51</b>
Biztonsági óvintézkedések	HU	<b>53</b>
Varnostni napotki	SL	<b>54</b>
BBezpečnostné pokyny	SK	<b>56</b>
Ohutusjuhised	ET	<b>57</b>
Drošības norādījumi	LV	<b>59</b>
Saugos nurodymai	LT	<b>60</b>
Указания за безопасност	BG	<b>62</b>


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5


Fig. 6


Fig. 7


Fig. 8


Fig. 9


Fig. 10


Fig. 11


Fig. 12


Fig. 13


Fig. 14


Fig. 17


Fig. 20


Fig. 23


Fig. 26


Fig. 29


Fig. 15


Fig. 18


Fig. 21


Fig. 24


Fig. 27


Fig. 16


Fig. 19


Fig. 22


Fig. 25


Fig. 28

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## Contents

	Page
EU Declaration of conformity	14
Safety instructions	14
Intended use	15
Unpacking	16
Commissioning	16
Connecting the pH-probe	16
Replacing extension cable and temperature sensor	17
Factory setting	17
Configuration	17
- Language	17
- Operating mode	17
- Maximum temperature	18
- Target temperature	18
- Temperature controller	18
- Time out	18
- pH-Measurement	18
- Reset	19
- Exit	19
- pH-Calibration	19
- PT 1000-Calibration	19
Error codes	20
Accessories	21
Warranty	21
Technical Data	21

## EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61326-1, EN 60529 and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at [sales@ika.com](mailto:sales@ika.com).

## For your protection

- Read the operating instructions in its entirety before using the device and follow the safety instructions. If this accessory is used with another device, observe also its operating instructions.
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the media to be processed.
- **Caution!** The supply voltage (Limited Power Supply) for the **ETS-D6** must be provided using an IEC 61558 compliant transformer with double or reinforced insulation. Alternatively, laboratory equipment that meets the requirements of IEC 61010 or NEC Class 2 may be used.
- **Caution!** In unsupervised and safe operation, this equipment must only be used with media that have a flash point higher than the temperature safety limit set for the magnet agitator and/or heating system used. The safety temperature limit must always be set to at least 25 °C lower than the flash point of the media used. (acc. to EN 61010-2-010)
- **Burn hazard!** Take care when touching the temperature sensor! The temperature sensor may heat up to 450°C when inserted in the media!
- Do not touch the temperature sensor while measurements are being taken. This will prevent incorrect results.
- Ensure that the external temperature sensor (PT 1000, ETS-D ...) is inserted in the media to a depth of at least 20 mm, see **Fig. 2**.
- Please ensure that the spiral cable does not come into contact with the heating plate.
- The stainless steel temperature sensor must not be used with aggressive media such as acids, caustic solutions or distilled water, due to the risk of corrosion. The glass sensor should be

used in such cases.

- Only use glass encapsulated temperature sensors for electrolysis procedures.
- Always use the extension cable H 70 when the media being processed produces vapour. This ensures that the control unit does not come into contact with the vapour.
- Beware of hazards due to:
  - flammable materials
  - combustible media with a low boiling temperature
- Do not use the device in explosive atmospheres, it is not EXprotected.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Check the instrument and accessories beforehand for damage each time you use them. Do not use damaged components.
- Safe operation is guaranteed only with the use of original IKA accessories.
- Accessories must be securely attached to the device and cannot come off by themselves. The centre of gravity of the assembly must lie within the surface on which it is set up.
- **Caution!** The user must check that all conditions have been met for the safe and correct operation of the equipment before operating the **ETS-D6** with a heater that does have a DIN 12878 compliant connector, and therefore does not provide cut-out protection. (The cut-out protection is a mechanism that ensures that the electrical heating circuit is switched off should the contact thermometer (**ETS-D6**) be broken or if the connector is removed.).

## For protection of the equipment

- Protect the instrument and accessories from bumps and impacts.

Use of "H+P" magnetic stirring systems (no break fuse protection included)

- When the supply current for the **ETS-D6** is interrupted, the electronics for the heatable magnetic agitator system will take over the temperature control. This means that the magnetic agitator system will heat up to the temperature set (e.g. 300 °C). Readjust the set-point temperature for the magnetic agitator system or switch the heating off.

## Use of "Corning" magnetic stirrer PC-400, PC-420, PC-600, PC-620

- When using the ETS-D6 in conjunction with a "Corning" Magnetic Agitator it is absolutely essential that the rotary "heat" dial on the magnetic agitator is set to the "O" (off) position. The **ETS-D6** will only control the set temperature when it is in this position.

**Warning:** The heating plate temperature limiter is not active in this operating mode.

- If the rotary "heat" dial on the magnetic agitator is not in the "O" position then the magnetic agitator will heat up to the temperature set, regardless of the settings on the **ETS-D6**.

**Warning:** **ETS-D6** is inactive.

The safety advices given above represent the current status of the known risks. Nonetheless, the user must check equipment for correct and safe operation before operating any system comprising the **ETS-D6** and a magnetic stirring/heating system.

## Intended use

## Use

- Precise measurement and control of temperature  
The instrument can be connected to any magnetic agitator or heating plate that has a contact thermometer connection, provided this conforms with the requirements listed in the "Technical Data", see chapter "Technical Data".

These are amongst others **IKA** magnetic stirrers and Heidolph magnetic stirrers. You can connect Corning magnetic stirrers by the adapter AD-C1.

- pH- measurement  
An optional pH-probe can be connected using the built in BNC connector.

When measuring temperature and pH at the same time, the glass sheathed temperature sensor should be used in order to avoid any interference with the pH measurement. To obtain useful results for pH measurements when using earthed heating plates (stainless steel, silumin) an electrically isolated container (glass beaker) must be used.

It is not possible to carry out pH measurements in metal containers.

The **ETS-D6** is not suitable for diagnostic measurements for medical purposes.

#### Area of use

- Indoor environments similar to that a laboratory of research, teaching, trade or industry area.
- The safety of the user cannot be guaranteed:
  - if the device is operated with accessories that are not supplied or recommended by the manufacturer,
  - if the device is operated improperly or contrary to the manufacture's specifications,
  - if the device or the printed circuit board are modified by third parties.

## Unpacking

#### Unpacking

- Unpack the device carefully. Any damage should immediately be reported to the carrier (mail, rail or freight forwarding company).

#### Scope of delivery

- ETS-D6
- Operating instructions.

## Commissioning

- ☞ Plug the instrument into the contact thermometer jack on the heating magnetic stirrer or heating device , e.g.

**IKA:** RCT basic, RH basic, RH digital  
**Heidolph:** MR3001K, MR3002  
**VWR:** VMS-D, VWR VMS-A  
**H+P** Variomag

Please observe all safety instructions in this respect.

Connecting a "Corning" magnetic stirrer:

- ☞ Turn the rotary "heat" dial on the magnetic agitator to the "O" (off) position. Next, connect the **ETS-D6** to the magnetic agitator using an accessory AD-C1 adapter.

Please observe all safety instructions in this respect.

- The device will be activated when the magnetic agitator is switched on.

If the supply voltage at the contact thermometer connector is not between 8 Vdc and 16 Vdc then it will be necessary to power the device using the H 52 power supply set (H 50 double connector adapter and H 51 plug-in power supply unit). The power supply and/or double connector adapter will provide the necessary supply voltage for the device. An automatic self-test will be performed when the device is activated. An introductory screen will be shown on the LCD display, see **Fig. 5**.

## Connecting the pH-probe

- ☞ Connect the BNC-connector for the pH-probe to the BNC-socket on the rear of the unit and lock in place using the bayonet fitting, see **Fig. 3**.

## Replacing extension cable and temperature sensor

- ☞ Pull the sensor and protective cap downwards, see **Fig. 4**
- The plastic snap fastener is released
- ☞ Push the sensor or extension cable and protective cap over the fitting on the instrument
- The connection is re-established
- LCD display see **Fig. 28**

## Factory setting

The **ETS-D6** is supplied with the following factory settings:

Maximum temperature	450 °C
Mode	A
pH-measurement	disabled
Controller	PID
pH-Calibration	OK
PT 1000-Calibration	-
Extension cable	-
Time Out Er 5	3 min
Language	English

## Configuration

The unit is configured using the Setup menu:

- ☞ Switch on the device
- LCD-Display: 1. Intro, see **Fig. 5**
  2. Information about the operating mode/ Maximum temperature and the request "OK to enter Setup" (Duration: 2 seconds), see **Fig. 6**
- ☞ Press the SET button during the two second period. Use the ▲ and ▼ buttons to navigate through the menu and select the desired menu item.

The SET button is used to activate a selected menu item and confirm changes, see **Figs. 7, 8, 9**.

The currently active settings are shown under the corresponding menu item.

Select the menu item "Exit" to leave the menu.

To accept and store changes select the "Yes" option when prompted with "Confirm changes" on the next screen.  
 ➢ LCD display see **Fig. 10**

## • Language

The following languages can be selected under the menu item "Language":

- English
- French
- Chinese
- Korean
- German
- Spanish
- Japanese

- LCD display see **Fig. 11**

Accept and store the selected language:

- ☞ Press button SET

## • Operating modes

The operating mode can be set under the menu item "Mode": The **ETS-D6** can be run in three different operating modes (A, B, C). The mode and maximum temperature currently set will be displayed for two seconds after power on, see **Fig. 6**.

#### Mode A

The set-point temperature is always set to -10°C.

The maximum temperature can only be adjusted if Mode A is active.

- ☞ Use the ▲ and ▼ buttons to change the target temperature.

#### Mode B

The last target temperature selected will be used.

- ☞ Use the ▲ and ▼ buttons to change the target temperature.

#### Mode C

The last values selected will be used.

☞ It is not possible to adjust the target temperature.

> LCD display see **Fig. 12**

Accept and store the selected mode:

☞ Press button SET

Change the mode

☞ Hold down the ▲ button when turning on the unit. The system will switch to the next operating mode.

Sequence: A-B-C-A-B-C-A ...

*Factory setting: Operating mode A*

#### • Maximum temperatur

The maximum temperature (maximum allowed target temperature) can be set using the "Max. Temp" menu item.

☞ Use the▲and▼buttons to change the maximum temperature.

(0 - 450 °C)

> LCD display see **Fig. 13**

Accept and store the selected maximum temperature:

☞ Press button SET

#### • Target temperature

The target temperature is set directly rather than using the Setup menu.

0,1 K steps

☞ Press buttons ▲ and ▼

1 K and 10 K steps respectively

☞ Press buttons ▲ and ▼ longer than 5 seconds

> LCD display see **Fig. 14, 15:**

- The instantaneous actual temperature with physical units

- The target temperature below the actual temperature

- The symbol  indicates heating is switched on

- The pH-value is shown at the bottom left if a pH-probe is connected

#### • Temperature controller

The type of control used can be set under menu item

"Controller":

PID: Gives best results with minimum overshoot

On/Off: Operates with maximum heating speed

☞ Select Controller

> LCD display see **Fig. 16**

Accept and store the selected controller:

☞ Press button SET

#### • Time out

"Time out" is the time taken to recognise that the sensor has not been placed in the medium when the heating is switched on.

**Caution!** This function will only be active if:

- Sensor temperature < 50 °C

- Difference target temperature/sensor temperature >5 K

The "Time out" period can be adjusted using the "Time-out Er5" menu item:

☞ Press buttons ▲ and ▼ (1 - 30 min)

> LCD display: "Time-out Er5" period see **Fig. 17**

Accept and store the selected "Time-out" period:

☞ Press button SET

#### • pH-Measurement

The pH measuring function can be activated using the "pH-measurement" menu item.

> LCD display see **Fig. 18**

Accept and store:

☞ Press button SET

#### • Reset (Restoring parameter to factory settings)

☞ Switch off the device

☞ Press and hold down the buttons ▲ and ▼, and switch on the equipment or use the "Reset" menu item

> LCD display see **Fig. 19**

Accept and store:

☞ Press button SET

#### • Exit

☞ Select the menu item "Exit" to leave the Setup menu.

> LCD display see **Fig. 9**

Any changes made can be accepted or rejected.

#### • pH-Calibration

The two-point calibration process can be started under menu item "ph calibration".

> LCD display see **Fig. 20**

☞ Calibration can be cancelled at any time using the ▲ and ▼ buttons.

Calibration is used to adjust the unit for the measurement chain being used. The process determines the zero point and the pH gradient for the measurement chain.

The calibration should be carried out using pH 4 and pH 7 buffer solution as per DIN 19266.

Other buffer solutions, including those with identical rated values, may exhibit different temperature characteristics, which will lead to measurement errors.

• First, immerse the measurement chain and the temperature sensor in the pH 7 buffer solution. The rated value for the required buffer solution will be displayed, see **Fig. 20**.

• Press the SET button to start the calibration process, see **Fig. 21** (duration up to 3 min.).

• Remove the measurement chain and the temperature sensor from the first buffer solution as soon as the rated value for the

second buffer solution appears on display (see **Fig. 22**). Rinse both items thoroughly with distilled water and dry using a soft, lint free paper towel. Do not rub!

• Next, immerse the measurement chain and the temperature sensor in the pH 4 buffer solution.

• Press the SET button to restart the calibration process, see **Fig. 23**.

The gradient for the measurement chain in % (mV/pH) will be displayed at the end of the calibration.

100%=>59,16 mV

The calibration process was successful if the gradient is between 90%-105%, see **Fig. 24**. The calibration values will be saved automatically.

• Confirm the setting using the SET button

If the gradient is outside this range then the calibration values will be rejected. The measurement chain should be swapped in this case.

#### • PT 1000-Calibration

The **ETS-D6** is pre-calibrated without a sensor before shipping. Sensor tolerances in accordance with DIN IEC 751 Class A are not taken into account.

The complete measurement chain can be calibrated in order to eliminate the effects of the contact resistance at the plug connector and the sensor tolerances (inclusive of temperature drift).

The calibration process is started using the "PT 1000 calibration" menu.

The following items are required:

- An additional temperature measuring device with an accuracy of + 0.05°C

- Two stable, constant temperature sources (media at different, controlled temperatures)

- Immerse both sensors in the first medium. Wait until the display on the temperature measurement device has stabilised.
- Enter the measured temperature on the **ETS-D6** using the ▲ and ▼ buttons.
- Confirm the value entered using the SET button, see **Fig. 25**.
- Immerse both sensors in the second medium. Wait until the

- display on the temperature measurement device has stabilised.
  - Enter the measured temperature on the **ETS-D6** using the ▲ and ▼ buttons.
  - Confirm the value entered using the SET button, see **Fig. 26**.
- The calibration is now successfully completed, see **Fig. 27**.

## Error codes

Any malfunctions during operation will be identified by an error message on the display, see e.g. fig. 29.

Proceed as follows in such cases:

- ☞ Disconnect the power supply for at least 10 sec.
- ☞ Carry out corrective measures
- ☞ Restart device

Error code	Cause	Effect	Correction
<b>Er 1</b>	Sensor not connected or defective (Sensor rupture)	Device off	- Connect temperature sensor - Replace temperature sensor
<b>Er 3</b>	Temperature inside device higher than 76 °C	Device off	- Switch off the device and allow to cool down
<b>Er 5</b>	Sensor not in media (no change in temperature of the media after heater active for 3 minutes)	Device off	- Insert the temperature sensor in the media to a depth of at least 20 mm. - Increase the "Time-out" period
<b>Er 6</b>	Temperature sensor short-circuit	Device off	- Switch off the device and replace temperature sensor
<b>Er 7</b>	The temperature of the medium is higher than the maximum target temperature that can be set	Device off	- Switch off the device and allow the media to cool down - Increase the maximum target temperature that can be set
<b>Er 8</b>	Calibration error	Device off	- Switch off the device - Re-calibration by authorised service personnel

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

## Accessories

- **H 38** Holding rod
- **H 16 V** Support rod
- **H 44** Boss head clamp
- **H 62.51** Standard sensor, Ø3 mm, length 260 mm, stainless steel
- **H 66.51** Glass-coated sensor, Ø 6 mm, length 260 mm
- **H 70** Extension cable, length 1m

See more accessories on [www.ika.com](http://www.ika.com).

## Warranty

In accordance with **IKA** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

## Technical Data

### Temperature

Range of temperature measurement	°C	-50 – 450
Resolution	K	0.01
Measuring accuracy	K	± 0.05 + PT 1000 Sensor tolerance, DIN IEC 751 Class A, ± ( 0,15 + 0,002 x  t  ), t in °C
Setting accuracy	K	0.1
Control deviation <small>(500 ml water in 600 ml glass beaker, magnetic stirring bar 40, 600 rpm, 50 °C)</small>	K	± 0.2

### pH-measurement

Range of measurement	pH	0...14
Accuracy	pH	± 0.1
Resolution	pH	± 0.01
pH-connection		BNC jack
Supply voltage	Vdc	8 – 16 (Limited Power Supply)
Current consumption	mA	30 (at 9V)
Perm. ambient temperature	°C	0 – 60
Perm. relative humidity	%	80
Perm. duration of operation	%	100
Protection type acc. to DIN EN 60529		IP 54
Length spiral cable	mm	350
Male connector		6-pin DIN 45322
Socket		DIN 12878 class 2
Dimensions (W x D x H)	mm	100 x 45 x 95 (without sensor)
Weight	kg	0.2
Operation at a terrestrial altitude	m	max. 2000

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**IKA Vietnam Company Limited**
**IKA India Private Limited**