OPERATING MANUAL

Cole-Parmer®

SP Overhead Mixer

Mixer Controller Model Number 55008-06

Mixer Motor Model Number 55008-08

SP Overhead Mixer Set

Model Number 55008-02

SP Overhead Mixer System

Model Number 55008-03

ENGLISH

Instruction Manual

pages 3 - 19

FRANÇAIS

Mode d'emploi

pages 23 - 40

ESPAÑOL

Manual de instrucciones

paginas 43 - 60

ACCESSORY TABLE

pages 62 - 63

CONTENTS

SAFETY PRECAUTIONS	4
GENERAL INSTRUCTIONS	7
INTRODUCTION	7
Packing List	8
Cord Sets	8
Motor	8
CONTROLLER AND LED FUNCTIONS	9
INSTALLATION AND SETUP PROCEDURES	10
Connecting Motor Cable to Controller	10
Installation of Adjustable Tilt Mounting Bracket	10
Installation Using a Mixer Stand	11
Installation of Pivot Mount	12
Installation Using a Tank Clamp	12
Impeller Attachment	13
Chuck Removal and Attachment	13
MIXER CONTROLLER OPERATION	13
FWD-OFF-REV Switch	13
SPEED Control	14
Overload Indicator	14
SPECIFICATIONS	14
Dimensions	14
Mixer Motor	15
Controller	15
Torque Capabilities	16
Torque Curve	16
Compliance	16
MAINTENANCE AND TROUBLESHOOTING	17
Replacement Parts	17
Troubleshooting	18
WARRANTY	19
PRODUCT RETURN	19
TECHNICAL ASSISTANCE	19

SAFETY PRECAUTIONS

1. All operators must be familiar with the mixer and read this entire manual. If the mixer is used with accessories not provided or recommended by the manufacturer or used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.
2. The mixer must be securely fixed to a stable support. Mount to a Cole-Parmer stand using a Cole-Parmer double (rod) clamp. If other stands or clamps are used, the mixer must be secured such that it will not fall if tipped.
3. **CAUTION:** This is not an explosion proof mixer. Do not use with highly flammable or explosive materials. Do not operate the mixer in explosive atmospheres, with hazardous substances or under water.
4. Spinning impellers can cause severe personal injuries. Operators must use extreme care and good judgment when mixing at any speed.
5. All mixing paddles and impellers must be in good condition with straight shafts. If the mixer vibrates at high speeds check the paddle shaft for damage and repair or replace if needed.
6. Extreme care must be taken when mixing chemicals to ensure that no chemicals are splashed outside the mixing vessel. Care must be taken when changing to faster mixing speeds.
7. Ensure that the mixing impeller does not contact the containment vessel.
8. Do not operate while standing in water. Keep the mixer dry and do not immerse any part, except the mixing paddle, into any liquids. Protect from splashing.
9. Ensure that no loose clothing, jewelry, or hair can become entangled in any rotating parts. A fast spinning chuck can cause injury to operator.
10. To stop the mixer in an emergency, turn the SPEED control knob counter-clockwise to the "0" position, place the FWD-OFF-REV switch in the OFF position, push the rear-panel POWER switch to the OFF position, or remove the power plug disconnect the mains cord.
11. Ensure the SPEED control knob, POWER switch, and mains power cord of the mixer can be accessed immediately, directly and without risk at any time, and above all, in case of fault.
12. The mixer must be used with the supplied power cord. The cord and grounded plug type are important safety elements. Always connect the mixer to a grounded outlet. Do not alter the plug.
13. The mixer must be supplied with the voltage specified on the serial label. See serial label.
14. Excluding procedures outlined under maintenance and troubleshooting, service must be done by an authorized service centre.
15. Be sure the SPEED control knob is turned fully counter clockwise before the FWD-OFF-REV switch is placed in the FWD or REV positions. Gradually increase speed to desired level.
16. Wear appropriate personal protective equipment such as eye protection and any clothing suitable for the application.
17. When the mixer is ready for use to mix liquids and liquid/solid solutions, the operator must ensure that all safety instructions are followed and must determine appropriate impeller and container specifications.

18. The mixer can operate continuously if the torque demand falls below the continuous torque operation line (see Torque Capabilities).
19. If the mixer is used for an extended period of time above the continuous torque operation line (see Torque Capabilities), then it may shut down. The mixer may be started after it has cooled down.

GENERAL INSTRUCTIONS

Read Safety Precautions. Unpack the product carefully. Ensure the product is not damaged and verify the contents against the packing list. If the product is damaged or the contents do not correspond with the packing list, notify your dealer immediately.

Read the manual carefully! Make sure that all users read this manual thoroughly before operating the mixer. For further questions or inquiries, please contact:

Cole-Parmer

625 East Bunker Court
Vernon Hills, IL 60061 United States

Tel: 847-549-7600
Toll Free: 1-800-323-4340

Website: www.coleparmer.com
E-mail: sales@antylia.com

INTRODUCTION

SP Overhead Mixer covered in this manual is specifically designed for professional use in the laboratory and for use in industry to stir numerous liquids with various properties and applications. It consists of a mixer motor and an analog controller.

SP Overhead Mixer rotates at 50-2500 rpm with a maximum torque of 91 in-oz (64 Ncm).

The mixer uses a brushless DC motor. The speed is kept constant with continuous feedback to the motor. The mixer has integrated overload protection which ensures no damage can occur to the electronics or when an overload situation occurs, such as too great a torque demand or mixer overheating.

The mixer is designed for reliability and durability. There are no wear parts that require user service during the normal life of the mixer.

Packing List

Compare the contents of the shipment to the items on the following lists to ensure all parts are received. Do not discard the packing materials until all parts are accounted for. It is recommended that the original packing materials be kept for transportation and storage.

SP Overhead Mixer Set 55008-02 covered in this manual consists of:

- Mixer Motor Model 55008-08 including:
 - Bracket with 6 in (152 mm) mounting rod
 - Hardware
- SP Mixer Controller Model 55008-06 including:
 - North American cord set

SP Overhead Mixer System 55008-03 consists of the set above plus:

- Mixer stand
- Double (rod) clamp
- Stainless steel shaft for propeller, 18 in x 3/8 in (457 mm x 10 mm)
- Stainless steel propeller, 3 in (78 mm)

Cord Sets

North American cord set included. RoHS compliant. For international cord sets contact your dealer.

Motor

The Cole-Parmer SP and SD mixer motor, Figure 1, has through-shaft capability and features a user replaceable keyless chuck. Bracket hardware is preinstalled to enable installation of the adjustable tilt mounting bracket.


Figure 1: Mixer Motor with Cable

CONTROLLER AND LED FUNCTIONS

SP Mixer Controller 55008-06

The variable-speed controller, Figures 2 and 3, is used to turn on and set the speed and direction of the mixer motor.


Figure 2: SP Controller Functions

The controller ensures accurate, dependable motor speed (mixing) control, with “soft-start” and smooth, gradual adjustment to the speed setting. Constant speed is maintained electronically via automatic compensation for changes in load torque demands at both high and low speeds. Built-in motor current compensation minimizes speed drift with changes in load. Motor current limiting protects the system in case of a locked rotor condition. Output circuit protection protects against inadvertent output shorting. Mixer rotation is reversible.

When LED is green it indicates the controller is powered on. When LED flashes, it indicates a limit condition has been reached.

The Forward-Off-Reverse Switch is used to set the direction of rotation or to stop motor rotation in the OFF position. The motor and controller are still powered on (as indicated by the LED) if the switch is in the OFF position.

The SPEED Control knob sets the speed of the mixer motor.


Figure 3: SP Controller Back Panel

POWER Switch controls power to controller and motor.

C14 Connector is for power cable to mains.

INSTALLATION AND SETUP PROCEDURES

Connecting Motor Cable to Controller

The Motor Cable connector has a locking mechanism when cable is connected. This is to prevent unintended disconnection in the event of accidental force on the motor cable.

To ensure proper connection:

1. Orient the motor cable connector so the unlock icon and arrows are facing upward.
2. Push the connector into the controller socket. The collar should rotate and return to its initial position, locking in the connection.
3. To disconnect the motor cable, rotate the blue collar counter-clockwise (left) to release the lock. The connector can now be pulled out of the socket, as shown in Figure 4.


Figure 4: Connecting and disconnecting Motor Cable to Controller

Installation of Adjustable Tilt Mounting Bracket

The mixer motor has an adjustable tilt mounting bracket included. The mounting bracket is designed to attach to a mixer stand or tank clamp.

Perform the following steps:

1. The motor comes with two black Phillips screws installed in the housing. Loosen the screws without removing in preparation for installing the mounting bracket.
2. Align the holes on the mounting assembly with the loosened screws. Place openings over the screws and slide the mounting assembly into the installed position.
3. Tighten the screws to secure the mounting assembly in the installed position, as shown in Figure 5.


Figure 5: Installing Adjustable Tilt Mounting Bracket

The mixer is now ready to be attached to a mixer stand (supplied with system) or to a tank clamp. Recommended stands, tank clamps, and a double (rod) clamp are available through your dealer. Part numbers are listed in the ACCESSORIES section.

Installation Using a Mixer Stand

1. Set up mixer stand in accordance with the manufacturer's instructions.
2. Attach a double (rod) clamp to the vertical rod on the mixer stand and tighten clamp screw.
3. Slide rod of adjustable tilt mounting bracket into clamp and secure with clamp screw.
4. Place mixer stand legs around vessel. Tilt angle is adjustable from 0 to 30 degrees.
5. Position the height of the mixer motor by moving the clamp up or down on the mixer stand. Once installed, the impeller can also be adjusted to the desired height.


Figure 6: Typical Mixer Setup using a Mixer Stand and Pivot Mount

Note: The controller can be placed on any flat surface within 6 ft (183 cm) of the mixer motor or can be mounted on a pivot mount (not supplied). Refer to the ACCESSORIES section for part numbers.

To attach the pivot mount, proceed to the section below. If controller is placed elsewhere proceed as follows:

6. Connect the 6 ft (183 cm) cable between the mixer motor and the controller rear panel.
7. Check that the controller is not powered. Ensure that the rear panel POWER switch is in the OFF position.
8. Connect the 6 ft (183 cm) line cord from the controller to the appropriate 3-wire, grounded power source.

Installation of Pivot Mount

1. Remove the two Phillips screws and lock washers on the bottom of controller.
2. Position Pivot Mount bracket over the bracket mounting holes, re-install the two Phillips screws and lock washers and tighten.
3. Attach a double (rod) clamp to the vertical support rod on the stand and position so that pivot mount bracket will be clear of mixer motor.
4. Slide the mounting rod of the pivot mount bracket into the double (rod) clamp or lower hole on the tank clamp. Tighten clamp.
5. Adjust position of controller by first loosening knob on pivot mount bracket, adjust position and then retighten knob.
6. Connect the 6 ft (183 cm) cable between the mixer motor and the controller rear panel.
7. Check that the controller is not powered. Ensure that the rear panel POWER switch is in the OFF position.
8. Connect the 6 ft (183 cm) line cord from the controller to the appropriate 3-wire, grounded power source.
9. Excess cordage from either the mixer motor or controller can be wrapped around the ears of the Pivot Mount to prevent entanglement with the rotating mixer shaft.

Installation Using a Tank Clamp

The tank clamp is attached to the side of a tank and supports both the mixer motor and the controller (see Figure 7).

1. Position the tank clamp over the edge of the tank with large knob to the outside of tank.
2. Allow clamp to rest on edge of tank for added stability and tighten knob.
3. Position mixer motor over material to be mixed and insert rod of adjustable tilt mounting bracket into upper hole in tank clamp. Tighten locking knob securely.
4. Select an impeller and shaft for the required operation and install as instructed below. Contact your dealer for obtaining the proper shaft size. Refer to the ACCESSORIES section for the part numbers.

Note: The controller can be placed on any flat surface within 6 ft (183 cm) of the mixer motor or can be mounted on a pivot mount (not supplied). Refer to the ACCESSORIES section for part numbers. To attach the pivot mount, proceed to the section above.


Figure 7: Typical Mixer Setup using a Tank Clamp.

Impeller Attachment

Select an impeller and shaft for the required operation and install as follows:

1. Place impeller on shaft and tighten set screw.
2. Open the chuck. Insert the impeller shaft into the open chuck and adjust the impeller to the desired height. The hollow through-shaft allows the impeller shaft to come through the top of the mixer motor if necessary. Insert the mixer shaft into the chuck so it extends into the through-shaft by at least 3 inches to ensure a secure hold on the shaft. Keep mixer shaft as short as possible to reduce resonant vibrations caused by shaft flexing at higher speeds.
3. Tighten the chuck. A clicking sound will be heard as the chuck approaches being fully tightened.

Chuck Removal and Attachment

1. Slide chuck onto motor shaft, aligning the set screw with the through hole on the shaft.
2. Tighten the set screw using the provided hex key.

MIXER CONTROLLER OPERATION

Powering On

Make sure that the motor cable is plugged and locked into the controller and the controller is connected to the mains supply with provided power cable. Place the rear-panel POWER switch to the ON position. If the controller is connected to mains, the LED will illuminate green indicating power. This is regardless of the position of FWD-OFF-REV switch.

FWD-OFF-REV Switch

In the FWD position, the shaft will rotate in a clockwise direction when viewed from the motor down the mixing shaft.

In the center or OFF position the shaft rotation is stopped but power is still connected to the controller and motor.

In the REV position, the shaft will rotate in a counter-clockwise direction when viewed from the motor down the mixing shaft.

If the following conditions apply:

- The FWD-OFF-REV switch is in FWD or REV positions
- The SPEED Control knob is not at the "0" position
- The controller remains plugged into AC power

The mixer will resume mixing after a power failure, or when the power has been removed by control equipment (i.e., PLC), or when the AC power has been cut (i.e., the mixer unplugged). When power is restored, the LED will flash 5 times and then the mixer will automatically ramp up to the last speed. This function cannot be enabled or disabled by the end user; it is always enabled.

SPEED Control

The SPEED Control precisely sets the output speed of the shaft. Speed can be as low as 50 rpm or as high as 2500 rpm.

Slowly adjust the SPEED Control from the “0” position to a setting that provides the desired mixing speed. The controller will ramp up the mixing speed to match the setting of the control.

The table below outlines the typical speeds for each speed setting

Setting	3	4	5	6	7	8	9	10
Speed rpm	500	750	1000	1500	1750	2000	2300	2500

Always return the SPEED Control to the “0” position when mixing is complete.

Note: With new mixers, allow approximately 20 minutes run-in time to allow warm-up drift to become minimal.

Overload Indicator

Motor overload can result from too much load (torque) on the mixer, a locked motor condition, or operating conditions causing overheating in the motor or controller. When the mixer is approaching an overload condition, the speed of the motor may automatically slow in order to reduce the load. If the overload situation continues, the motor will stop rotation and the green LED indicator will flash.

Resolve the condition causing too much load or a locked motor.

The mixer will remain stopped until it is reset by turning power off and back on again. If the shutdown occurred due to overheating, the mixer will not turn back on until sufficient time has passed for the mixer to cool. Turn mixer off, wait approximately 30 minutes, then turn mixer back on to resume mixing.

SPECIFICATIONS

Dimensions


Figure 8: Mixer Motor and SP Controller Dimensions

Mixer Motor

Direction of Rotation	Clockwise or counter-clockwise controlled with the FWD-OFF-REV switch
Speed Range	50 to 2500 rpm
Torque Output max At 500 rpm	54 in-oz (38 Ncm) Continuous 91 in-oz (64 Ncm) Intermittent
Maximum Volume	25 gal (95 L) [250 cP]
Maximum Viscosity	1000 cP [< 1 gal (4 L)]
Chuck	Three-jawed type chuck accepts up to 3/8 in (10 mm) diameter shafts
Through Shaft diameter	3/8 in (10 mm) maximum
Mechanical Output Power	1/15 hp (50 W) brushless DC Motor
Input Voltage max	24 V DC
Input Current max	8.4 Amps
Installation Category	Installation Category I per IEC 664
Motor Cord Length	6 ft (183 cm) - DC
Weight	6.6 lb (3.0 kg)

Controller


Speed Accuracy	± 2% Full Speed
Voltage	100-240 VAC, 50/60 Hz +/- 10% of the nominal voltage
Maximum Current Draw	2.6 Amps
Electrical Output to Motor	DC 24 Volts, 8.40 Amps, 201.6 Watts
Installation Category	Installation Category II per IEC 664
NA Power Cord Length	6 ft (183 cm) - AC
Weight	3.7 lb (1.7 kg)

Torque Capabilities

	Continuous	Intermittent
Model 55008-02	54 in-oz (38 Ncm)	91 in-oz (64 Ncm)

Torque Conversion Table		
in-lb to in-oz	Multiply by 16	A in-lb = B in-oz = C Ncm
in-lb to Ncm	Multiply by 11.36	
Ncm to in-oz	Multiply by 1.416	

Torque Curve

Compliance

Compliance of mixer system includes both the mixer motor and the controller.

Enclosure Rating (IP)	IP44 per IEC 60529
Chemical Resistance	Exposed material is anodized aluminum, ABS, and polyester.
Compliance	CAN/CSA C22.2 No. 61010-1:2019 CAN/CSA-C22.2 No. 61010-2-051:2015 UL 61010-1:2019 UL 61010-2-051:2015 IEC 61010-1:2010, IEC 61010-2-051:2015 EN 61010-1:2010, EN 61010-2-051:2015 FCC Part 15 Subpart B, Class A ICES-003:2020 IEC 61326-1:2012 CISPR 11:2017 EN 61326-1:2013, EN 55011:2016

Environment	<p>Indoors</p> <ol style="list-style-type: none"> 1. Altitudes up to 6500 ft (2000 m) 2. Temperatures from 41°F to 104°F (5°C to 40°C) 3. Maximum relative humidity 80% for temperatures up to 88°F (31°C) decreasing linearly to 50% relative humidity at 104°F (40°C) 4. Pollution Degree 2 in accordance with IEC 664
--------------------	--

RoHS	Restriction of Hazardous Substances (EU)
WEEE	Waste Electrical and Electronic Equipment (EU)
CSA	Canadian Standards Association
CE	Communauté Européenne (European Community)
UKCA	UK Conformity Assessed
UL	Underwriters Laboratories
FCC	Federal Communications Commission
IEC	The International Electrotechnical Commission
ICES	Interference Causing Equipment Standards
EMC	Electromagnetic Compatibility
EMI	Electromagnetic Interference

MAINTENANCE AND TROUBLESHOOTING

Maintenance is limited to cleaning the unit. If a fault is detected when troubleshooting and cannot be corrected by the suggested remedy, contact your dealer for service.

Cleaning the Mixer/Controller System

Keep the mixer housing and controller enclosure clean with mild detergents. Never immerse or use excessive fluid.

Replacement Parts

No user-serviceable parts are inside of this instrument. Refer servicing to your dealer.

The following parts are user-replaceable:

Description	Part Number
Mixer Keyless Chuck replacement - SP/SD	55008-10
Adjustable Tilt Mounting Bracket	55008-11

Troubleshooting

Symptom	Cause	Remedy
<p>POWER switch is ON.</p> <p>Mixer shaft stops or does not rotate when SPEED Control is turned past "0".</p> <p>LED Indicator is not illuminated.</p>	<ul style="list-style-type: none"> • Disconnected AC line cord. • Defective line cord. • Defective internal power supply. • Defective control. 	<ol style="list-style-type: none"> 1. Check controller is plugged into a live line. 2. Check cord connections to the mains supply and controller. 3. Contact dealer for service.
<p>POWER switch is ON.</p> <p>Mixer shaft stops or does not rotate when SPEED Control is turned past "0".</p> <p>LED Indicator is not illuminated.</p>	<ul style="list-style-type: none"> • Controller overload. 	<ol style="list-style-type: none"> 1. Turn controller POWER OFF, then ON again to reset. 2. If controller LED is not illuminated after reset, turn controller POWER OFF. Reduce environmental temperature if applicable. 3. Wait approximately 30 minutes to allow the controller to cool down. 4. Turn controller POWER ON again to resume mixing.
<p>POWER switch is ON.</p> <p>Mixer shaft stops or does not rotate when SPEED Control is turned past "0".</p> <p>LED Indicator is illuminated.</p>	<ul style="list-style-type: none"> • Defective motor cable. • Defective mixer motor or controller. 	<ol style="list-style-type: none"> 1. Check motor cable is properly connected to mixer controller. 2. Contact dealer for service.
<p>POWER switch is ON.</p> <p>Mixer shaft stops or does not rotate when SPEED Control is turned past "0".</p> <p>LED Indicator is flashing.</p>	<ul style="list-style-type: none"> • Disconnected motor cable. • Defective motor cable. • Defective mixer motor or controller. • Motor overheat caused by long running time at high loads. 	<ol style="list-style-type: none"> 1. Check motor cable is properly connected to mixer controller. 2. Turn controller POWER OFF, then ON again to reset. 3. If Mixer did not reset, turn controller POWER OFF again. Reduce motor loads and/or environmental temperature. 4. Wait approximately 30 minutes to allow the mixer to cool down. 5. Turn controller POWER ON again to resume mixing. 6. If LED Indicator continues to flash, contact dealer for service.
<p>POWER switch is ON.</p> <p>Mixer shaft stops or does not rotate when SPEED Control is turned past "0".</p> <p>LED Indicator is flashing.</p>	<ul style="list-style-type: none"> • Motor overload caused by too much load (torque) on the motor or by a locked motor condition. 	<ol style="list-style-type: none"> 1. Turn controller POWER OFF. 2. Resolve the condition causing too much load or a locked motor. 3. Turn controller POWER ON again to resume mixing.

WARRANTY

The manufacturer warrants this product to be free from significant deviations from published specifications. If repair or adjustment is necessary within the 2 year warranty period, the problem will be corrected at no charge if it is not due to misuse or abuse as determined by the manufacturer. Repair costs outside the warranty period, or those resulting from product misuse or abuse, may be invoiced. To register the warranty online, go to the Cole-Parmer website address below:

www.coleparmer.com/product-warranty-registration

PRODUCT RETURN

To limit charges and delays, contact your dealer for authorization and shipping instructions before returning the product, either within or outside of the warranty period. Clean the motor and controller so that it is free from any materials which may be harmful. Provide an SDS or other safety data sheets for all chemicals used in conjunction with this system where appropriate. When returning the product, please state the reason for the return. For your protection, pack the product carefully, add labels and follow shipping directions provided. Insure it against possible damage or loss. Any damages resulting from improper packaging are your responsibility.

TECHNICAL ASSISTANCE

If you have any questions about the use of this product, contact your local dealer.

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Cole-Parmer®

ACCESSORIES

The accessories listed below, used for mounting the mixer head and controller, are available from your dealer. Many other accessories such as shafts and impellers are also available. Contact your dealer for the full line of available accessories.

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	Part Number	Description
	55008-06	Cole-Parmer SP Controller 91 in-oz (64 Ncm) max torque 25 gal (95 L) [250 cP] 1000 cP [< 1 gal (4 L)]
	55008-07	Cole-Parmer SV Controller 113 in-oz (80 Ncm) max torque 31 gal (120 L) [250 cP] 2000 cP [< 1 gal (4 L)]
	55008-08	Cole-Parmer SP/SV Mixer Motor 50-2500 rpm Keyless chuck 6 in (152 mm) mounting rod with tilt bracket
	08041-22	Double (rod) Clamp
	04552-70	Pivot Mount for Controller, includes clamp
	04552-80	Mixer Stand, 19-1/4 x 11 in (489 x 279 mm) base with 5/8 in dia x 28 in (16 x 711 mm) stainless steel rod
	04552-65	Tank Clamp
	04552-25	Shaft, Stainless Steel, 3/8 in dia x 12 in (10 x 305 mm) • Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers
	04552-35	Shaft, Stainless Steel 3/8 in dia x 18 in (10 x 457 mm) • Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers
	04552-00	Shaft, Stainless Steel, 3/8 in dia x 24 in (10 x 610 mm) • Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers
	04552-05	Shaft, Stainless Steel, 3/8 in dia x 30 in (10 x 762 mm) • Accepts 5/16 in or 3/8 in (8 mm or 10 mm) bore propellers

	04552-20	Shaft, Stainless Steel 5/16 in dia x 12 in (8 x 305 mm)
	04552-30	Shaft, Stainless Steel, 5/16 in dia x 18 in (8 x 457 mm)
	04552-40	Three-blade Propeller, Stainless Steel, 1 in dia x 5/16 in bore dia (25 x 8 mm)
	04552-45	Three-blade Propeller, Stainless Steel, 1.5 in dia x 5/16 in bore dia (38 x 8 mm)
	04552-50	Three-blade Propeller, Stainless Steel, 2 in dia x 5/16 in bore dia (50 x 8 mm)
	04552-60	Three-blade Propeller, Stainless Steel, 3 in dia x 5/16 in bore dia (78 x 8 mm)
	04552-55	Three-blade Propeller, Stainless Steel, 2 in dia x 3/8 in bore dia (50 x 10 mm)
	04552-10	Turbine Propeller, Stainless Steel, 2 in dia x 5/16 in bore dia (50 x 8 mm)
	04552-15	Turbine Propeller, Stainless Steel, 2 in dia x 3/8 in bore dia (50 x 10 mm)
	04552-61	Straight Blade impeller only 1-9/16 in dia x 5/16 in bore dia (40 x 8 mm)
	04552-62	Dispersion Blade impeller only 1-7/8 in dia x 5/16 in bore dia (48 x 8 mm)
	50801-53	Square Blade impeller; 2-3/4 in (70 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)
	50801-56	Cross Blade impeller; 2 in (50 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)
	50801-58	Collapsible blade impeller; 3-1/2 in (90 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)
	50801-60	Straight blade impeller; 2 in (50 mm) dia with 5/16 x 16 in shaft (8 x 406 mm)
	50801-62	Pitched blade impeller; 2-1/2 in (64 mm) dia with 5/16 x 16 in shaft (64 x 406 mm)

Cole-Parmer®

625 E. Bunker Court
Vernon Hills, IL 60061 United States

847-549-7600

Toll Free: 1-800-323-4340 (USA)

sales@antylia.com

www.coleparmer.com

Canada

1-514-355-6100

info@antylia.ca

www.coleparmer.ca

United Kingdom

+44 (0) 1480 272279

uk.sales@antylia.com

www.coleparmer.co.uk

Germany

+49 (0) 9377 92030

de.sales@antylia.com

www.coleparmer.de

France

+33 (0) 1486 37800

fr.sales@antylia.com

www.coleparmer.fr

Italy

+39 (0) 2843 49215

it.sales@antylia.com

www.coleparmer.it

India

+9122 61394444

info@coleparmer.in

www.coleparmer.in

China

+1 847-549-7600

sales@antylia.com

www.coleparmer.com