
# Temperature Controller SCT1

## Instruction Manual Version 2.2

# English

Figure 1 - Front view


Figure 2 - Rear view


Figure 3 - Hotplate front view



Figure 4 - Hotplate rear view

# Temperature Controller

## SCT1 Instructions for use

### Introduction

Thank you for purchasing this piece of Stuart equipment. To get the best performance from the equipment and for your own safety, please read these instructions carefully before use. Before discarding the packaging check that all parts are present and correct.

### General Description

The Stuart SCT1 temperature controller is the ideal instrument for the accurate temperature control of aqueous and oil based samples in the laboratory. The SCT1 can be used in two different modes; as a precise temperature controller from 20 to 200°C or as a digital thermometer from -4 to 325°C. With digital setting and display of temperature, the SCT1 temperature controller is accurate to  $\pm 0.5^\circ\text{C}$ .

The SCT1 is compatible with the following Stuart hotplates and hotplate stirrers: US150, US152, US152D, UC150, UC152, UC152D, SC162.

### Safety Advice before use

**If the equipment is not used in the manner described in this manual and with accessories other than those recommended by Cole-Parmer the protection provided might be impaired.**

This equipment is designed to operate under the following conditions: -

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range  $+5^\circ\text{C}$  to  $+40^\circ\text{C}$
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Use with a minimum distance all around of 200mm from walls or other items.
- ❖ Never move or carry the unit when in use or connected to the hotplate.
- ❖ Do not use in a hazardous atmospheres or with hazardous materials

### Preparation for use

To use the SCT1 temperature controller with either a metal top or ceramic top hotplate, place the hotplate on a firm, level surface.

Open the jaws of the retort clamp on the rear of the SCT1 temperature controller by turning the clamp knob anticlockwise. The clamp can be used with retort rods from 9.5mm to 13mm diameter. Slide the SCT1 on to the retort rod and tighten by turning the clamp knob clockwise.

#### **NB: DO NOT OVER TIGHTEN**

Connect the SCT1 controller plug to the DIN probe socket at the rear of the hotplate (Figures 1 and 4).

Connect the probe plug to the probe socket at the rear of the controller (Figures 2 and 4). To ensure the best connection, tighten the locking collar on the plug by turning clockwise.

Connect the mains supply lead to the IEC socket at the rear of the hotplate and switch on at the mains. Where the unit has a mains on/off switch, switch the unit on using this switch.

To switch on the SCT1 temperature controller, turn the heat control knob on the hotplate clockwise (Figure 3).

#### **NB: The controller power is supplied by the hotplate**

The SCT1 temperature controller and hotplate are now ready for use. Correct communication between the SCT1 temperature controller and the hotplate is indicated by the illuminated amber contact thermometer indicator on the front panel of the hotplate.

**NB.: The hotplate and controller will perform an automatic self-test at switch on, before the current temperature of the probe is displayed. The factory default set temperature is 20°C.**

# Operation

## Controller positioning

The SCT1 temperature controller is supplied with a detachable temperature probe (Figure 2). This useful feature enables the SCT1 temperature controller to be positioned away from damaging fumes.

The temperature controller is factory set such that the retort clamp is positioned for mounting on a vertical retort rod. To mount the SCT1 temperature controller on to a horizontal retort rod, ease back the clamp situated at the rear of the controller sufficiently to clear the right-angled raised location. Rotate clockwise by 90° until the clamp snaps back into the location.

## Probe positioning

For optimum temperature control, ensure the end of the sensing probe is immersed at least 20mm deep into the medium being controlled. The temperature probe when not in use should be placed in the probe holder located on the left hand side of the SCT1 control panel. For safety and to avoid contamination, always wipe the probe clean prior to placing into holder.

## Using the SCT1 temperature controller

The SCT1 can be used in two different modes, either as a temperature controller ('control mode') or as a digital thermometer ('monitor mode').

### To use the SCT1 in control mode

Ensure that the SCT1 temperature controller is connected to the hotplate. Also ensure the temperature probe is attached to the controller and immersed in the medium to be controlled before switching on at the mains.

### US150, US152, UC150 and UC152


Turn the unit on using the On/Off switch at the side of the unit. The SCT1 will power up in monitor mode. To turn on the temperature controller in control mode, turn the heat control clockwise on the front panel of the hotplate.

### US152D and UC152D

Turn the unit on using the On/Off switch at the side of the unit. The SCT1 will power up in the last mode used. To change from monitor mode to control mode, press either the '+' or '-' key on

the SCT1. To turn on the temperature controller in control mode, turn the heat control clockwise on the front panel of the hotplate; the LED display will show "---".

The SCT1 temperature controller in control mode is indicated by the illuminated amber contact thermometer indicator on the front panel of the hotplate (Figure 3). The LED and display will be visible on the SCT1 temperature controller panel (Figure 5) and the current temperature of the medium will be displayed.

**Figure 5 - SCT1 temperature controller panel**

**NB: Turning the heat control clockwise on the hotplate turns the SCT1 temperature controller into control mode and cannot be used to set the hotplate temperature.**

### To set the temperature in controller mode

Press the + or - key once to instigate the temperature set point mode. This is indicated by a flashing LED on the SCT1 temperature controller display. To select the temperature, use the '+' or '-' keys to obtain the desired temperature. To set the temperature, release the key. The LED will cease flashing once the temperature has been set. With the probe in the sample, the thermometer LED will illuminate as follows:

- ❖ Green – temperature set point not reached
- ❖ Amber – liquid at set point
- ❖ Red – liquid hotter than set point

**NB: To obtain accurate temperature control, the minimum temperature set point must be at least 15°C above the surrounding ambient temperature.**

### To use the SCT1 in monitor mode *US150, US152, US152D, UC150, UC152 and UC152D*

To use the SCT1 controller in monitor mode together with the hotplate, ensure the unit is switched off using the mains control switch at the side. Hold down both the '+' and '-' keys simultaneously on the SCT1 temperature controller and turn the mains switch to 'ON'. Once the digital display illuminates, release the keys on the SCT1 temperature controller. Monitor mode will be indicated by a red flashing LED at the bottom right hand corner of the display of the SCT1 temperature controller. The amber control LED on the front panel of the hotplate will NOT be illuminated in monitor mode.

**NB: When in monitor mode, the SCT1 temperature controller operates as a digital thermometer only, and no longer controls the temperature of the medium in use. The hotplate temperature can now be manually adjusted to the desired setting.**

### To change from monitor mode to control mode

To change from monitor mode to controller mode, press either the '+' or '-' key on the SCT1 temperature controller once. The unit will revert to controller mod, and will be indicated by the illuminated amber contact thermometer indicator visible on the front panel of the hotplate. The red LED on the bottom right hand corner of the SCT1 temperature controller display will extinguish when the unit is in controller mode.

**NB: For safety reasons, when reverting from monitor to controller mode the SCT1 temperature controller will default to the set point temperature of 20°C.**

## Troubleshooting

In the event of malfunction the SCT1 temperature controller can self-diagnose certain fault conditions. In this case the display will show a message in the form of 'Err XX' where XX is either a single or double digit code, which identifies the fault.

---

### Error

code	Fault condition	Remedy
Err 1	Thermometer not calibrated	Return to factory
Err 2	Electronic fault	Return to factory
Err 4	Probe plug removed in control or monitor mode	Reconnect probe. Continuing error code indicates hotplate fault - return to factory
Err 8	Probe out of sample when in control mode	Replace probe in sample
Err 16	Ambient temperature inside hotplate too high (controller attached)	Switch off at 'mains' and allow to cool

---

## Servicing and Repair

**WARNING: Ensure the hotplate is disconnected from the mains electricity supply before attempting maintenance or servicing.**

This equipment does not require routine servicing. The only maintenance required is to clean external surfaces with a damp cloth and mild detergent. Spillages or splashes should be cleaned up straightaway after isolating the hotplate from the mains electricity supply.

### Repairs

In the event of product failure it is recommended that any repair is only undertaken by suitably qualified personnel.

**Only spare parts supplied by Cole-Parmer or its agent should be used. Fitting of non-approved parts may affect the performance of the safety features of the instrument.**

For a comprehensive list of parts required by service engineers conducting internal repairs, or to receive a service manual, please contact the Service Department of Cole-Parmer quoting both the model and serial number.

Cole-Parmer Ltd.  
Beacon Road,  
Stone, Staffordshire,  
ST15 0SA, United Kingdom  
Tel: +44 (0)1785 810475  
Email: cpservice@coleparmer.com  
Web: www.stuart-equipment.com

For any other technical enquiries, please contact the Technical Support Department:  
E-mail: cptechsupport@coleparmer.com  
Tel: +44 (0)1785 810433

## Accessories

A range of accessories detailed below is available for use with the SCT1 temperature controller.

Part code	Description
SCT1/1	Probe holder
SCT1/2	PTFE probe
SR1	Retort rod
SCT1/MP	Spare metal probe

## Warranty

Cole-Parmer Ltd. warrants this equipment to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of three (3) years. In the event of a justified claim, Cole-Parmer will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by installation, adaptation, modification, fitting of non-approved parts or repair by unauthorised personnel.

Cole-Parmer Ltd.  
Beacon Road,  
Stone, Staffordshire,  
ST15 0SA, United Kingdom  
Tel: +44 (0)1785 810475  
Email: cpservice@coleparmer.com  
Web: www.stuart-equipment.com

## Technical Specification

Probe	Stainless steel
Temperature range °C	20 to 200°C
Accuracy, °C	±0.5°C
Resolution	1°C
Dimensions, mm (w x d x h)	90 x 75 x 123
Net weight, kg	0.3 (inc. probe)
IP rating	54

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---

**CE** This product meets the applicable EC harmonized standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards

and so we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, contact the manufacturer.

## EU Declaration of Conformity

<b>Product</b>	Laboratory Equipment	<b>File Number</b>	P225
<b>Manufacturer</b>	Cole-Parmer Ltd Beacon Road Stone, Staffordshire ST15 0SA United Kingdom		

**This declaration of conformity is issued under the sole responsibility of the manufacturer**

**Object of Declaration** Melting Point Apparatus  
(reference the attached list of catalogue numbers)

The object of the declaration described above is in conformity with the relevant Union Harmonisation Legislation:

<b>Low Voltage Directive</b>	2006/95/EC	(until 19 April 2016)
	2014/35/EU	(from 20 April 2016)
<b>EMC Directive</b>	2004/108/EC	
<b>RoHS Directive</b>	2011/65/EC	

References to the relevant harmonised standards used or references to the other technical specifications in relation to which conformity is declared:

IEC/EN 61010-1:2010	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1: General requirements.
IEC/EN 61010-2-010:2014	Particular requirements for laboratory equipment for the heating of materials.
IEC/EN 61326-1:2006	Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1: General requirements (Class A).

**Signed for and on behalf of the above manufacturer**

**Additional Information** Year of CE Marking: 2008

**Place of Issue** Stone, Staffordshire, UK

**Date of Issue** June 2009. Revised 26 January 2017

**Authorised Representative** Carl Warren

**Title** Technical Manager

**Signature**

# INSPECTION REPORT

**MODEL** SCT1


## ELECTRICAL SAFETY

- 1. Earth continuity
- 2. Insulation
- 3. Flash test

## FUNCTIONAL

- 1. Indicators
- 2. Stirring function
- 3. Visual acceptance

**QUALITY CONTROL INSPECTOR**

**Cole-Parmer Ltd.**

Beacon Road,  
Stone,

Staffordshire,  
ST15 0SA,

United Kingdom

Tel: +44 (0)1785 812121

Email: [cpinfo@coleparmer.com](mailto:cpinfo@coleparmer.com)

Web: [www.stuart-equipment.com](http://www.stuart-equipment.com)