







INSTRUCTION FOR USE

Hygiene Sponge In Container

Samples should be kept cool until testing is performed.

<p>STEP 1</p> <p>Identify the site to be sampled. It is recommended to use a template for accurate results where possible (TS/15-T).</p> 	<p>STEP 2</p> <p>Unscrew the container cap and carefully remove the glove attempting not to touch the sponge or outer surface of glove fingers with bare skin.</p> 
<p>STEP 3</p> <p>Whilst wearing the glove, grip the sponge in the container and remove.</p> 	<p>STEP 4</p> <p>Use the sponge to collect the sample over the desired (or measured) area.</p> 
<p>STEP 5</p> <p>Place the sponge back in the container and screw the cap back on tight.</p> 	<p>STEP 6</p> <p>Store at 4 – 8°C (or cool & dry as possible) and send to the laboratory within 24 hours for analysis.</p> 

Have you considered? Sterile Sampling Sponge - TS/15-B



Technical Service
Consultants Ltd

www.tscswabs.co.uk
Tel: +44 (0)1706 620600

IFU_Sponge_Container_V4 – 23/05/2014