DOC027.98.00800

**POCKET COLORIMETER™ II -
Phosphate LCK 348, LCK 349, LCK 350
Phosphat LCK 348, LCK 349, LCK 350**

05/2011 Edition 3A

APPLICATION INSTRUCTION
PROGRAMMIERANLEITUNG
INSTRUCTIONS DE FONCTIONNEMENT
INSTRUCCIONES DE APLICACIÓN
ISTRUZIONI APPLICATIVE
GEBRUIKSAANWIJZING
INSTRUKTIONSHANDBOK

Table of contents

Section 1 General information	17
1.1 Safety information	17
1.2 General product information	17
Section 2 Installation	19
2.1 Packing list	19
2.2 Insert the adapter	19
2.3 Replace the instrument cap	19
Section 3 System start up	21
3.1 Instrument keys and display	21
3.2 Enter calibration data of LANGE Cuvette Tests	21
3.3 Programming procedure	23
3.4 Check entered data	26
3.5 Overwrite stored calibration data	27
Section 4 Operation	29
4.1 Analyse Phosphate-Phosphorus/Phosphate with LANGE Cuvette Tests ..	29
4.2 Quality assurance	29
4.3 Evaluate LANGE Phosphate-Phosphorus/ Phosphate Cuvette Tests	29
Section 5 Replacement Parts and Accessories	31
5.1 Replacement Parts	31
5.2 Accessories	31
Section 6 Warranty and liability	33

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Section 1 General information

1.1 Safety information

Please read this entire application instruction and the POCKET COLORIMETER™ II instruction manual before unpacking, setting up, or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

To ensure that the protection provided by this equipment is not impaired, do not use or install this equipment in any manner other than that specified in this manual.

1.1.1 Use of hazard information

DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

Important Note: Information that requires special emphasis.

Note: Information that supplements points in the main text.

1.2 General product information

This application instruction describes how to evaluate LANGE Cuvette Tests on a single-wavelength POCKET COLORIMETER II.

Section 2 Installation

Important Note: This application instruction is only valid in conjunction with the POCKET COLORIMETER™ II instruction manual and the applied LCK working procedure. Please read this entire application instruction and the entire POCKET COLORIMETER II instruction manual before unpacking, setting up or operating the instrument.

2.1 Packing list

- POCKET COLORIMETER II 655 nm, Cat. No. 5870065
- Instrument cap with HACH LANGE label, Cat. No. LZV 759
- Cuvette adapter 1 cm, Cat. No. 5954600
- Blank value cuvettes, Cat. No. LCW 919
- Application instruction DOC027.98.00800

2.2 Insert the adapter

Remove the instrument cap and carefully position the cuvette adapter into the cell holder with the arrow pointing downward. Fix it by pressing it down.

2.3 Replace the instrument cap

Use the instrument cap with the HACH LANGE label when evaluating LANGE cuvettes.

Section 3 System start up

3.1 Instrument keys and display


Figure 1 Instrument keys and display

1	POWER/BACKLIGHT	6	Measuring range indicator (2)
2	ZERO/SCROLL (BLUE key)	7	Menu indicator
3	MENU	8	Calibration adjusted indicator
4	Numeric display	9	Battery low indicator
5	Measuring range indicator (1)	10	READ/ENTER (GREEN key)

3.1.1 Key-Functions

ZERO/SCROLL (BLUE key) : In measurement mode, sets the instrument to zero. In menu mode, scrolls through menu options. Also scrolls numbers when entering or editing a value.

READ/ENTER (GREEN key) : In measurement mode, initializes sample measurement. In menu mode, selects a menu option. When entering numbers, moves one space to the right. Confirm the function when the entry is complete.

3.2 Enter calibration data of LANGE Cuvette Tests

3.2.1 Overview

Each POCKET COLORIMETER II has two channels in which measurements can be made. Each channel (measuring range) will accept user-defined calibration data.

This application instruction contains the calibration data of the LANGE Cuvette Tests for Phosphate-Phosphorus/Phosphate:

-
- LCK 348: 0.5-5 mg/l PO₄-P or 1.5-15 mg/l PO₄
 - LCK 349: 0.05-1.5 mg/l PO₄-P or 0.15-4.5 mg/l PO₄
 - LCK 350: 2-20 mg/l PO₄-P or 6-60 mg/l PO₄

Select the calibration data of maximum two measuring ranges which will be the suitable one.

Each programming procedure has only to be followed once prior to the first test evaluation.

For detailed information please read the *POCKET COLORIMETER II instruction manual*.

3.2.2 Data table

3.2.2.1 Evaluation as Phosphate-Phosphorus PO₄-P

Table 1 Evaluation as Phosphate-Phosphorus PO₄-P

	LCK 348	LCK 349	LCK 350
RES (resolution)	--.--	--.---	--.--
S0 (standard 0)	00.50	00.05	02.00
A0 (absorbance 0)	0.104	0.090	0.104
S1 (standard 1)	05.00	01.50	20.00
A1 (absorbance 1)	0.600	0.728	0.753

3.2.2.2 Evaluation as Phosphate PO₄

Table 2 Evaluation as Phosphate PO₄

	LCK 348	LCK 349	LCK 350
RES (resolution)	--.---	--.---	--.---
S0 (standard 0)	01.50	00.15	06.00
A0 (absorbance 0)	0.104	0.090	0.104
S1 (standard 1)	15.00	04.50	60.00
A1 (absorbance 1)	0.600	0.728	0.753

3.3 Programming procedure

Note: Each programming procedure has only to be followed once prior to the first test evaluation.

Important Note: The following programming procedures are valid for the evaluation form Phosphate-Phosphorus PO₄-P. To evaluate as Phosphate PO₄ please choose the data mentioned in table 3.2.2.2 Evaluation as Phosphate PO₄ on page 23.

Note: After programming the calibration data, it could be helpful to mark the channel(s) 1 and/or 2 accordingly with a label.

3.3.1 LCK 348 Phosphate-Phosphorus PO₄-P

1. Press **POWER** to switch on the instrument.
2. Press **MENU**. "SEL" is displayed. A flashing arrow indicates the current measuring range (1 or 2).
3. Select measuring range 1 or 2 for LCK 348 (0.5-5 mg/l PO₄-P).
Note: Press **GREEN** key to toggle between measuring ranges.
4. Press **MENU** for about five seconds until "USER" followed by "CAL" is displayed.
5. Press **BLUE** key until "Edit" is displayed. Press **GREEN** key. "RES" is displayed.
6. Press **GREEN** key. "-.---" is displayed.
7. Press **BLUE** key until "--.---" is displayed. Press **GREEN** key. "S0" is displayed.
8. Press **GREEN** key until "A0" is displayed.
9. Press **GREEN** key until "S1" is displayed.
10. Press **GREEN** key. Enter "05.00" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "A1" is displayed.
11. Press **GREEN** key. Enter "0.600" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "Add" is displayed.

12. Press **MENU** until "Edit" is displayed. Press **GREEN** key until „S0“ is displayed.
13. Press **GREEN** key. Enter "00.50" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "A0" is displayed.
14. Press **GREEN** key. Enter "0.104" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until „S1“ is displayed.
15. Press **MENU** until "Edit" is displayed to check the entered data. Press **GREEN** key. „RES“ is displayed.
16. Press **BLUE** key to scroll through the entries. Compare the entries with the data shown in the table below.

Table 3 Calibration data of LCK 348

RES	--.--
S0	0.50
A0	0.104
S1	5.00
A1	0.600

17. After comparing the entries press **MENU** two times to enter the measuring mode. See [4.3 Evaluate LANGE Phosphate-Phosphorus/ Phosphate Cuvette Tests on page 29](#).

Note: In case of any mistake retrieve the factory calibration first. See [3.5.1 Retrieve the factory calibration on page 27](#). After retrieving the factory calibration, start the programming procedure again. If the instrument is shut-off or loses power during data entry, all edits will be lost. Automatic shut-off is 60 minutes.

3.3.2 LCK 349 Phosphate-Phosphorus PO₄-P

1. Press **POWER** to switch on the instrument.
2. Press **MENU**. "SEL" is displayed. A flashing arrow indicates the current measuring range (1 or 2).
3. Select measuring range 1 or 2 for LCK 349 (0.05-1.5 mg/l PO₄-P).
Note: Press **GREEN** key to toggle between measuring ranges.
4. Press **MENU** for about five seconds until "USER" followed by "CAL" is displayed.
5. Press **BLUE** key until "Edit" is displayed. Press **GREEN** key. " RES" is displayed.
6. Press **GREEN** key. "-.---" is displayed.
7. Press **BLUE** key until "--.--" is displayed. Press **GREEN** key. "S0" is displayed.
8. Press **GREEN** key until "A0" is displayed.
9. Press **GREEN** key until "S1" is displayed.

10. Press **GREEN** key. Enter "01.50" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "A1" is displayed.
11. Press **GREEN** key. Enter "0.728" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "Add" is displayed.
12. Press **MENU** until "Edit" is displayed. Press **GREEN** key until „S0“ is displayed.
13. Press **GREEN** key. Enter "00.05" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "A0" is displayed.
14. Press **GREEN** key. Enter "0.090" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until „S1“ is displayed.
15. Press **MENU** until "Edit" is displayed to check the entered data. Press **GREEN** key. „RES“ is displayed.
16. Press **BLUE** key to scroll through the entries. Compare the entries with the data shown in the table below.

Table 4 Calibration data of LCK 349

RES	--.--
S0	0.05
A0	0.090
S1	1.50
A1	0.728

17. After comparing the entries press **MENU** two times to enter the measuring mode. See [4.3 Evaluate LANGE Phosphate-Phosphorus/ Phosphate Cuvette Tests on page 29](#).

Note: In case of any mistake retrieve the factory calibration first. See [3.5.1 Retrieve the factory calibration on page 27](#). After retrieving the factory calibration, start the programming procedure again. If the instrument is shut-off or loses power during data entry, all edits will be lost. Automatic shut-off is 60 minutes.

3.3.3 LCK 350 Phosphate-Phosphorus PO₄-P

1. Press **POWER** to switch on the instrument.
2. Press **MENU**. "SEL" is displayed. A flashing arrow indicates the current measuring range (1 or 2).
3. Select measuring range 1 or 2 for LCK 350 (2-20 mg/l PO₄-P).
Note: Press **GREEN** key to toggle between measuring ranges.
4. Press **MENU** for about five seconds until "USER" followed by "CAL" is displayed.
5. Press **BLUE** key until "Edit" is displayed. Press **GREEN** key. " RES" is displayed.
6. Press **GREEN** key. "-.---" is displayed.

7. Press **BLUE** key until "--.--" is displayed. Press **GREEN** key until "S0" is displayed.
8. Press **GREEN** key until "A0" is displayed.
9. Press **GREEN** key until "S1" is displayed.
10. Press **GREEN** key. Enter "20.00" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "A1" is displayed.
11. Press **GREEN** key. Enter "0.753" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "Add" is displayed.
12. Press **MENU** until "Edit" is displayed. Press **GREEN** key until „S0“ is displayed.
13. Press **GREEN** key. Enter "02.00" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until "A0" is displayed.
14. Press **GREEN** key. Enter "0.104" by using **BLUE** key to change value and **GREEN** key to move to the next space. Press **GREEN** key until „S1“ is displayed.
15. Press **MENU** until "Edit" is displayed to check the entered data. Press **GREEN** key. „RES“ is displayed.
16. Press **BLUE** key to scroll through the entries. Compare the entries with the data shown in the table below.

Table 5 Calibration data of LCK 350

RES	--.--
S0	2.00
A0	0.104
S1	20.00
A1	0.753

17. After comparing the entries press **MENU** two times to enter the measuring mode. See [4.3 Evaluate LANGE Phosphate-Phosphorus/ Phosphate Cuvette Tests on page 29](#).

Note: In case of any mistake retrieve the factory calibration first. See [3.5.1 Retrieve the factory calibration on page 27](#). After retrieving the factory calibration, start the programming procedure again. If the instrument is shut-off or loses power during data entry, all edits will be lost. Automatic shut-off is 60 minutes.

3.4 Check entered data

To avoid erroneous measurements it is recommended to check the data after completing the programme procedure.

1. Press **MENU** for about five seconds until "USER" followed by "CAL" is displayed.
2. Press **BLUE** key to scroll to "Edit". Press **GREEN** key.

-
3. "RES" is displayed. Press **BLUE** key to scroll through the entries. Compare the entries with the figures in the [3.2.2.1 Evaluation as Phosphate-Phosphorus PO4-P on page 22](#).

In case a value was not stored correctly, retrieve the factory calibration first. See [3.5.1 Retrieve the factory calibration on page 27](#). After retrieving the factory calibration, start the programming procedure again.

3.5 Overwrite stored calibration data

With the POCKET COLORIMETER II it is possible to overwrite stored calibration data. Before overwriting calibration data it is recommended to retrieve the factory calibration first. After retrieving the factory calibration, you can program new calibration data.

In case of any mistake retrieve the factory calibration first.

3.5.1 Retrieve the factory calibration

1. Press **MENU**. "SEL" is displayed. A flashing arrow indicates the current measuring range (1 or 2). Press **GREEN** key to toggle between measuring ranges. Select measuring range retrieving the factory calibration.
1. Press **MENU** for about five seconds until "USER" followed by "CAL" is displayed.
2. Press **BLUE** key to scroll to "dFL". Press **GREEN** key.

Section 4 Operation

4.1 Analyse Phosphate-Phosphorus/Phosphate with LANGE Cuvette Tests

DANGER

Handling chemical samples, standards and reagents can be dangerous. Users of this product are advised to familiarize themselves with safety procedures and the correct use of chemicals, and to carefully read all relevant material safety data sheets.

To analyse Phosphate-Phosphorus/Phosphate (ortho or total) follow the working procedure enclosed in the according Cuvette Test LCK 348, LCK 349 or LCK 350.

4.2 Quality assurance

The measurement results must be subjected to plausibility checks (dilute and/or spike the sample).

Regular checks with corresponding standard solutions (see [Section 5 Replacement Parts and Accessories on page 31](#)) ensure that your measurement system is functioning properly and is being correctly operated.

4.3 Evaluate LANGE Phosphate-Phosphorus/ Phosphate Cuvette Tests

Note: Make sure that the instrument is programmed with the corresponding calibration data, see [3.4 Check entered data on page 26](#).

1. Fill DI water into an empty cuvette (LCW 919) up to the lower edge of the barcode label.
2. Close the cuvette with the rubber stopper. This is the blank cuvette.
3. Thoroughly clean the outside of the cuvette.
4. Press **POWER** to switch on the instrument.
5. Make sure that the correct measuring range is selected! Otherwise press **MENU**. "SEL" is displayed. A flashing arrow indicates the current measuring range (1 or 2). Press **GREEN** key to select measuring range. Press **MENU** to get back to the measuring mode.

6. Insert the blank cuvette into the cell holder equipped with the cuvette adapter 1 cm.

Note: Make sure that the cuvette adapter is positioned correctly.

7. Place the instrument cap with the HACH LANGE label over the cell holder to cover the cuvette.
8. Press **BLUE** key. "0.00" is displayed.
9. Remove the instrument cap and insert the sample cuvette.
10. Place the instrument cap over the cell holder to cover the cuvette.
11. Press **GREEN** key.

12. The result is displayed in mg/l Phosphate-Phosphorus/Phosphate.

Note: *In case the Phosphate-Phosphorus/Phosphate concentration is out of the selected measuring range the value for the lowest or highest concentration is displayed and flashes. Repeat the test using the appropriate LCK Cuvette Test measuring range or dilute your sample.*

Note: *A new blank measurement should be performed prior to each individual test series.*

Section 5 Replacement Parts and Accessories

5.1 Replacement Parts

Description	Quantity	Cat. Number
POCKET COLORIMETER™ II 655 nm	1	5870065
Instrument cap with HACH LANGE label	1	LZV 759
Cuvette adapter 1 cm	1	5954600
Blank value cuvettes	1 pk/5 pc	LCW 919

5.2 Accessories

Description	Quantity	Cat. Number
Phosphate Cuvette Test (0.5-5 mg/l PO ₄ -P)	1 pk/25 pc	LCK 348
Phosphate Cuvette Test (0.05-1.5 mg/l PO ₄ -P)	1 pk/25 pc	LCK 349
Phosphate Cuvette Test (2-20 mg/l PO ₄ -P)	1 pk/25 pc	LCK 350
Standard and Ring Test solution for LCK 348 (ortho)	1	LCA 700
Standard and Ring Test solution for LCK 348 (total)	1	LCA 707
Standard and Ring Test solution for LCK 349 (ortho)	1	LCA 704
Standard and Ring Test solution for LCK 349 (total)	1	LCA 709
Standard and Ring Test solution for LCK 350 (ortho)	1	LCA 703
Standard and Ring Test solution for LCK 350 (total)	1	LCA 708
Pipette 0.2-1 ml	1	BBP 078
Pipette tips for BBP 078	1 pk/100 pc	BBP 079
Pipette 1-5 ml	1	BBP 065
Pipette tips for BBP 065	1 pk/75 pc	BBP 068
Cuvette Stand for LANGE cuvettes	1	LYW 915
Thermostat LT200	1	LTV082.99.10002

Section 6 Warranty and liability

The manufacturer warrants that the product supplied is free of material and manufacturing defects and undertakes the obligation to repair or replace any defective parts at zero cost.

The warranty period for instruments is 24 months. If a service contract is taken out within 6 months of purchase, the warranty period is extended to 60 months.

With the exclusion of the further claims, the supplier is liable for defects including the lack of assured properties as follows: all those parts that, within the warranty period calculated from the day of the transfer of risk, can be demonstrated to have become unusable or that can only be used with significant limitations due to a situation present prior to the transfer of risk, in particular due to incorrect design, poor materials or inadequate finish will be improved or replaced, at the supplier's discretion. The identification of such defects must be notified to the supplier in writing without delay, however at the latest 7 days after the identification of the fault. If the customer fails to notify the supplier, the product is considered approved despite the defect. Further liability for any direct or indirect damages is not accepted.

If instrument-specific maintenance and servicing work defined by the supplier is to be performed within the warranty period by the customer (maintenance) or by the supplier (servicing) and these requirements are not met, claims for damages due to the failure to comply with the requirements are rendered void.

Any further claims, in particular claims for consequential damages cannot be made.

Consumables and damage caused by improper handling, poor installation or incorrect use are excluded from this clause.

The manufacturer process instruments are of proven reliability in many applications and are therefore often used in automatic control loops to provide the most economical possible operation of the related process.

To avoid or limit consequential damage, it is therefore recommended to design the control loop such that a malfunction in an instrument results in an automatic change over to the backup control system; this is the safest operating state for the environment and the process.

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