
Cole-Parmer®

FJT-200 Series

6 Position Flocculator

Instruction Manual

SW1061(ES) Version 1.7

Cole-Parmer®
essentials

English


Figure 1: FJT-200D-6


Figure 2: Control panel

Thank you for purchasing this Cole-Parmer product. To get the best performance from the equipment, and for your personal safety, please read these instructions carefully before use.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuations not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

Electrical Installation

THIS EQUIPMENT MUST BE EARTHED

Before connection please ensure that the line supply corresponds to that stated on the rating label.

Power requirements:

Model	Voltage & Wattage
FJT-200D-6	230V, 50Hz, ~ , 200W

There is an IEC socket at the side of the instrument for connection to the mains supply. The unit is supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Schuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. See the enclosed instruction sheet for advice on how to carry out this procedure.

Should the mains lead require replacement a cable of 1mm² of harmonised code H05W-F connected to an IEC 320 plug should be used.

IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

The mains lead should be connected to the instrument BEFORE connection to the mains supply.

Safety Advice Before Use

- ❖ Never move or carry the unit when in use or connected to the mains electricity supply.
- ❖ Do not mix flammable liquids or use the equipment in hazardous atmospheres.
- ❖ In the case of mains interruption, a fault or mechanical failure, the unit will NOT continue to operate on removal of fault.
- ❖ Mechanical energy can lead to breakage of glass vessels. Use with care.

General Description

The Cole-Parmer FJT-200D-6 Flocculator is designed to ensure repeatable conditions between up to six samples during water and sewage flocculation tests.

The speed control has electronic feedback to accurately maintain the rotational speed of the stainless steel paddles. Digital displays clearly indicate the speed, which is variable between 25 and 250 rpm and time, which can be set to count down from 1 to 99 minutes. Alternatively, the use of two adjustable pre-set programs (P1 and P2) allows quick and accurate repetition of pre-selected speeds and times.

All paddles can be easily lifted out of the samples when required. Observation and monitoring of samples is assisted by diffused base illumination and a black background.

A coagulant dispenser is available as an accessory.

Operation (see Fig 1 & 2)

Place the unit on a firm, level surface and connect to the mains supply. Leave sufficient space to access the side IEC socket.

Turn the unit ON at the side mains On / Off switch and observe the fluorescent lamp illuminate.

To position a beaker, lift the paddle up to its maximum height using the paddle lift and place the vessel over the illuminated base. Once all the samples are in position, lower the paddles into the beakers to the required depth.

Selecting a speed / time

Press the <MAN> button to selected the manual mode. To set the speed, hold down the <SET> key under the RPM display and use the <UP> and <DOWN> arrow keys until the desired speed appears on the digital display.

To set a running time, hold down the <SET> key under the TIME display and use the <UP> and <DOWN> arrow keys until the desired running time in minutes appears on the digital display.

IF A RUNNING TIME IS NOT REQUIRED, SET 00 IN THE TIME DISPLAY.

To start the run press the <START / STOP> key. The paddles will turn at the selected speed. If a time has been set, the timer will count down.

When time zero is reached, the unit will beep and the paddles stop rotating. The timer display will then flash and count up, in minutes, to indicate time elapsed since the paddles stopped.

The program can be halted at any time by pressing the <START / STOP> key. To change the settings, stop rotation and repeat the above procedure.

Pre-set programs

In order to input a pre-set speed and time, press and hold the <P1> key for 5 seconds. The green P1 light will flash. To set the speed, hold down the <SET> key under the RPM display and use the <UP> and <DOWN> arrow keys until the desired speed appears on the digital display. To set a running time, hold down the <SET> key under the TIME display and use the <UP> and <DOWN> arrow keys until the desired time in minutes appears on the digital display.

IF A RUNNING TIME IS NOT REQUIRED, SET 00 IN THE TIME DISPLAY.

Once the correct speed and time have been set, press the <P1> key to confirm. The settings will be held in the memory.

Pressing the <P1> key at any time during operation will select the preset values. A green LED shows that the P1 program has been selected.

If another preset speed and time is required press and hold down the <P2> key and repeat the above setting procedure.

New P1 and P2 values can be inputted at any time by repeating the above procedure.

Maintenance & Servicing

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

The following spares and accessories are available from your laboratory dealer.

Description	Catalogue Number
Coagulant Dispenser	SW6/2

For a comprehensive list of parts required by service engineers conducting internal repairs and a service manual, please contact the Service Department at Antylia Scientific quoting both the model and serial number.

Only spare parts supplied or specified by Antylia Scientific or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

If in any doubt, please contact the Service Department of Antylia Scientific or the point of sale.

Technical Specification

FJT-200D-6

No. of samples	6
Speed range	25-250rpm
Timer	0 to 99 minutes
Pre-set programs	2
Digital displays	LED
Dimensions	750 x 210 x 460mm
Net weight	17kg
Power	200W


Figure 1: FJT-200D-6


Figure 2: Tableau de commande

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Notes

CE UK CA This product meets the applicable CE Directives and UKCA Legislation for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in its vicinity will meet these standards and so we cannot guarantee

that interference will not occur in practise. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advise before use, contact the manufacturer.

Declaration of Conformity is available to view online at www.coleparmer.com

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Ordering Information

Order No.	Series	Model	Legacy SKU
99551-05	FJT-200	FJT-200D-6	SW6

Warranty Registration

UK

T: +44 (0) 1480 272279
E: uk.sales@antylia.com
W: coleparmer.co.uk

India

T: +9122 61394444
E: info@coleparmer.in
W: coleparmer.in

Germany

T: +49 (0) 9377 92030
E: de.sales@antylia.com
W: coleparmer.de

China

T: +1 847 549 7600
E: sales@antylia.com
W: coleparmer.com

France

T: +33 (0) 1486 37800
E: fr.sales@antylia.com
W: coleparmer.fr

USA

T: +1 847 549 7600
E: sales@antylia.com
W: coleparmer.com

Italy

T: +39 (0) 284349215
E: it.sales@antylia.com
W: coleparmer.it

Canada

T: +514 355 6100
E: info@antylia.ca
W: coleparmer.ca

Other

T: +1 847 549 7600

Cole-Parmer®
essentials

Antylia Scientific Ltd.
Beacon Road,
Stone,
Staffordshire,
ST15 0SA,
United Kingdom