

Validation Package

BRAIDMOTION

Braided Silicone Pharma Tubing
made by RAUMEDIC

| | |
|---|-----------|
| 1. Introduction | 04 |
| 1.1 About RAUMEDIC | 04 |
| 1.2 Security of Supply | 05 |
| 1.3 cGMP Conformity | 05 |
| 1.4 Quality Management System | 05 |
| 2. Technical Specifications | 06 |
| 2.1 Tubing dimensions, Type and Part Numbers Overview | 06 |
| 2.2 Material Information | 06 |
| 2.3 Material Compliance | 07 |
| 2.4 Physical Properties | 08 |
| 2.5 Physico-Chemical Properties | 08 |
| 2.6 Biocompatibility | 09 |
| 2.7 Tubing Printing | 09 |
| 2.8 Extractable Profile | 09 |
| 2.9 Sterilization Compatibility | 10 |
| 2.10 General Properties, Shelf Life and Storage | 10 |
| 3. Pressure Testing | 11 |
| A. Additional Documents | 12 |
| A.1 Certificates – ISO 13485:2016 | 12 |
| A.2 Certificates – ISO 15378:2018 | 14 |
| A.3 Certificates – ISO 50001:2018 | 15 |
| A.4 Documentation – Material Safety Data Sheet | 17 |



For us the issue of the customer's approval is the proof of the functional efficiency respectively of the fitness of our articles, which the customer determined for the specification application within the scope of testing. The obligated quality of goods has to conform to the expressly stipulated performance features in written form only (e.g. Specifications, Technical Delivery Specification, drawings, markings and other data).

We only assume a warranty exceeding this stipulation of obligated quality for a specific application or a specific fitness, period of application or life

period after passing of risk, if this is expressly agreed upon in written form; otherwise the risk of application and fitness is solely with the purchaser.

RAUMEDIC products and medical devices manufactured from these may not be sold or used in applications in the human body lasting more than 24 hours. No other guarantees are or will be given by RAUMEDIC.

This document is valid until revoked or until new information is provided. If you require further information, please contact your responsible sales contact.

1. Introduction

1.1 About RAUMEDIC

As a partner of the international medical technology and pharmaceutical industry, RAUMEDIC develops and produces components for customers, including tubing, catheters, and molded parts as well as complex groups of components and systems for a wide range of diagnostic and therapeutic uses. For the clinical areas of neuro-monitoring and traumatology, RAUMEDIC produces high-precision pressure-measurement sensors. The company is 100 percent family-owned.



The information in this document is subject to change without notice and should not be construed as a commitment by RAUMEDIC.

1.2 Security of Supply

The RAUMEDIC manufacturing plants in Europe and the US operate under ISO 7 clean room conditions in order to offer flexibility and reliability to our customers worldwide. Consistent process performance is ensured by the ongoing qualification of all manufacturing processes and personnel.

1.3 cGMP Conformity

Consistent high quality of silicone tubing is assured by careful selection of the raw material, well planned and validated production technologies and an efficient Quality Assurance Department, all of which result in high batch-to-batch reproducibility.

1.4 Quality Management System

RAUMEDIC is certified according to ISO 13485. The respective certificate can be found attached to this document in appendix A.1 Certificates.

2.0 Technical Specifications

2.1 Tubing Dimensions, Type and Part Numbers Overview

BRAIDMOTION high-pressure Silicone tubing coils made of Silicone SIK8649 and polyester yarn are of dimensions between 1/8" x 0,355" (ID x OD) and 1" x 1-3/8" (ID x OD).

The color of the tubing is natural without addition of color pigments (RAUMEDIC color no. 20900). All tubing coils are provided in double PE bags packaged in a cardboard box.

| Dimension ID x OD | Tolerances | | Coil length [m] | Printed / non-printed | Article code |
|------------------------------------|------------------------|------------------------|--------------------|--------------------------|--------------|
| | Inner diameter [mm] | Wall thickness [mm] | | | |
| 1/8" x 0,355" 3,2 mm x 9,0 mm | 2,70 - 3,2 | 2,70 - 3,2 | 25 | printed | 807710-800 |
| 1/4" x 1/2" 6,4 mm x 12,7 mm | 5,85 - 6,45 | 2,90 - 3,50 | 15 | printed | 807711-800 |
| 3/8" x 5/8" 9,5 mm x 15,9 mm | 9,05 - 9,75 | 2,90 - 3,50 | 15 | printed | 807712-800 |
| 1/2" x 7/8" 12,7 mm x 22,2 mm | 12,20 - 13,00 | 4,30 - 5,10 | 15 | printed | 807713-800 |
| 3/4" x 1-1/8" 19,1 mm x 28,6 mm | 18,30 - 19,50 | 4,30 - 5,30 | 10 | printed | 807714-800 |
| 1" x 1-3/8" 25,4 mm x 34,9 mm | 24,70 - 26,10 | 4,30 - 5,30 | 7,5 | printed | 807715-800 |

2.2 Material Information

BRAIDMOTION silicone pharma tubing is made of the material grade Silicone SIK8649, a high purity silicone rubber, and a medical-grade multi-filament polyester yarn, suitable for pharma applications. Chemically Silicone SIK8649 is an addition cross-linked hot vulcanisate based on vinyl methyl dimethyl polysiloxane using silicic acid fillers and a platinum catalyst.

BRAIDMOTION silicone tubing can be used continuously in a temperature range from -20°C to +135°C without losing its integrity or deterioration of chemical/ physiological properties. A short-term use up to 200°C is possible after validation by the customer.

The material Silicone SIK8649 shows good resistance to water up to 100 °C and to low-pressure steam up to 135 °C. Steam from higher temperatures leads to degradation, especially after prolonged exposure.

Silicone SIK8649 has good resistance to weak acids and alkalis. However, the material will degrade by strong acids and alkalis, which will be promoted at higher temperatures.

The resistance is strongly dependent on the polar or non-polar character of the contact medium. While there is almost no swelling in polar contact media (e.g. water/alcohol), non-polar contact media (e.g. petrol/oil) cause medium to strong but reversible swelling.

The surface is coated in a plasma process. This coating provides a less sticky surface of these braided silicone tubing (Low Tack - see also 2.7) in comparison with common non-coated silicone products.

Braided silicone tubing is colourless transparent or translucent.

2.3 Material Compliance

Based on our experience and knowledge about the materials and production processes the base material formulation complies with / does not contain:

- Bisphenol A (BPA)
- Latex
- Phthalates
- Endocrine disruptors
- Ozone depleting substances (Regulation (EC) No 1005/2009)
- Perfluorinated chemicals (PFOA, PFOS)
- Persistent Bioaccumulative Toxins (PBTs)
- Persistent Organic Pollutants (POPs) (Regulation (EU) 2019/1021)
- Genetically Modified Organisms (GMO)
- Animal derived materials (Materials which present BSE/TSE risks)
- Conflict minerals (Dodd-Frank Act, Section 1502)
- RoHS (Directive 2011/65/EU, Directive (EU) 2015/863)
- Heavy Metals (CONEG, Directive 94/62/EC)
- CMR substances (Regulation (EC) No 1272/2008)
- Substances of Very High Concern (SVHC) as identified under article 33 of Regulation (EC) No 1907/2006 (REACH or substances listed in Annex XIV and Annex XVII)
- Nanomaterials (acc. to the definition Medical Device Regulation (EU) 2017/745) as intentionally added components above threshold limits stated in aforementioned Regulations, Directives or Laws

as intentionally added components above threshold limits stated in aforementioned Regulations, Directives or Laws.

However, the raw materials used, and our final product(s) are not analyzed for the presence of traces of the above-mentioned substances.

If a provision of this Certificate should be or become invalid or loophole, the validity of the other provisions shall not be affected thereby.

If you require further information, please contact your responsible sales contact.

2.4 Physical Properties

| Property | Norm | SIK8649 |
|-----------------------------------|------------------|---------|
| Shore hardness A | ISO 48-4 | 60 ± 5 |
| Density (g/cm ³) | ISO 1183 | 1.15 |
| Tear strength (MPa) | ISO 527 | ≥ 8.0 |
| Elongation at break (%) | ISO 527 | ≥ 500 |
| Compression Set (%; 22h / 175 °C) | ISO 815-1 type B | ≤ 15 % |

The physical properties of the raw material were determined on standard test specimen. The specified values are approximate values.

2.5 Physico-Chemical Properties

European Pharmacopoeia: 3.1.9 & FDA regulation 21CFR, § 177.2600 & § 177.2800.

Test results | BRAIDMOTION silicone tubing in silicone rubber SIK8649 and polyester yarn meets the requirements of the European Pharmacopoeia 3.1.9. and regulation 21 CFR, § 177.2600, § 177.2800. The test methods, limits and results are those described by the E.P. monograph and listed in the table below.

| Test Description | E.P. 3.1.9 Limits | Results on tubing: Pass or Fail |
|------------------------------|---------------------|------------------------------------|
| Appearance of Solution | Colourless | Pass |
| Acidity | ≤ 2.5 mL NaOH 0.01M | Pass |
| Alkalinity | ≤ 1.0 mL HCl 0.01M | Pass |
| Reducing Substances | ≤ 1 mL | Pass |
| Substances Soluble in Hexane | ≤ 15 mg | Pass |
| Volatile Matter | < 2 % | Pass |
| Mineral Oils | < 1 ppm | Pass |
| Platinum | < 30 ppm | Pass |

USP <661> – Containers, Physico-chemical Tests – Plastic

Test results | BRAIDMOTION silicone tubing meets the USP <661> requirements when sterilized at 50 kGy with and without aging conditions corresponding to a shelf life of 5 years.

| Test Description | USP <661> Limits | Results on tubing: Pass or Fail |
|----------------------|------------------|------------------------------------|
| Non Volatile Residue | < 15 mg | Pass |
| Residue on Ignition | < 5 mg | Pass |
| Heavy Metals | < 1 ppm | Pass |
| Buffering Capacity | < 10 mL | Pass |

2.6 Biocompatibility

Samples of BRAIDMOTION silicone tubing, consisting of SIK8649 (inner/ outer tube layer) and polyester yarn (middle layer), has passed the following biocompatibility tests:

| Testing | Standard | Result |
|---|----------|--------|
| Biological reactivity, in vivo (class VI) | USP <88> | Pass |
| Biological reactivity, in vitro | USP <87> | Pass |
| Bacterial endotoxins | USP <85> | Pass |

2.7 Tubing Printing

The tubing is available with a standardized printing.

The standardized printing text is defined as follows:

„BRAIDMOTION A“ x B“ MADE BY RAUMEDIC LOT XXXXXXXX“

A - Inner diameter in inch

B - Outer diameter in inch

BRAIDMOTION silicone tubing features the patented Low-Tack surface treatment providing enhanced technical and handling properties.

The biocompatibility of the ink was proven by successfully performed ISO 10993-5 and USP <88> testing.

2.8 Extractable Profile

The study on the extractables profile for RAUMEDIC non-braided silicone pharma tubing was performed according to BPOG “BioPhorum Best Practices Guide for Extractables Testing of Polymeric Single-Use Components Used in biopharmaceutical Manufacturing” released in April 2020. RAUMEDIC has no indication of a significant migration effect caused by the polyester braid.

The assessment on gamma-irradiated and autoclaved non-braided silicone pharma tubing will be provided upon request.

2.9 Sterilization Compatibility

It is the responsibility of the user to validate a sterilization process with autoclave for BRAIDMOTION silicone tubing.

BRAIDMOTION silicone tubing is produced for single use. The product will be delivered non-sterile.

Sterilization with gas (ethylene oxide, ETO), steam (up to 135 °C), gamma or X-rays (Dose: max. 50 kGy) is generally recommended. The suitability of the material must be evaluated by the user. Possible material damage depends on the sterilization conditions (temperature, pressure, energy, time, packaging unit and distance to the radiation source).

Multiple sterilizations in the scope of CIP/ SIP cycles are technically possible, but the corresponding process validation lies within the responsibility of the customer.

RAUMEDIC is not able to anticipate implications that may be caused by the intended

application conditions and repeated use (e.g., media contact, temperature, pressure, time, etc.).

However, sterilization of BRAIDMOTION silicone tubing with multiple CIP/SIP-cycles at 135 °C does not impair fit, form or function of the products. The amount of the repetitions and the process conditions used are within the validation responsibility of the user. The user is obliged to examine the resterilizability under application-related conditions, since the number of resterilization cycles may be influenced by the flow medium/ cleaning agent and/ or mechanical stress applied throughout application lifetime.

2.10 General Properties, Shelf Life and Storage

According to the current state-of-the-art individual fisheyes due to raw material and processing, foreign material, dirt inclusions and air bubbles as well as contamination on the tubing surface, like intrinsic particles and fluff, cannot be completely excluded.

The following shelf life is valid for the non-sterile, originally packed products.

A slight discoloration of the material might appear over storage time which may vary from batch to batch in intensity. This material related aging effect is common and does neither affect functionality nor material properties of Silicone SIK8649 in the scope of specified physical properties as well as shelf-life.

The shelf life of the tubing made of RAUMEDIC Silicone SIK8649 in non-sterile condition and originally packed is 5 years. There is no substantial change of the specified physical, chemical and physiological parameters of the tubing material after that period of time.

The following storage conditions shall apply:

- Storage in original packaging
- Dry (Standard conditions 10 – 75% relative humidity)
- Temperature 15 – 25°C with short-term deviations 5 – 35°C
- Goods protected from moisture
- Goods protected from directly exposure to UV-light
- Not in the environment and in presence of chemicals, e.g. solvent and disinfectant

3. Pressure Testing

The pressure ratings represent typical reference values for selected dimensions. The tests were conducted on a test bench in accordance with ISO 1402 at ambient conditions.

The values are applicable for non-sterile as well as sterilized tubing.

| Tube size ID (") x OD (") | Failure Pressure [bar] | Max. recommended working pressure [bar] |
|------------------------------|------------------------|--|
| 1/8 x 0,355 | ≥ 40 | 17 |
| 1/4 x 1/2 | ≥ 35 | 12 |
| 3/8 x 5/8 | ≥ 30 | 12 |
| 1/2 x 7/8 | ≥ 25 | 9 |
| 3/4 x 1-1/8 | ≥ 15 | 6 |
| 1 x 1-3/8 | ≥ 10 | 4 |

The specified pressures were determined on tubes in the condition as they are delivered by RAUMEDIC. With regards to pressure testing of silicone tubing RAUMEDIC does not use the common terminology of burst but failure pressures. At this specific point, the tube is no longer able build up pressure but inflates uncontrollably due to defect patterns generated by the applied test conditions. Failure pressures and working pressures of silicone tubes can be significantly influenced by several determinants like flow media, operating temperature and time, shore hardness, connection technology and sterilization conditions. For this reason, the failure pressure is to be considered as typical reference value while representing the minimum requirement for each dimension. The maximum working pressures as recommended by RAUMEDIC represent values that have been determined considering the influencing factors mentioned above to the best of our knowledge. The failure and working pressures stated by RAUMEDIC do not release the customer from process validation under application conditions since it is not possible to test all variables within the scope of the type tests performed.

Additional Documents

A.1 Certificates – ISO 13485:2016

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證書 ◆ CERTIFICADO ◆ CERTIFIKAT ◆ CERTIFICATE



Product Service

Certificate

No. Q5 053268 0081 Rev. 01

Holder of Certificate: RAUMEDIC AG
Hermann-Staudinger-Strasse 2
95233 Helmbrechts
GERMANY

Certification Mark:



Scope of Certificate: Design, development, production and sales of
- systems, components and semi-finished products made from polymer materials for medical devices and medical accessories, based on extrusion, injection moulding and assembly techniques
- precision measurement catheters and accessories
- compounds for the manufacturing of products for medical applications

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsud.com/ps-cert?q=cert:Q5 053268 0081 Rev. 01

Report No.: 713204491

Valid from: 2022-04-06

Valid until: 2025-03-31

Date, 2022-04-06

Christoph Dicks
Head of Certification/Notified Body



Certificate

No. Q5 053268 0081 Rev. 01

Applied Standard(s): EN ISO 13485:2016
 Medical devices - Quality management systems -
 Requirements for regulatory purposes
 (ISO 13485:2016)
 DIN EN ISO 13485:2016

Facility(ies): RAUMEDIC AG
 Hermann-Staudinger-Strasse 2, 95233 Helmbrechts, GERMANY

 See Scope of Certificate

 RAUMEDIC AG
 Crailsheimer Strasse 34, 91555 Feuchtwangen, GERMANY

 See Scope of Certificate

 RAUMEDIC AG
 Am Mühlgraben 10, 08297 Zwönitz, GERMANY

 See Scope of Certificate

Parameters: -

A.2 Certificates – ISO 15378:2018



CERTIFICATE



It is hereby certified, that the company

Raumedic AG

Hermann-Staudinger-Straße 2
95233 Helmbrechts, Germany

has implemented and applied a Management System in accordance to the
Standard

DIN EN ISO 15378:2018

Certified Scope:

Design, development and production of systems, components, semi finished products made from polymer materials as primary packaging materials for medicinal products based on extrusion, injection moulding and assembly techniques.

This certificate is valid from 01-15-2021 until 01-14-2024.

certificate-no.: K1068/PMA/10.20

Evidence has been provided with audit report no. K1068/10054.

bavaria certification GmbH
Oberschneiding, 01-15-2021



bavaria certification GmbH
Straubinger Straße 19
94363 Oberschneiding
www.bavaria-cert.com

A.4 Documentation – Material Safety Data Sheet

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9) Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878
Issue date: 12/5/2023 Revision date: 12/5/2023 Supersedes version of: 1/15/2020 Version: 2.00

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name : SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)
Type of product : Polymer

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Use of the substance/mixture : Extrusion, Injection molding

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

RAUMEDIC AG
Hermann Staudinger-Straße 2
DE– 95233 Helmbrechts
Germany
T +49-(0) 92 52 / 3 59-0 - F +49-(0) 92 52 / 3 59-10 00
info@RAUMEDIC.com - www.RAUMEDIC.com

Email competent person

sds@kft.de

1.4. Emergency telephone number

Emergency number : GIZ-Nord, Göttingen
Germany
+49 551 19240

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

To our knowledge, this product does not present any particular risk, provided it is handled in accordance with good occupational hygiene and safety practice.

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

2.3. Other hazards

PBT: not relevant – no registration required

vPvB: not relevant – no registration required

Contains no PBT and/or vPvB substances $\geq 0.1\%$ assessed in accordance with REACH Annex XIII

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

3.2. Mixtures

Comments : Vinyl methyl dimethyl polysiloxane using silicic acid fillers
Safety data sheet is for the processed (tempered) product.
This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of REACH Annex II

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general : In all cases of doubt, or when symptoms persist, seek medical attention.
First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing. Inhalation of product is unlikely.
First-aid measures after skin contact : Wash skin with plenty of water. Get medical advice if skin irritation persists.
First-aid measures after eye contact : Eye contact is unlikely. If particles get into the eyes, remove as customary for foreign bodies. Rinse eyes with water as a precaution.
First-aid measures after ingestion : Spit. Rinse mouth out with water. Do not induce vomiting. Call a poison center or a doctor if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire. Water spray. Dry powder. Foam. Carbon dioxide.
Unsuitable extinguishing media : Strong water jet.

5.2. Special hazards arising from the substance or mixture

Reactivity in case of fire : Decompose in hydrogen which is fire and explosion generator.
Hazardous decomposition products in case of fire : Toxic fumes may be released. Carbon monoxide. Carbon dioxide. Silicon oxide. hydrogen.

5.3. Advice for firefighters

Firefighting instructions : If there is a fire close by, use suitable extinguishing agents. Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire fighting water from entering the environment.
Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.
Other information : Do not allow run-off from fire fighting to enter drains or water courses. Disposal must be done according to official regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Wear breathing apparatus if exposed to vapours/dusts/aerosols.

6.1.1. For non-emergency personnel

Emergency procedures : Prohibit unauthorized persons. Ventilate spillage area. Do not breathe mist, vapours, spray. Keep away from sources of ignition. Spill area may be slippery.

6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

6.2. Environmental precautions

Avoid release to the environment. Avoid sub-soil penetration. Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up liquid spill into absorbent material. Take up mechanically (sweeping, shovelling) and collect in suitable container for disposal. Clean contaminated surfaces with an excess of water.

Other information : Disposal must be done according to official regulations.

6.4. Reference to other sections

Information for safe handling. See section 7. Concerning personal protective equipment to use, see section 8. For further information refer to section 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Additional hazards when processed : Product may release hydrogen gas. Increased storage temperatures will accelerate this process.

Precautions for safe handling : Ensure good ventilation of the work station. Wear personal protective equipment. Do not breathe mist, vapours, spray. Keep away from sources of ignition - No smoking. Avoid the build-up of electrostatic charge.

Hygiene measures : Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in original container. Store in a well-ventilated place. Keep cool. Protect from moisture.

Incompatible materials : Do not store in brand new glass containers with an alkaline surface.

Heat and ignition sources : Keep away from heat and direct sunlight.

Information about storage in one common storage facility : Keep away from food, drink and animal feeding stuffs.

7.3. Specific end use(s)

No additional information available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1 National occupational exposure and biological limit values

No additional information available

8.1.2. Recommended monitoring procedures

No additional information available

8.1.3. Air contaminants formed

No additional information available

8.1.4. DNEL and PNEC

No additional information available

8.1.5. Control banding

No additional information available

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

8.2. Exposure controls

8.2.1. Appropriate engineering controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

8.2.2. Personal protection equipment

8.2.2.1. Eye and face protection

Eye protection:

Use splash goggles when eye contact due to splashing is possible. ISO 16321-1

8.2.2.2. Skin protection

Skin and body protection:

Wear suitable protective clothing. EN 13034. EN ISO 13688

Hand protection:

In case of repeated or prolonged contact wear gloves. Nitrile rubber. Butyl rubber. ISO 374-1. Choosing the proper glove is a decision that depends not only on the type of material, but also on other quality features, which differ for each manufacturer. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer. Gloves must be replaced after each use and whenever signs of wear or perforation appear

8.2.2.3. Respiratory protection

Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Breathing apparatus with filter. P2. EN 143. Short term exposure. Breathing equipment is only to be used in order to handle the residual risk of short term jobs if all other risk minimizing measures have been carried out e.g. retention and/or local exhaust.

8.2.2.4. Thermal hazards

No additional information available

8.2.3. Environmental exposure controls

Environmental exposure controls:

Avoid release to the environment.

Other information:

Avoid contact with skin and eyes. Do not eat, drink or smoke when using this product. Wash hands before breaks and after work.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

| | |
|---|-----------------------------|
| Physical state | : Liquid |
| Colour | : Diverse. |
| Appearance | : Paste. |
| Odour | : weak. |
| Odour threshold | : Not available |
| Melting point | : Not applicable |
| Freezing point | : Not available |
| Boiling point | : Not available |
| Flammability | : Not applicable |
| Explosive properties | : Product is not explosive. |
| Oxidising properties | : Non oxidizing. |
| Lower explosion limit | : Not available |
| Upper explosion limit | : Not available |
| Flash point | : > 150 °C |
| Auto-ignition temperature | : Not available |
| Decomposition temperature | : Not available |
| pH | : Not available |
| Viscosity, kinematic | : Not available |
| Solubility | : Water: Insoluble |
| Partition coefficient n-octanol/water (Log Kow) | : Not available |
| Vapour pressure | : Not available |

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

| | |
|---------------------------------|------------------|
| Vapour pressure at 50°C | : Not available |
| Density | : Not available |
| Relative density | : Not available |
| Relative vapour density at 20°C | : Not available |
| Particle characteristics | : Not applicable |

9.2. Other information

9.2.1. Information with regard to physical hazard classes

No additional information available

9.2.2. Other safety characteristics

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use. Contact with metals produces hydrogen gas which may form explosive mixtures with air.

10.4. Conditions to avoid

Protect from moisture. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

10.5. Incompatible materials

Strong acids. alkalis. Oxidizing agent. Alcohol. Water, humidity.

10.6. Hazardous decomposition products

Hydrogen. Formaldehyde.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

| | |
|-----------------------------------|---|
| Acute toxicity (oral) | : Not classified (Based on available data, the classification criteria are not met) |
| Acute toxicity (dermal) | : Not classified (Based on available data, the classification criteria are not met) |
| Acute toxicity (inhalation) | : Not classified (Based on available data, the classification criteria are not met) |
| Skin corrosion/irritation | : Not classified (Based on available data, the classification criteria are not met) |
| Serious eye damage/irritation | : Not classified (Based on available data, the classification criteria are not met) |
| Respiratory or skin sensitisation | : Not classified (Based on available data, the classification criteria are not met) |
| Germ cell mutagenicity | : Not classified (Based on available data, the classification criteria are not met) |
| Carcinogenicity | : Not classified (Based on available data, the classification criteria are not met) |
| Reproductive toxicity | : Not classified (Based on available data, the classification criteria are not met) |
| STOT-single exposure | : Not classified (Based on available data, the classification criteria are not met) |
| STOT-repeated exposure | : Not classified (Based on available data, the classification criteria are not met) |
| Aspiration hazard | : Not classified (Based on available data, the classification criteria are not met) |

11.2. Information on other hazards

No additional information available

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

SECTION 12: Ecological information

12.1. Toxicity

Hazardous to the aquatic environment, short-term (acute) : Not classified (Based on available data, the classification criteria are not met)

Hazardous to the aquatic environment, long-term (chronic) : Not classified (Based on available data, the classification criteria are not met)

12.2. Persistence and degradability

| SIK5XXX, SIK6XXX, SIK8XXX (X=0-9) | |
|-----------------------------------|---|
| Persistence and degradability | Product is practically not biodegradable. |

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

| SIK5XXX, SIK6XXX, SIK8XXX (X=0-9) | |
|---|--|
| PBT: not relevant – no registration required | |
| vPvB: not relevant – no registration required | |

12.6. Endocrine disrupting properties

No additional information available

12.7. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Disposal must be done according to official regulations. Do not dispose of with domestic waste. Do not discharge into drains or the environment. European waste catalogue.
Product/Packaging disposal recommendations : Recycle or dispose of in compliance with current legislation.
Additional information : Packaging must be completely emptied.

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA / ADN / RID

| ADR | IMDG | IATA | ADN | RID |
|---|---------------|---------------|---------------|---------------|
| 14.1. UN number or ID number | | | | |
| Not regulated for transport | | | | |
| 14.2. UN proper shipping name | | | | |
| Not regulated | Not regulated | Not regulated | Not regulated | Not regulated |
| 14.3. Transport hazard class(es) | | | | |
| Not regulated | Not regulated | Not regulated | Not regulated | Not regulated |
| 14.4. Packing group | | | | |
| Not regulated | Not regulated | Not regulated | Not regulated | Not regulated |
| 14.5. Environmental hazards | | | | |
| Not regulated | Not regulated | Not regulated | Not regulated | Not regulated |
| No supplementary information available | | | | |

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

14.6. Special precautions for user

Overland transport

Not regulated

Transport by sea

Not regulated

Air transport

Not regulated

Inland waterway transport

Not regulated

Rail transport

Not regulated

14.7. Maritime transport in bulk according to IMO instruments

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Other information, restriction and prohibition regulations : A safety data sheet is not required for this product under Article 31 of REACH. This Product Safety Information Sheet has been created on a voluntary basis.

REACH Annex XVII (Restriction List)

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

REACH Annex XIV (Authorisation List)

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

REACH Candidate List (SVHC)

Contains no substance(s) listed on the REACH Candidate List

PIC Regulation (Prior Informed Consent)

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

POP Regulation (Persistent Organic Pollutants)

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

Explosives Precursors Regulation (2019/1148)

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

Drug Precursors Regulation (273/2004)

Contains no substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

15.1.2. National regulations

Germany

National Rules and Recommendations : TRGS 510: Storage of hazardous substances in non-stationary containers.
Water hazard class (WGK) : WGK 1, Slightly hazardous to water (Classification according to AwSV, Annex 1).
Storage class (LGK, TRGS 510) : LGK 10 - Combustible liquids.
Hazardous Incident Ordinance (12. BImSchV) : Is not subject of the Hazardous Incident Ordinance (12. BImSchV)

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

SECTION 16: Other information

Indication of changes

| Section | Changed item | Change | Comments |
|---------|------------------|--------|----------|
| | General revision | | |

Abbreviations and acronyms:

| | |
|---------|---|
| ADN | European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways |
| ADR | European Agreement concerning the International Carriage of Dangerous Goods by Road |
| ATE | Acute Toxicity Estimate |
| BCF | Bioconcentration factor |
| CLP | Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008 |
| DMEL | Derived Minimal Effect level |
| DNEL | Derived-No Effect Level |
| EC50 | Median effective concentration |
| IARC | International Agency for Research on Cancer |
| IATA | International Air Transport Association |
| IMDG | International Maritime Dangerous Goods |
| LC50 | Median lethal concentration |
| LD50 | Median lethal dose |
| LOAEL | Lowest Observed Adverse Effect Level |
| NOAEC | No-Observed Adverse Effect Concentration |
| NOAEL | No-Observed Adverse Effect Level |
| NOEC | No-Observed Effect Concentration |
| OECD | Organisation for Economic Co-operation and Development |
| PBT | Persistent Bioaccumulative Toxic |
| PNEC | Predicted No-Effect Concentration |
| REACH | Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006 |
| RID | Regulations concerning the International Carriage of Dangerous Goods by Rail |
| SDS | Safety Data Sheet |
| STP | Sewage treatment plant |
| TLM | Median Tolerance Limit |
| vPvB | Very Persistent and Very Bioaccumulative |
| CAS-No. | Chemical Abstract Service number |

Data sources : European Chemicals Agency, <http://echa.europa.eu/>. Information provided by the manufacturer.

Department issuing data specification sheet: : KFT Chemieservice GmbH
Im Leuschnerpark 3
D-64347 Griesheim

Phone: +49 6155-8981-400
Fax: +49 6155 8981-500
SDS Service: +49 6155 8981-522

Contact person : Dr. Christian Rank

KFT SDS EU 00 - Version 23.1

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

RAUMEDIC biopharma tubings offer significant advantages



RAUMEDIC AG

Hermann-Staudinger-Str. 2
95233 Helmbrechts
+49 9252 359-0
contact@raumedic.com

Item No. C86-N0044/11.2023

Validation Package for BRAIDMOTION
(TLV_002273_EN, Rev. 3.0)

Court of Registration: Amtsgericht Hof: HRB 3643

Executive Board: Stefan Seufferling (Chairman) – Thomas Knechtel – Klaus Schabert – Dr. Robert Schilling

Chairman of the Supervisory Board: Jürgen Werner