

# Validation Package

SILMOTION

Silicone Pharma Tubing  
made by RAUMEDIC

<b>1. Introduction</b>	<b>04</b>
1.1 About RAUMEDIC	04
1.2 Security of Supply	05
1.3 cGMP Conformity	05
1.4 Quality Management System	05
<b>2. Technical Specifications</b>	<b>06</b>
2.1 Tubing Dimensions, Type and Part Numbers Overview	06
2.2 Material Information	07
2.3 Material Compliance	07
2.4 Physical Properties	08
2.5 Physico-Chemical Properties	08
2.6 Biocompatibility	09
2.7 Tubing Printing	10
2.8 Extractable Profile	10
2.9 Sterilization Compatibility	10
2.10 General Properties, Shelf Life and Storage	11
<b>3. Functional Tests</b>	<b>12</b>
3.1 Pressure Testing	12
3.2 Pumping Lifetime	13
3.3 Flow Rate Data	13
<b>A. Additional Documents</b>	<b>14</b>
A.1 Certificates – ISO 13485:2016	14
A.2 Certificates – ISO 15378:2018	16
A.3 Certificates – ISO 50001:2018	17
A.4 Documentation – Material Safety Data Sheet	19



For us the issue of the customer's approval is the proof of the functional efficiency respectively of the fitness of our articles, which the customer determined for the specification application within the scope of testing. The obligated quality of goods has to conform to the expressly stipulated performance features in written form only (e.g. Specifications, Technical Delivery Specification, drawings, markings and other data).

We only assume a warranty exceeding this stipulation of obligated quality for a specific application or a specific fitness, period of application or life

period after passing of risk, if this is expressly agreed upon in written form; otherwise the risk of application and fitness is solely with the purchaser.

RAUMEDIC products and medical devices manufactured from these may not be sold or used in applications in the human body lasting more than 24 hours. No other guarantees are or will be given by RAUMEDIC.

This document is valid until revoked or until new information is provided. If you require further information, please contact your responsible sales contact.

# 1. Introduction

## 1.1 About RAUMEDIC

As a partner of the international medical technology and pharmaceutical industry, RAUMEDIC develops and produces components for customers, including tubing, catheters, and molded parts as well as complex groups of components and systems for a wide range of diagnostic and therapeutic uses. For the clinical areas of neuro-monitoring and traumatology, RAUMEDIC produces high-precision pressure-measurement sensors. The company is 100 percent family-owned.

### RAUMEDIC Group

**RAUMEDIC AG**



**The information in this document is subject to change without notice and should not be construed as a commitment by RAUMEDIC.**

## **1.2 Security of Supply**

The RAUMEDIC manufacturing plants operate under ISO 7 clean room conditions in order to offer flexibility and reliability to our customers worldwide. Consistent process performance is ensured by the ongoing qualification of all manufacturing processes and personnel.

RAUMEDIC pursues a multiple sourcing strategy for the Silicone Pharma Grade Silicone SIK8649. Different suppliers are used for this grade to ensure the security of supply. RAUMEDIC has validated the mechanical, physical and chemical equivalence and will use the sources in parallel according to the market situation. The used sources can be traced back via batch traceability at any time and fulfill the agreed specification.

## **1.3 cGMP Conformity**

Consistent high quality of silicone tubing is assured by careful selection of the raw material, well planned and validated production technologies and an efficient Quality Assurance Department, all of which result in high batch-to-batch reproducibility.

## **1.4 Quality Management System**

RAUMEDIC is certified according to ISO 13485. The respective certificate can be found attached to this document in appendix A.1 Certificates.

# 2.0 Technical Specifications

## 2.1 Tubing Dimensions, Type and Part Numbers Overview

SILMOTION tubing coils made of Silicone SIK8649 are of dimensions between 1/32" x 5/32" (ID x OD) and 1" x 1 3/8" (ID x OD). The color of the tubing is natural without addition of color pigments (RAUMEDIC color no. 20900). All tubing coils are provided in double PE bags packaged in a cardboard box. Further dimensions upon request.

Dimension ID x OD	Tolerances		Coil length [m]	Printed / non-printed	Article code
	ID [mm]	WT (*OD) [mm]			
1/32" x 5/32" 0,8 mm x 4,0 mm	0,70 - 0,80	1,50 - 1,70	100	non-printed	818830-800
1/16" x 3/16" 1,6 mm x 4,8 mm	1,45 - 1,75	1,45 - 1,75	15	non-printed	818832-323
			15	printed	818832-324
			100	non-printed	818832-800
			100	printed	818832-801
3/32" x 7/32" 2,4 mm x 5,6 mm	2,20 - 2,60	1,40 - 1,80	100	non-printed	818816-800
1/8" x 1/4" 3,2 mm x 6,4 mm	3,00 - 3,40	1,40 - 1,80	15	non-printed	818834-323
			15	printed	818834-324
			100	non-printed	818834-800
			100	printed	818834-801
1/8" x 5/16" 3,2 mm x 7,9 mm	3,00 - 3,40	2,10 - 2,70	50	non-printed	818819-800
3/16" x 5/16" 4,8 mm x 7,9 mm	4,55 - 4,75	* 7,50 - 7,92	50	non-printed	818847-800
3/16" x 3/8" 4,8 mm x 9,5 mm	4,60 - 5,00	2,10 - 2,70	50	non-printed	819040-800
1/4" x 3/8" 6,4 mm x 9,5 mm	6,10 - 6,70	1,40 - 1,80	15	non-printed	818835-323
			15	printed	818835-324
			50	non-printed	818835-800
			50	printed	818835-801
1/4" x 7/16" 6,4 mm x 11,1 mm	6,10 - 6,70	2,10 - 2,70	15	non-printed	818841-323
			15	printed	818841-324
			50	non-printed	818841-800
			50	printed	818841-801
5/16" x 1/2" 7,9 mm x 12,7 mm	7,62 - 7,92	2,20 - 2,60	15	non-printed	818035-323
			15	printed	818035-324
			30	non-printed	818035-800
			30	printed	818035-801
5/16" x 3/4" 7,9 mm x 19,1 mm	7,62 - 7,92	5,30 - 5,90	25	non-printed	818806-800
3/8" x 1/2" 9,5 mm x 12,7 mm	9,40 - 9,80	1,40 - 1,80	50	non-printed	818880-800
3/8" x 9/16" 9,5 mm x 14,3 mm	9,40 - 9,80	2,10 - 2,70	25	non-printed	818762-800

Dimension ID x OD	Tolerances		Coil length [m]	Printed / non-printed	Article code
	ID [mm]	WT (*OD) [mm]			
3/8" x 5/8" 9,5 mm x 15,9 mm	9,40 - 9,80	2,90 - 3,30	15	non-printed	818842-323
			15	printed	818842-324
			25	non-printed	818842-800
			25	printed	818842-801
1/2" x 3/4" 12,7 mm x 19,1 mm	12,50 - 13,50	2,70 - 3,30	15	non-printed	818844-323
			15	printed	818844-324
			25	non-printed	818844-800
			25	printed	818844-801
1/2" x 7/8" 12,7 mm x 22,2 mm	12,20 - 12,70	4,45 - 4,75	20	non-printed	818828-800
3/4" x 1" 19,1 mm x 25,4 mm	18,50 - 19,70	2,85 - 3,45	15	non-printed	818839-323
			15	non-printed	818839-800
			15	printed	818839-801
3/4" x 1 1/8" 19,1 mm x 28,6 mm	18,50 - 19,70	4,45 - 5,05	15	non-printed	818840-323
			15	non-printed	818840-800
			15	printed	818840-801
1" x 1 3/8" 25,4 mm x 34,9 mm	24,60 - 26,20	4,40 - 5,20	10	non-printed	818846-323
			10	non-printed	818846-800
			10	printed	818846-801

## 2.2 Material Information

SILMOTION silicone pharma tubing is made of the material grade Silicone SIK8649, a high purity silicone rubber suitable for pharma applications. Chemically Silicone SIK8649 is an addition cross-linked hot vulcanisate based on vinyl methyl dimethyl polysiloxane using silicic acid fillers and a platinum catalyst.

SILMOTION tubing can be used continuously in a temperature range from -60 °C (-76 °F) to +200 °C (+392 °F) without losing its integrity or deterioration of its chemical | physiological properties

The material Silicone SIK8649 shows good resistance to water up to 100 °C and to low-pressure steam up to 135 °C. Steam from higher temperatures leads to degradation, especially after prolonged exposure.

Silicone SIK8649 has good resistance to weak acids and alkalis. However, the material will degrade by strong acids and alkalis, which will be promoted at higher temperatures. The resistance is strongly dependent on the polar or non-polar character of the contact medium. While there is almost no swelling in polar contact media (e.g. water/alcohol), non-polar contact media (e.g. petrol/oil) cause medium to strong but reversible swelling.

The surface is coated in a plasma process. This coating provides a less sticky surface of these silicone tubing (Low Tack - see also 2.7) in comparison with common non-coated silicone products. Silicone tubing is colourless transparent or translucent.

## 2.3 Material Compliance

Based on our experience and knowledge about the materials and production processes the material formulation complies with / does not contain:

- Bisphenol A (BPA)
- Latex
- Phthalates
- Endocrine disruptors
- Ozone depleting substances (Regulation (EC) No 1005/2009)
- Perfluorinated chemicals (PFOA, PFOS)
- Persistent Bioaccumulative Toxins (PBTs)
- Persistent Organic Pollutants (POPs) (Regulation (EU) 2019/1021)
- Genetically Modified Organisms (GMO)
- Animal derived materials (Materials which present BSE/TSE risks)
- Conflict minerals (Dodd-Frank Act, Section 1502)
- RoHS (Directive 2011/65/EU, Directive (EU) 2015/863)
- Heavy Metals (CONEG, Directive 94/62/EC)
- CMR substances (Regulation (EC) No 1272/2008)
- Substances of Very High Concern (SVHC) as identified under article 33 of Regulation (EC) No 1907/2006 (REACH or substances listed in Annex XIV and Annex XVII)
- Nanomaterials (acc. to the definition Medical Device Regulation (EU) 2017/745) as intentionally added components above threshold limits stated in aforementioned Regulations, Directives or Laws

However, the raw materials used, and our final product(s) are not analyzed for the presence of traces of the above-mentioned substances.

**This statement is valid until revoked or until new information is provided.**

If you require further information, please contact your responsible sales contact.

## 2.4 Physical Properties

Property	Norm	SIK8649
Shore hardness A	ISO 48-4	60 ± 5
Density (g/cm <sup>3</sup> )	ISO 1183	1.15
Tear strength (MPa)	ISO 527	≥ 8.0
Elongation at break (%)	ISO 527	≥ 500
Compression Set (% , 22h / 175 °C)	ISO 815-1 type B	≤ 15 %

The physical properties of the raw material were determined on standard test specimen. The specified values represent guide values.



## 2.5 Physico-Chemical Properties

European Pharmacopoeia: 3.1.9 & FDA regulation 21CFR, § 177.2600

**Test results** | SILMOTION tubing in silicone rubber meets the requirements of the European Pharmacopoeia 3.1.9. and regulation 21 CFR, § 177.2600. The test methods, limits and results are those described by the E.P. monograph and listed in the table below.

Test Description	E.P. 3.1.9 Limits	Results on silicone tubing: Pass or Fail
Appearance of Solution	Colourless	Pass
Acidity	≤ 2.5 mL NaOH 0.01M	Pass
Alkalinity	≤ 1.0 mL HCl 0.01M	Pass
Reducing Substances	≤ 1 mL	Pass
Substances Soluble in Hexane	≤ 15 mg	Pass
Volatile Matter	< 2 %	Pass
Mineral Oils	< 1 ppm	Pass
Platinum	< 30 ppm	Pass

## USP <661> - Containers, Physico-chemical Tests – Plastic

**Test results** | The SILMOTION tubing meets the USP <661> requirements when sterilized at 50 kGy with and without aging conditions corresponding to a shelf life of 5 years.

Test Description	USP <661> Limits	Results on silicone tubing: Pass or Fail
Non Volatile Residue	< 15 mg	Pass
Residue on Ignition	< 5 mg	Pass
Heavy Metals	< 1 ppm	Pass
Buffering Capacity	< 10 mL	Pass

## 2.6 Biocompatibility

Samples of SILMOTION pharma tubing were tested in accordance with the below mentioned standards to assess state-of-the-art biocompatibility requirements. This includes ISO 10993 testing as well as USP class VI tests to evaluate toxicity systemically, intracutaneously, and through implantation.

Testing	Standard	Result
Hemolysis	ISO 10993-4	Pass
In-vitro Cytotoxicity	ISO 10993-5	Pass
Implantation	ISO 10993-6	Pass
Sensitization	ISO 10993-10	Pass
Pyrogenicity	Eur. Ph.	Pass
USP Class VI – Intracutaneous Reactivity (Irritation Rabbit) <sup>1)</sup>	USP <88>	Pass
USP Class VI – Systemic Toxicity (Mice) <sup>1)</sup>	USP <88>	Pass
USP Class VI – Implantation 7 days (Rabbit) <sup>1)</sup>	USP <88>	Pass

<sup>1)</sup>Tests were performed on gamma irradiated samples (50 kGy)

## 2.7 Tubing Printing

The tubing is available unprinted, with a standardized printing (refer table point 2.1) or a customized printing.

The standardized printing text is defined as follows:

„SILMOTION A“ x B“ MADE BY RAUMEDIC LOT XXXXXXXX“

A – Inner diameter in inch

B – Outer diameter in inch

RAUMEDIC biopharma silicone tubing features the patented Low-Tack surface treatment providing enhanced technical and handling properties.

The biocompatibility of the ink was proven by successfully performed ISO 10993-5 and USP <88> testing.

## 2.8 Extractable Profile

The study on the extractables profile for SILMOTION pharma tubing was performed according to BPOG “BioPhorum Best Practices Guide for Extractables Testing of Polymeric Single-Use Components Used in biopharmaceutical Manufacturing” released in April 2020. The assessment on gamma-irradiated and autoclaved tubing will be provided upon request.

## 2.9 Sterilization Compatibility

It is the responsibility of the user to validate a sterilization process with autoclave for SILMOTION tubing.

SILMOTION tubing is produced for single use. The product will be delivered non-sterile.

Sterilization with gas (ethylene oxide, ETO), steam (up to 135 °C), gamma or X-rays (Dose: max. 50 kGy) is generally recommended. The suitability of the material must be evaluated by the user. Possible material damage depends on the sterilization conditions (temperature, pressure, energy, time, packaging unit and distance to the radiation source).

Multiple sterilization is not recommended, neither with different nor the same sterilization method.

However, sterilization of SILMOTION tubing with multiple CIP/SIP-cycles at 135 °C does not impair fit, form or function of the products. The amount of the repetitions and the process conditions used are within the validation responsibility of the user. The user is obliged to examine the re-sterilizability under application-related conditions, since the number of re-sterilization cycles may be influenced by the flow medium/ cleaning agent and/ or mechanical stress applied throughout application lifetime.

## **2.10 General Properties, Shelf Life and Storage**

According to the current state-of-the-art individual fisheyes due to raw material and processing, foreign material, dirt inclusions and air bubbles as well as contamination on the tubing surface, like intrinsic particles and fluff, cannot be completely excluded.

The following shelf life is valid for the non-sterile, originally packed products.

A slight discoloration of the material might appear over storage time which may vary from batch to batch in intensity. This material related aging effect is common and does neither affect functionality nor material properties of Silicone SIK8649 in the scope of specified physical properties as well as shelf-life.

The shelf life of the tubing made of RAUMEDIC Silicone SIK8649 in non-sterile condition and originally packed is 5 years. There is no substantial change of the specified physical, chemical and physiological parameters of the tubing material after that period of time.

The following storage conditions shall apply:

- Storage in original packaging
- Dry (Standard conditions 10 – 75% relative humidity)
- Temperature 15 - 25°C with short-term deviations 5 – 35°C
- Goods protected from moisture
- Goods protected from directly exposure to UV-light
- Not in the environment and in presence of chemicals, e.g. solvent and disinfectant

# 3. Functional Tests

## 3.1 Pressure Testing

The pressure ratings represent typical reference values for selected dimensions. The tests were conducted on a test bench in accordance with ISO 1402 at ambient conditions. The values are applicable for non-sterile as well as sterilized tubing.

Tube size	Failure Pressure	Max. recommended working pressure
ID (") x OD (")	[bar]	[bar]
1/16 x 3/16	≥ 6,0	≤ 1,5
1/8 x 1/4	≥ 5,0	
1/4 x 3/8	≥ 3,0	
1/4 x 7/16	≥ 4,0	
3/8 x 5/8	≥ 4,5	
1/2 x 3/4	≥ 4,5	
3/4 x 1	≥ 3,0	
3/4 x 1 1/8	≥ 3,0	
1 x 1 3/8	≥ 3,0	

The specified pressures were determined on tubes in the condition as they are delivered by RAUMEDIC. With regards to pressure testing of silicone tubing RAUMEDIC does not use the common terminology of burst but failure pressures. The reason lies within the elastomeric behavior of silicone, which causes a balloon-like inflation of the tube after reaching the failure pressure. At this specific point, the tube is no longer able build up pressure but inflates uncontrollably. Failure pressures and working pressures of silicone tubes can be significantly influenced by several determinants like flow media, operating temperature and time, shore hardness, connection technology and sterilization conditions. For this reason, the failure pressure is to be considered as typical reference value while representing the minimum requirement for each dimension. The maximum working pressures as recommended by RAUMEDIC represent values that have been determined considering the influencing factors mentioned above to the best of our knowledge. The failure and working pressures stated by RAUMEDIC do not release the customer from process validation under application conditions since it is not possible to test all variables within the scope of the type tests performed.

## 3.2 Pumping Lifetime

Reference pumping lifetime for selected dimensions to assess the mechanical resistance of the tubing under pumping conditions. The stated values are based on external test data. The tests were conducted in a peristaltic roller pump. The tubing was pumping water at ambient temperature between 2 tanks mimicking recirculation conditions. The test was stopped and time measured at tubing break resulting in leak.

Tube size ID (") x OD (")	Pump Speed [rpm]	Pumping Lifetime [h]
1/4 x 7/16	210	> 70
1/2 x 3/4	310	
3/4 x 1 1/8	310	

## 3.3 Flow Rate Data

The objective of this test was to assess the maximum flow rate of selected dimensions by means of an peristaltic roller pump. The stated values are based on external test data.

Tube size ID (") x OD (")	Pump Speed [rpm]	Flow Rate [L/min] *
1/4 x 7/16	310	> 1,2
1/2 x 3/4		> 8,9
3/4 x 1 1/8		> 20,0

\*The measured values given are empirically determined guide values .

# Additional Documents

## A.1 Certificates – ISO 13485:2016

ZERTIFIKAT ◆ CERTIFICATE ◆ 認證書 ◆ CERTIFICADO ◆ CERTIFIKAT ◆ CERTIFICADO ◆ CERTIFICAT



Product Service

### Certificate

No. Q5 053268 0081 Rev. 01

**Holder of Certificate:** **RAUMEDIC AG**  
Hermann-Staudinger-Strasse 2  
95233 Helmbrechts  
GERMANY

**Certification Mark:**



**Scope of Certificate:** Design, development, production and sales of  
- systems, components and semi-finished products made from polymer materials for medical devices and medical accessories, based on extrusion, injection moulding and assembly techniques  
- precision measurement catheters and accessories  
- compounds for the manufacturing of products for medical applications

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: [www.tuvsud.com/ps-cert?q=cert:Q5 053268 0081 Rev. 01](http://www.tuvsud.com/ps-cert?q=cert:Q5_053268_0081_Rev_01)

**Report No.:** 713204491

**Valid from:** 2022-04-06

**Valid until:** 2025-03-31

**Date,** 2022-04-06

Christoph Dicks  
Head of Certification/Notified Body



# Certificate

No. Q5 053268 0081 Rev. 01

**Applied Standard(s):** EN ISO 13485:2016  
 Medical devices - Quality management systems -  
 Requirements for regulatory purposes  
 (ISO 13485:2016)  
 DIN EN ISO 13485:2016

**Facility(ies):** RAUMEDIC AG  
 Hermann-Staudinger-Strasse 2, 95233 Helmbrechts, GERMANY  
  
 See Scope of Certificate  
  
 RAUMEDIC AG  
 Crailsheimer Strasse 34, 91555 Feuchtwangen, GERMANY  
  
 See Scope of Certificate  
  
 RAUMEDIC AG  
 Am Mühlgraben 10, 08297 Zwönitz, GERMANY  
  
 See Scope of Certificate

**Parameters:** -



## A.2 Certificates – ISO 15378:2018



# CERTIFICATE



It is hereby certified, that the company

## Raumedic AG

Hermann-Staudinger-Straße 2  
95233 Helmbrechts, Germany

has implemented and applied a Management System in accordance to the  
Standard

## DIN EN ISO 15378:2018

Certified Scope:

Design, development and production of systems, components, semi finished products made from polymer materials as primary packaging materials for medicinal products based on extrusion, injection moulding and assembly techniques.

This certificate is valid from 01-15-2021 until 01-14-2024.

certificate-no.: K1068/PMA/10.20

Evidence has been provided with audit report no. K1068/10054.

bavaria certification GmbH  
Oberschneiding, 01-15-2021



bavaria certification GmbH  
Straubinger Straße 19  
94363 Oberschneiding  
www.bavaria-cert.com









## A.4 Documentation – Material Safety Data Sheet

### SIK5XXX, SIK6XXX, SIK8XXX (X=0-9) Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878  
Issue date: 12/5/2023 Revision date: 12/5/2023 Supersedes version of: 1/15/2020 Version: 2.00

#### SECTION 1: Identification of the substance/mixture and of the company/undertaking

##### 1.1. Product identifier

Trade name : SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)  
Type of product : Polymer

##### 1.2. Relevant identified uses of the substance or mixture and uses advised against

###### 1.2.1. Relevant identified uses

Use of the substance/mixture : Extrusion, Injection molding

###### 1.2.2. Uses advised against

No additional information available

##### 1.3. Details of the supplier of the safety data sheet

###### Supplier

RAUMEDIC AG  
Hermann Staudinger-Straße 2  
DE– 95233 Helmbrechts  
Germany  
T +49-(0) 92 52 / 3 59-0 - F +49-(0) 92 52 / 3 59-10 00  
[info@RAUMEDIC.com](mailto:info@RAUMEDIC.com) - [www.RAUMEDIC.com](http://www.RAUMEDIC.com)

###### Email competent person

sds@kft.de

##### 1.4. Emergency telephone number

Emergency number : GIZ-Nord, Göttingen  
Germany  
+49 551 19240

#### SECTION 2: Hazards identification

##### 2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

##### Adverse physicochemical, human health and environmental effects

To our knowledge, this product does not present any particular risk, provided it is handled in accordance with good occupational hygiene and safety practice.

##### 2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

##### 2.3. Other hazards

PBT: not relevant – no registration required

vPvB: not relevant – no registration required

Contains no PBT and/or vPvB substances  $\geq 0.1\%$  assessed in accordance with REACH Annex XIII

#### SECTION 3: Composition/information on ingredients

##### 3.1. Substances

Not applicable

# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

### 3.2. Mixtures

Comments : Vinyl methyl dimethyl polysiloxane using silicic acid fillers  
Safety data sheet is for the processed (tempered) product.  
This mixture does not contain any substances to be mentioned according to the criteria of section 3.2 of REACH Annex II

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures general : In all cases of doubt, or when symptoms persist, seek medical attention.  
First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing. Inhalation of product is unlikely.  
First-aid measures after skin contact : Wash skin with plenty of water. Get medical advice if skin irritation persists.  
First-aid measures after eye contact : Eye contact is unlikely. If particles get into the eyes, remove as customary for foreign bodies. Rinse eyes with water as a precaution.  
First-aid measures after ingestion : Spit. Rinse mouth out with water. Do not induce vomiting. Call a poison center or a doctor if you feel unwell.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire. Water spray. Dry powder. Foam. Carbon dioxide.  
Unsuitable extinguishing media : Strong water jet.

### 5.2. Special hazards arising from the substance or mixture

Reactivity in case of fire : Decompose in hydrogen which is fire and explosion generator.  
Hazardous decomposition products in case of fire : Toxic fumes may be released. Carbon monoxide. Carbon dioxide. Silicon oxide. hydrogen.

### 5.3. Advice for firefighters

Firefighting instructions : If there is a fire close by, use suitable extinguishing agents. Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire fighting water from entering the environment.  
Protection during firefighting : Do not attempt to take action without suitable protective equipment. Self-contained breathing apparatus. Complete protective clothing.  
Other information : Do not allow run-off from fire fighting to enter drains or water courses. Disposal must be done according to official regulations.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

General measures : Wear breathing apparatus if exposed to vapours/dusts/aerosols.

#### 6.1.1. For non-emergency personnel

Emergency procedures : Prohibit unauthorized persons. Ventilate spillage area. Do not breathe mist, vapours, spray. Keep away from sources of ignition. Spill area may be slippery.

#### 6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information refer to section 8: "Exposure controls/personal protection".

# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

### 6.2. Environmental precautions

Avoid release to the environment. Avoid sub-soil penetration. Prevent entry to sewers and public waters.

### 6.3. Methods and material for containment and cleaning up

- Methods for cleaning up : Take up liquid spill into absorbent material. Take up mechanically (sweeping, shovelling) and collect in suitable container for disposal. Clean contaminated surfaces with an excess of water.
- Other information : Disposal must be done according to official regulations.

### 6.4. Reference to other sections

Information for safe handling. See section 7. Concerning personal protective equipment to use, see section 8. For further information refer to section 13.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

- Additional hazards when processed : Product may release hydrogen gas. Increased storage temperatures will accelerate this process.
- Precautions for safe handling : Ensure good ventilation of the work station. Wear personal protective equipment. Do not breathe mist, vapours, spray. Keep away from sources of ignition - No smoking. Avoid the build-up of electrostatic charge.
- Hygiene measures : Do not eat, drink or smoke when using this product. Always wash hands after handling the product.

### 7.2. Conditions for safe storage, including any incompatibilities

- Storage conditions : Store in original container. Store in a well-ventilated place. Keep cool. Protect from moisture.
- Incompatible materials : Do not store in brand new glass containers with an alkaline surface.
- Heat and ignition sources : Keep away from heat and direct sunlight.
- Information about storage in one common storage facility : Keep away from food, drink and animal feeding stuffs.

### 7.3. Specific end use(s)

No additional information available

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### 8.1.1 National occupational exposure and biological limit values

No additional information available

#### 8.1.2. Recommended monitoring procedures

No additional information available

#### 8.1.3. Air contaminants formed

No additional information available

#### 8.1.4. DNEL and PNEC

No additional information available

#### 8.1.5. Control banding

No additional information available

# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

### 8.2. Exposure controls

#### 8.2.1. Appropriate engineering controls

**Appropriate engineering controls:**

Ensure good ventilation of the work station.

#### 8.2.2. Personal protection equipment

##### 8.2.2.1. Eye and face protection

**Eye protection:**

Use splash goggles when eye contact due to splashing is possible. ISO 16321-1

##### 8.2.2.2. Skin protection

**Skin and body protection:**

Wear suitable protective clothing. EN 13034. EN ISO 13688

**Hand protection:**

In case of repeated or prolonged contact wear gloves. Nitrile rubber. Butyl rubber. ISO 374-1. Choosing the proper glove is a decision that depends not only on the type of material, but also on other quality features, which differ for each manufacturer. Please follow the instructions related to the permeability and the penetration time provided by the manufacturer. Gloves must be replaced after each use and whenever signs of wear or perforation appear

##### 8.2.2.3. Respiratory protection

**Respiratory protection:**

In case of insufficient ventilation, wear suitable respiratory equipment. Breathing apparatus with filter. P2. EN 143. Short term exposure. Breathing equipment is only to be used in order to handle the residual risk of short term jobs if all other risk minimizing measures have been carried out e.g. retention and/or local exhaust.

##### 8.2.2.4. Thermal hazards

No additional information available

### 8.2.3. Environmental exposure controls

**Environmental exposure controls:**

Avoid release to the environment.

**Other information:**

Avoid contact with skin and eyes. Do not eat, drink or smoke when using this product. Wash hands before breaks and after work.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Colour	: Diverse.
Appearance	: Paste.
Odour	: weak.
Odour threshold	: Not available
Melting point	: Not applicable
Freezing point	: Not available
Boiling point	: Not available
Flammability	: Not applicable
Explosive properties	: Product is not explosive.
Oxidising properties	: Non oxidizing.
Lower explosion limit	: Not available
Upper explosion limit	: Not available
Flash point	: > 150 °C
Auto-ignition temperature	: Not available
Decomposition temperature	: Not available
pH	: Not available
Viscosity, kinematic	: Not available
Solubility	: Water: Insoluble
Partition coefficient n-octanol/water (Log Kow)	: Not available
Vapour pressure	: Not available



# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

Vapour pressure at 50°C	: Not available
Density	: Not available
Relative density	: Not available
Relative vapour density at 20°C	: Not available
Particle characteristics	: Not applicable

### 9.2. Other information

#### 9.2.1. Information with regard to physical hazard classes

No additional information available

#### 9.2.2. Other safety characteristics

No additional information available

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

### 10.2. Chemical stability

Stable under normal conditions.

### 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use. Contact with metals produces hydrogen gas which may form explosive mixtures with air.

### 10.4. Conditions to avoid

Protect from moisture. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

### 10.5. Incompatible materials

Strong acids. alkalis. Oxidizing agent. Alcohol. Water, humidity.

### 10.6. Hazardous decomposition products

Hydrogen. Formaldehyde.

## SECTION 11: Toxicological information

### 11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity (oral)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (dermal)	: Not classified (Based on available data, the classification criteria are not met)
Acute toxicity (inhalation)	: Not classified (Based on available data, the classification criteria are not met)
Skin corrosion/irritation	: Not classified (Based on available data, the classification criteria are not met)
Serious eye damage/irritation	: Not classified (Based on available data, the classification criteria are not met)
Respiratory or skin sensitisation	: Not classified (Based on available data, the classification criteria are not met)
Germ cell mutagenicity	: Not classified (Based on available data, the classification criteria are not met)
Carcinogenicity	: Not classified (Based on available data, the classification criteria are not met)
Reproductive toxicity	: Not classified (Based on available data, the classification criteria are not met)
STOT-single exposure	: Not classified (Based on available data, the classification criteria are not met)
STOT-repeated exposure	: Not classified (Based on available data, the classification criteria are not met)
Aspiration hazard	: Not classified (Based on available data, the classification criteria are not met)

### 11.2. Information on other hazards

No additional information available

# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

### SECTION 12: Ecological information

#### 12.1. Toxicity

Hazardous to the aquatic environment, short-term (acute) : Not classified (Based on available data, the classification criteria are not met)

Hazardous to the aquatic environment, long-term (chronic) : Not classified (Based on available data, the classification criteria are not met)

#### 12.2. Persistence and degradability

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)	
Persistence and degradability	Product is practically not biodegradable.

#### 12.3. Bioaccumulative potential

No additional information available

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)	
PBT: not relevant – no registration required	
vPvB: not relevant – no registration required	

#### 12.6. Endocrine disrupting properties

No additional information available

#### 12.7. Other adverse effects

No additional information available

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

Waste treatment methods : Disposal must be done according to official regulations. Do not dispose of with domestic waste. Do not discharge into drains or the environment. European waste catalogue.

Product/Packaging disposal recommendations : Recycle or dispose of in compliance with current legislation.

Additional information : Packaging must be completely emptied.

### SECTION 14: Transport information

In accordance with ADR / IMDG / IATA / ADN / RID

ADR	IMDG	IATA	ADN	RID
<b>14.1. UN number or ID number</b>				
Not regulated for transport				
<b>14.2. UN proper shipping name</b>				
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
<b>14.3. Transport hazard class(es)</b>				
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
<b>14.4. Packing group</b>				
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
<b>14.5. Environmental hazards</b>				
Not regulated	Not regulated	Not regulated	Not regulated	Not regulated
No supplementary information available				



# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

### 14.6. Special precautions for user

#### Overland transport

Not regulated

#### Transport by sea

Not regulated

#### Air transport

Not regulated

#### Inland waterway transport

Not regulated

#### Rail transport

Not regulated

### 14.7. Maritime transport in bulk according to IMO instruments

Not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Other information, restriction and prohibition regulations : A safety data sheet is not required for this product under Article 31 of REACH. This Product Safety Information Sheet has been created on a voluntary basis.

#### REACH Annex XVII (Restriction List)

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

#### REACH Annex XIV (Authorisation List)

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

#### REACH Candidate List (SVHC)

Contains no substance(s) listed on the REACH Candidate List

#### PIC Regulation (Prior Informed Consent)

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

#### POP Regulation (Persistent Organic Pollutants)

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

#### Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

#### Explosives Precursors Regulation (2019/1148)

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

#### Drug Precursors Regulation (273/2004)

Contains no substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

#### 15.1.2. National regulations

##### Germany

National Rules and Recommendations : TRGS 510: Storage of hazardous substances in non-stationary containers.  
Water hazard class (WGK) : WGK 1, Slightly hazardous to water (Classification according to AwSV, Annex 1).  
Storage class (LGK, TRGS 510) : LGK 10 - Combustible liquids.  
Hazardous Incident Ordinance (12. BImSchV) : Is not subject of the Hazardous Incident Ordinance (12. BImSchV)

### 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

# SIK5XXX, SIK6XXX, SIK8XXX (X=0-9)

## Safety Data Sheet

according to the REACH Regulation (EC) 1907/2006 amended by Regulation (EU) 2020/878

### SECTION 16: Other information

#### Indication of changes

Section	Changed item	Change	Comments
	General revision		

#### Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ATE	Acute Toxicity Estimate
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DMEL	Derived Minimal Effect level
DNEL	Derived-No Effect Level
EC50	Median effective concentration
IARC	International Agency for Research on Cancer
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
NOAEL	No-Observed Adverse Effect Level
NOEC	No-Observed Effect Concentration
OECD	Organisation for Economic Co-operation and Development
PBT	Persistent Bioaccumulative Toxic
PNEC	Predicted No-Effect Concentration
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
RID	Regulations concerning the International Carriage of Dangerous Goods by Rail
SDS	Safety Data Sheet
STP	Sewage treatment plant
TLM	Median Tolerance Limit
vPvB	Very Persistent and Very Bioaccumulative
CAS-No.	Chemical Abstract Service number

Data sources : European Chemicals Agency, <http://echa.europa.eu/>. Information provided by the manufacturer.

Department issuing data specification sheet: : KFT Chemieservice GmbH  
Im Leuschnerpark 3  
D-64347 Griesheim

Phone: +49 6155-8981-400  
Fax: +49 6155 8981-500  
SDS Service: +49 6155 8981-522

Contact person : Dr. Christian Rank

KFT SDS EU 00 - Version 23.1

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

RAUMEDIC biopharma tubings offer significant advantages



**RAUMEDIC AG**

Hermann-Staudinger-Str. 2  
95233 Helmbrechts  
+49 9252 359-0  
[contact@raumedic.com](mailto:contact@raumedic.com)

Item No. C86-N0044/12.2023

Validation Package for SILMOTION  
(TLV\_002255\_EN Rev. 6.0)

**Court of Registration:** Amtsgericht Hof: HRB 3643

**Executive Board:** Stefan Seufferling (Chairman) – Thomas Knechtel – Klaus Schabert – Dr. Robert Schilling

**Chairman of the Supervisory Board:** Jürgen Werner