

Ultraclean SA Technical Data Sheet



Description

Ultraclean SA is mild alkaline detergent for general use in ultrasonic baths:

- Aqueous based product
- General purpose cleaning
- Improves surface tension & enhances ultrasonic performance
- Contains de-foaming agent
- Suitable for use on aluminium

Applications

Used in a wide range of applications including instrument decontamination in healthcare and cleaning different components in manufacturing and processing industries. Ultraclean SA can be used to remove contaminants from plastics, glass, metals and other materials.

Directions for Use

Add 20-50ml of Ultraclean SA per 1 Litre of water in the ultrasonic cleaner. Do not mix with other solutions. The required dosage will depend on your cleaning application and the level of contamination. Items should be rinsed with water and dried immediately after processing.

Technical Data

These figures are typical and should not be interpreted as specifications, refer to the product safety data sheet for details

Characteristic	Value/Description
Appearance	Green liquid
Odour	Mild
pH (@ 1%)	9.8
Relative Density at 20°C	1.112

Supplied as:

- 6 x 1Litre Bottles (SKU: F50025)
- 25 Litre Drum (SKU: F50024)

Disposal

When diluted to concentration of 5% or less, Ultraclean SA is not hazardous to aquatic organisms and may be safely disposed of down the drain. Refer to regional or national regulations.

Storage

Store in a cool, dry place protected from frost and away from acids and strong oxidising agents. Store upright in original containers. Shelf-life up to 2years from date of manufacture.

Specifications are subject to change for continual improvements