

Thermo Scientific Finnpiquette Performance Specifications

Salla Ruosaari, Product Manager, Thermo Fisher Scientific, Vantaa, Finland.

Key words

- **ISO 8655** - Calibration standard for the piston operated volumetric devices
- **Factory Specifications** - Often more stringent requirements than the above ISO 8655, these are set by manufacturers in the development phase of a new product
- **Systematic error** - The difference between the dispensed volume and the selected volume of a pipette
- **Random error** - Error referring to the repeatability of pipetting
- **Pipetting Systems** - In liquid handling, pipettes and tips always form a pipetting system. Using a pipette with the pipette tip it was designed, developed and manufactured to fit to allow for ideal pipetting performance
- **ISO 9001** - Standard for Quality management systems
- **ISO 13485** - Quality management standard for medical devices

ISO 8655 Specifications

The ISO 8655 standard forms a guideline for industry and independent service companies concerning the entire calibration process of piston operated volumetric devices. The acceptance limits stated in the ISO standard are designed to be used by manufacturers, service companies as well as end users. The standard criteria are agreed to be sufficient for the use with the devices described therein.

The standard gives the maximum permissible limits for systematic and random error, also referred to as inaccuracy and imprecision, respectively. Systematic error is the difference between the dispensed volume and the selected volume of a pipette. The error is repeated in each measurement taken but can be eliminated if the error source, such as imperfect calibration or a change in the environment where the testing is performed, is identified. Random error, on the other hand, refers to the repeatability of the pipetting. It is caused by random factors including features of the pipette, laboratory practice and user experience.



Factory Specifications

All manufacturers, Thermo Fisher Scientific included, use their own factory performance specifications to describe the actual functionality of their pipettes. As the factory specifications are tighter than those required by the ISO 8655 standard, this is an acceptable practice in the industry. These specifications are provided for pipetting systems, which comprises both pipette and tip. Performance achievements may be different if using non-system tips. The Finnpiquette performance has been tested using their system tip which are Thermo Scientific Finntips.

Optimal Pipetting Performance with the Thermo Scientific Systems

The main difference between the factory and ISO 8655 specifications is that the factory error limits are testing volume dependent for variable volume pipettes. The ISO 8655 standard, on the other hand, gives the same error limits for the whole volume range of each pipette. The difference is especially noticeable in the lower volume error limits ensured by Thermo Scientific. For the mid and minimum volumes, the factory specifications are remarkably tighter than what is required by the ISO standard. With more stringent specifications, pipetting performance is guaranteed with the Finnpiettes.

To demonstrate the difference, ISO 8655 specifies a systematic error of 8 µL for a 1000 µL pipette regardless of the selected volume. Thermo Scientific F1 and F2 100-1000 µL pipettes are however given a systematic error of 1 µL at minimum and 6 µL at maximum volume. Thus, when pipetting a sample of 100 µL, the maximum permissible error according to Thermo Scientific factory specifications will be only 1% but 8% according to ISO 8655 (Table 1). The same applies for all Finnpiette models as demonstrated in Table 2

Quality matters

Your tip choice has a big impact on the pipetting performance, which is why Finnpiettes are designed for use with the Finntips. Both Finnpiettes and Finntips are manufactured from high quality raw materials, tested thoroughly and certified according to strict quality standards. Specifically, Thermo Fisher Scientific Oy is a company certified to the standards ISO 9001 and ISO 13485, taking quality matters seriously. Further, the pipetting systems comprising Finnpiettes and Finntips are CE-marked according to IVD Directive in Europe and intended for use in in-vitro diagnostics for transferring of liquids.

For further information on calibration of pipettes, refer to The Good Laboratory Pipetting Guide.

Volume	Thermo Scientific Factory Specifications		ISO 8655 Specifications	
	Inaccuracy µl (%)	Imprecision µl (%)	Inaccuracy µl (%)	Imprecision µl (%)
1000	±6,0 (±0,6%)	2,0 (0,2%)	±8,0 (±0,8%)	3,0 (0,3%)
500	±4,0 (±0,8%)	1,5 (0,3%)	±8,0 (±1,6%)	3,0 (0,6%)
100	±1,0 (±1,0%)	0,6 (0,6%)	±8,0 (±8,0%)	3,0 (3,0%)

Table 1. Comparison of the Thermo Scientific Factory Specifications, and ISO 8655 limits for Finnpiette F1 and F2 100-1000 µl models. Notice how remarkably the mid and min volume error limits differ between the Thermo Scientific Factory Specifications and ISO 8655.

Table 2. Comparison of the Thermo Scientific Factory Specifications and ISO 8655 limits.**Thermo Scientific Finnpiquette F1 Pipettes****Finnpiquette F1 Variable Volume Single Channel Pipettes**

Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4641010	0.2-2 µL	2	±0,050	0,040	±0,080	0,040
		1	±0,040	0,035	±0,080	0,040
		0,2	±0,024	0,020	±0,080	0,040
4641020	0.5-5 µL	5	±0,075	0,050	±0,125	0,075
		2,5	±0,063	0,040	±0,125	0,075
		0,5	±0,030	0,025	±0,125	0,075
4641030	1-10 µL	10	±0,100	0,050	±0,120	0,080
		5	±0,075	0,040	±0,120	0,080
		1	±0,025	0,020	±0,120	0,080
4641040	1-10 µL	10	±0,100	0,080	±0,120	0,080
		5	±0,075	0,040	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4641050	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,06	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4641060	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,06	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4641130	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4641140	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4641070	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4641080	20-200 µL	200	±1,2	0,4	±1,60	0,60
		100	±1,0	0,4	±1,60	0,60
		20	±0,36	0,14	±1,60	0,60
4641090	30-300 µL	300	±1,8	0,6	±4,0	1,5
		150	±1,5	0,6	±4,0	1,5
		30	±0,45	0,18	±4,0	1,5
4641100	100-1000 µL	1000	±6,0	2,0	±8,0	3,0
		500	±4,0	1,5	±8,0	3,0
		100	±1,0	0,6	±8,0	3,0
4641110	0,5-5 mL	5000	±25,0	10,0	±40,0	15,0
		2500	±17,5	7,5	±40,0	15,0
		500	±10,0	4,0	±40,0	15,0
4641120	1-10 mL	10000	±50,0	20,0	±60,0	30,0
		5000	±40,0	15,0	±60,0	30,0
		1000	±20,0	8,0	±60,0	30,0

Finnpipette F1 Fixed Volume Single Channel Pipettes

Cat. No.	Range	Volume (μL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (μL)	Imprecision (μL)	Inaccuracy (μL)	Imprecision (μL)
4651000	1 μL	1	$\pm 0,040$	0,040	$\pm 0,050$	0,050
4651010	5 μL	5	$\pm 0,070$	0,070	$\pm 0,125$	0,075
4651020	10 μL	10	$\pm 0,090$	0,080	$\pm 0,120$	0,080
4651130	20 μL	20	$\pm 0,14$	0,10	$\pm 0,20$	0,10
4651030	25 μL	25	$\pm 0,15$	0,125	$\pm 0,50$	0,20
4651040	50 μL	50	$\pm 0,30$	0,20	$\pm 0,50$	0,20
4651050	100 μL	100	$\pm 0,40$	0,30	$\pm 0,80$	0,30
4651140	200 μL	200	$\pm 0,80$	0,60	$\pm 1,60$	0,60
4651060	250 μL	250	$\pm 1,0$	0,75	$\pm 4,0$	1,5
4651070	500 μL	500	$\pm 1,5$	1,5	$\pm 4,0$	1,5
4651080	1000 μL	1000	$\pm 3,0$	3,0	$\pm 8,0$	3,0
4651090	2000 μL	2000	$\pm 6,0$	4,0	$\pm 16,0$	6,0
4651100	3000 μL	3000	$\pm 9,0$	6,0	$\pm 40,0$	15,0
4651150	4000 μL	4000	$\pm 12,0$	8,0	$\pm 40,0$	15,0
4651110	5000 μL	5000	$\pm 15,0$	10,0	$\pm 40,0$	15,0
4651120	10000 μL	10000	$\pm 30,0$	20,0	$\pm 60,0$	30,0



Finnpipette F1 Multichannel Pipettes

Cat. No.	Range	Channels	Volume (μL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
				Inaccuracy (μL)	Imprecision (μL)	Inaccuracy (μL)	Imprecision (μL)
4661000	1.0-10 μL	8	10	$\pm 0,240$	0,160	$\pm 0,240$	0,160
			5	$\pm 0,200$	0,150	$\pm 0,240$	0,160
			1	$\pm 0,120$	0,080	$\pm 0,240$	0,160
4661010	5-50 μL	8	50	$\pm 0,75$	0,35	$\pm 1,00$	0,40
			25	$\pm 0,625$	0,30	$\pm 1,00$	0,40
			5	$\pm 0,25$	0,10	$\pm 1,00$	0,40
4661020	10-100 μL	8	100	$\pm 1,30$	0,50	$\pm 1,60$	0,60
			50	$\pm 1,25$	0,60	$\pm 1,60$	0,60
			10	$\pm 0,50$	0,20	$\pm 1,60$	0,60
4661030	30-300 μL	8	300	$\pm 3,0$	0,9	$\pm 8,0$	3,0
			150	$\pm 2,25$	0,75	$\pm 8,0$	3,0
			30	$\pm 1,5$	0,6	$\pm 8,0$	3,0
4661040	1.0-10 μL	12	10	$\pm 0,240$	0,160	$\pm 0,240$	0,160
			5	$\pm 0,200$	0,150	$\pm 0,240$	0,160
			1	$\pm 0,120$	0,080	$\pm 0,240$	0,160
4661050	5-50 μL	12	50	$\pm 0,75$	0,35	$\pm 1,00$	0,40
			25	$\pm 0,625$	0,30	$\pm 1,00$	0,40
			5	$\pm 0,25$	0,10	$\pm 1,00$	0,40
4661060	10-100 μL	12	100	$\pm 1,30$	0,50	$\pm 1,60$	0,60
			50	$\pm 1,25$	0,60	$\pm 1,60$	0,60
			10	$\pm 0,50$	0,20	$\pm 1,60$	0,60
4661070	30-300 μL	12	300	$\pm 3,0$	0,9	$\pm 8,0$	3,0
			150	$\pm 2,25$	0,75	$\pm 8,0$	3,0
			30	$\pm 1,5$	0,6	$\pm 8,0$	3,0
4661080	1.0-10 μL	16	10	$\pm 0,240$	0,160	$\pm 0,240$	0,160
			5	$\pm 0,200$	0,150	$\pm 0,240$	0,160
			1	$\pm 0,120$	0,080	$\pm 0,240$	0,160
4661090	5-50 μL	16	50	$\pm 0,75$	0,35	$\pm 1,00$	0,40
			25	$\pm 0,625$	0,30	$\pm 1,00$	0,40
			5	$\pm 0,25$	0,10	$\pm 1,00$	0,40



Thermo Scientific Finnpipette F2 Pipettes

Finnpipette F2 Variable Volume Single Channel Pipettes



Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4642010	0.2-2 µL	2	±0,050	0,040	±0,080	0,040
		1	±0,040	0,035	±0,080	0,040
		0,2	±0,024	0,020	±0,080	0,040
4642020	0.5-5 µL	5	±0,075	0,050	±0,125	0,075
		2,5	±0,063	0,038	±0,125	0,075
		0,5	±0,030	0,025	±0,125	0,075
4642030	1-10 µL	10	±0,100	0,050	±0,120	0,080
		5	±0,075	0,040	±0,120	0,080
		1	±0,025	0,020	±0,120	0,080
4642040	1-10 µL	10	±0,100	0,080	±0,120	0,080
		5	±0,075	0,040	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4642050	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,06	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4642060	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,06	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4642120	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4642130	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4642070	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4642080	20-200 µL	200	±1,2	0,4	±1,60	0,60
		100	±1,0	0,4	±1,60	0,60
		20	±0,36	0,14	±1,60	0,60
4642090	100-1000 µL	1000	±6,0	2,0	±8,0	3,0
		500	±4,0	1,5	±8,0	3,0
		100	±1,0	0,6	±8,0	3,0
4642100	0,5-5 mL	5000	±25,0	10,0	±40,0	15,0
		2500	±17,5	7,5	±40,0	15,0
		500	±10,0	4,0	±40,0	15,0
4642110	1-10 mL	10000	±50,0	20,0	±60,0	30,0
		5000	±40,0	15,0	±60,0	30,0
		1000	±20,0	8,0	±60,0	30,0

Finnpipette F2 Fixed Volume Single Channel Pipettes

Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4652000	1 µL	1	±0,040	0,040	±0,050	0,050
4652010	5 µL	5	±0,070	0,070	±0,125	0,075
4652020	10 µL	10	±0,090	0,080	±0,120	0,080
4652130	20 µL	20	±0,14	0,10	±0,20	0,10
4652030	25 µL	25	±0,15	0,125	±0,50	0,20
4652040	50 µL	50	±0,30	0,20	±0,50	0,20
4652050	100 µL	100	±0,40	0,30	±0,80	0,30
4652140	200 µL	200	±0,80	0,60	±1,60	0,60
4652060	250 µL	250	±1,0	0,75	±4,0	1,5
4652070	500 µL	500	±1,5	1,5	±4,0	1,5
4652080	1000 µL	1000	±3,0	3,0	±8,0	3,0
4652090	2000 µL	2000	±6,0	4,0	±16,0	6,0
4652100	3000 µL	3000	±9,0	6,0	±40,0	15,0
4652150	4000 µL	4000	±12,0	8,0	±40,0	15,0
4652110	5000 µL	5000	±15,0	10,0	±40,0	15,0
4652120	10000 µL	10000	±30,0	20,0	±60,0	30,0



Finnpipette F2 Multichannel Pipettes

Cat. No.	Range	Channels	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
				Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4662000	1.0-10 µL	8	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4662010	5-50 µL	8	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4662020	10-100 µL	8	100	±1,30	0,50	±1,60	0,60
			50	±1,25	0,60	±1,60	0,60
			10	±0,50	0,20	±1,60	0,60
4662030	30-300 µL	8	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0
4662040	1.0-10 µL	12	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4662050	5-50 µL	12	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4662060	10-100 µL	12	100	±1,30	0,50	±1,60	0,60
			50	±1,25	0,60	±1,60	0,60
			10	±0,50	0,20	±1,60	0,60
4662070	30-300 µL	12	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0
4662080	1.0-10 µL	16	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4662090	5-50 µL	16	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40



Thermo Scientific Finnpipette F3 Pipettes*

Finnpipette F3 Variable Volume Single Channel Pipettes



Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4640000	1-10 µL	10	±0,100	0,050	±0,120	0,080
		5	±0,100	0,050	±0,120	0,080
		1	±0,025	0,020	±0,120	0,080
4640010	1-10 µL	10	±0,100	0,080	±0,120	0,080
		5	±0,100	0,065	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4640020	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,07	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4640030	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,07	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4640080	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4640090	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,13	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4640040	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4640050	20-200 µL	200	±1,2	0,6	±1,60	0,60
		100	±1,0	0,5	±1,60	0,60
		20	±0,6	0,3	±1,60	0,60
4640060	100-1000 µL	1000	±5,0	2,0	±8,0	3,0
		500	±5,0	2,0	±8,0	3,0
		100	±1,5	0,6	±8,0	3,0
4640100	0,5-5 mL	5000	±25,0	10,0	±40,0	15,0
		2500	±17,5	7,5	±40,0	15,0
		500	±10,0	4,0	±40,0	15,0
4640070	1-10 mL	10000	±50,0	20,0	±60,0	30,0
		5000	±40,0	15,0	±60,0	30,0
		1000	±20,0	8,0	±60,0	30,0

Finnpipette F3 Fixed Volume Single Channel Pipettes



Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4650000	1 µL	1	±0,040	0,040	±0,050	0,050
4650010	5 µL	5	±0,070	0,070	±0,125	0,075
4650020	10 µL	10	±0,090	0,080	±0,120	0,080
4650090	20 µL	20	±0,14	0,10	±0,20	0,10
4650030	25 µL	25	±0,15	0,125	±0,50	0,20
4650040	50 µL	50	±0,30	0,20	±0,50	0,20
4650050	100 µL	100	±0,40	0,30	±0,80	0,30
4650100	200 µL	200	±0,80	0,60	±1,60	0,60
4650060	250 µL	250	±1,0	0,8	±4,0	1,5
4650070	500 µL	500	±1,5	1,5	±4,0	1,5
4650080	1000 µL	1000	±3,0	3,0	±8,0	3,0

* Please note: Availability may vary by country.

Finnpipette F3 Multichannel Pipettes

Cat. No.	Range	Channels	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
				Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4660010	5-50 µL	8	50	±0,75	0,35	±1,00	0,40
			25	±0,63	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4660020	30-300 µL	8	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0
4660040	5-50 µL	12	50	±0,75	0,35	±1,00	0,40
			25	±0,63	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4660050	30-300 µL	12	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0



Thermo Scientific Finn timer Digital Pipettes*

Finn timer Digital Variable Volume Single Channel Pipettes

Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4500000	0.2-2 µL, micro	2	±0,050	0,040	±0,080	0,040
		1	±0,035	0,035	±0,080	0,040
		0,2	±0,024	0,020	±0,080	0,040
4500010	0.5-10 µL, micro	10	±0,100	0,050	±0,120	0,080
		5	±0,100	0,050	±0,120	0,080
		1	±0,025	0,020	±0,120	0,080
4500020	0.5-10 µL	10	±0,100	0,080	±0,120	0,080
		5	±0,100	0,050	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4500080	2-20 µL	20	±0,20	0,08	±0,20	0,10
		10	±0,15	0,07	±0,20	0,10
		2	±0,06	0,05	±0,20	0,10
4500100	5-50 µL	50	±0,30	0,15	±0,50	0,20
		25	±0,25	0,125	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4500110	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4500090	20-200 µL	200	±1,2	0,4	±1,60	0,60
		100	±1,0	0,4	±1,60	0,60
		20	±0,36	0,14	±1,60	0,60
4500120	100-1000 µL	1000	±6,0	2,0	±8,0	3,0
		500	±4,0	1,5	±8,0	3,0
		100	±1,0	0,6	±8,0	3,0
4500060	1-5 mL	5000	±25,0	10,0	±40,0	15,0
		2500	±17,5	7,5	±40,0	15,0
		1000	±8,0	3,0	±40,0	15,0
4500070	2-10 mL	10000	±50,0	20,0	±60,0	30,0
		5000	±40,0	15,0	±60,0	30,0
		2000	±20,0	6,0	±60,0	30,0



* Please note: Availability may vary by country.



Finnpipette Digital Multichannel Pipettes

Cat. No.	Range	Channels	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
				Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4510000	0.5-10 µL	8	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4510020	5-50 µL	8	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4510030	50-300 µL	8	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			50	±2,3	0,75	±8,0	3,0
4510010	0.5-10 µL	12	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4510040	5-50 µL	12	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4510050	50-300 µL	12	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			50	±2,3	0,75	±8,0	3,0
4510070	5-50 µL	16	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40

Thermo Scientific Finnpipette Novus Pipettes



Finnpipette Novus Single Channel Electronic Pipettes

Cat. No.	Range	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
			Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4620000	1-10 µL, micro	10	±0,100	0,050	±0,120	0,080
		5	±0,075	0,040	±0,120	0,080
		1	±0,035	0,030	±0,120	0,080
4620010	1-10 µL	10	±0,100	0,050	±0,120	0,080
		5	±0,075	0,040	±0,120	0,080
		1	±0,070	0,060	±0,120	0,080
4620020	5-50 µL, micro	50	±0,40	0,15	±0,50	0,20
		25	±0,25	0,125	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4620030	5-50 µL	50	±0,40	0,15	±0,50	0,20
		25	±0,25	0,125	±0,50	0,20
		5	±0,15	0,125	±0,50	0,20
4620040	10-100 µL	100	±0,80	0,20	±0,80	0,30
		50	±0,60	0,20	±0,80	0,30
		10	±0,30	0,10	±0,80	0,30
4620050	30-300 µL	300	±1,8	0,6	±4,0	1,5
		150	±1,5	0,6	±4,0	1,5
		30	±0,9	0,21	±4,0	1,5
4620060	100-1000 µL	1000	±6,0	2,0	±8,0	3,0
		500	±4,0	1,5	±8,0	3,0
		100	±3,0	0,6	±8,0	3,0
4620070	0.5-5 mL	5000	±30,0	10,0	±40,0	15,0
		2500	±20,0	10,0	±40,0	15,0
		500	±15,0	4,0	±40,0	15,0
4620080	1-10 mL	10000	±60,0	20,0	±60,0	30,0
		5000	±40,0	20,0	±60,0	30,0
		1000	±30,0	8,0	±60,0	30,0

Finnpipette Novus Multichannel Electronic Pipettes

Cat. No.	Range	Channels	Volume (µL)	Thermo Scientific Factory Calibration Specifications		ISO 8655 Calibration Specifications	
				Inaccuracy (µL)	Imprecision (µL)	Inaccuracy (µL)	Imprecision (µL)
4630000	1.0-10 µL	8	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4630020	5-50 µL	8	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4630040	30-300 µL	8	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0
4630080	100-1200 µL	8	1200	±12,0	2,4	±32,00	12,00
			600	±12,0	2,4	±32,00	12,00
			100	±3,0	0,9	±32,00	12,00
4630010	1.0-10 µL	12	10	±0,240	0,160	±0,240	0,160
			5	±0,200	0,150	±0,240	0,160
			1	±0,120	0,080	±0,240	0,160
4630030	5-50 µL	12	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40
4630050	30-300 µL	12	300	±3,0	0,9	±8,0	3,0
			150	±2,25	0,75	±8,0	3,0
			30	±1,5	0,6	±8,0	3,0
4630070	5-50 µL	16	50	±0,75	0,35	±1,00	0,40
			25	±0,625	0,30	±1,00	0,40
			5	±0,25	0,10	±1,00	0,40



* Please note: Availability may vary by country.

thermoscientific.com/liquidhandling

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: +1 800 522 7763
Europe:
Austria +43 1 801 40 0
Belgium +32 2 482 30 30
Finland/Nordic +358 9 329 100
France +33 2 28 03 20 00
Germany National Toll Free 08001-536 376

Germany International +49 6184 90 6940
Italy +39 02 95059 1
Netherlands +31 76 571 4440
Russia/CIS +7 (812) 703 42 15
Spain/Portugal +34 93 223 3154
Switzerland +41 44 454 12 12
UK/Ireland +44 870 609 9203

Asia:
India +91 22 5542 9494
Japan +81 45 453 9220
China +86 21 6865 4588 or +86 10 5850 3588
Other Asian countries +852 2885 4613
Countries not listed:
+49 6184 90 6940 or +33 2 28 03 20 00

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific