




art & science of
amazing protection



*Superior **Grip**,
Heavy Duty Protection*

Series No.:
21590OR X-Tra Thick Pro Orange T-Grip Nitrile Powder Free Disposable Gloves



BENEFITS OF X-TRA THICK T-GRIP NITRILE POWDER FREE GLOVES

Heavy duty protection:

Our ASAP X-Tra Thick Orange T-Grip Nitrile is one of our toughest ever gloves! Providing superior heavy duty protection for some of the most demanding tasks.

Secure Grip:

The T-Grip pattern provides the "best in class" grip, ensuring a secure hold and minimizing accidental slips or drops, especially in wet and oily conditions.

Versatile:

Unleash versatility with ASAP's X-Tra Thick Orange T-Grip Nitrile Gloves. Perfect for automotive, manufacturing, food processing, and more. Experience the power of gloves that effortlessly elevate your performance in any professional environment.

Colourful option:

Spark productivity with the vibrant **ORANGE** colour of ASAP's X-Tra Thick Orange T-Grip Nitrile Gloves, ensuring high visibility and empowering you to work with confidence and precision.

Application:

ASAP x-Tra Thick T-Grip Orange Powder Free Disposable Gloves are recommended to be used in industries for:

- Automotive repair shops and dealerships
- Food processing and manufacturing facilities
- Industrial manufacturing and processing plants
- Cleaning and janitorial services
- Plumbing and HVAC service providers
- Painters and decorators
- Agriculture and farming operations

Regulations

- PPE Regulation (EU) 2016/425
- Food Contact Regulation (EU) 2020/1245 of Regulation (EU) No 10/2011

Harmonised Standards

- EN ISO 21420:2020
- EN 374-1:2016+A1:2018
- EN 374-5:2016

Quality Assurance

- ISO 9001:2015
- ISO 13485:2016
- ISO 14001:2015



At ASAP, we are committed to hygiene control and quality assurance. Proper hygiene standard is practiced throughout the development of all ASAP products from raw materials handling, processing, production, inspection, to our finished product to deliver high quality products while limiting risk of cross-contamination.

Look for the Hygiene Matters™ logo, quality and hygiene you can trust.



X-TRA THICK T-GRIP

NITRILE POWDER FREE DISPOSABLE GLOVES



Color Option:



Automotive



Safety



Food & Beverage



Series No.:
21590OR X-Tra Thick Pro Orange T-Grip Nitrile Powder Free Disposable Gloves

Series Size Codes

Medium	Large	Extra Large	Extra Extra Large
M, 8	L, 9	XL, 10	XXL, 11
21593OR	21594OR	21595OR	21596OR

Product Specifications

Design	Ambidextrous, Diamond Grip Surface Beaded Cuff
Colour	Orange
Acceptance Quality Level (AQL)	1.5
Packing Mode	50 pcs per box, 10 boxes per carton

Dimension Specifications

Glove Size	Item Codes	Palm Width (mm)	Length (mm)		Thickness Single Wall (mm)		
			EN 455	ASTM	Cuff (25±5 from bead)	Palm (centre of palm)	Finger (13±3 from tip)
M / 8	21593OR	95 ± 5	Min. 240	Min. 230	0.09± 0.03	0.22 ± 0.03	0.25 ± 0.02
L / 9	21594OR	105 ± 5					
XL / 10	21595OR	115 ± 5					
XXL / 11	21596OR	125 ± 5					

Physical Properties Specifications

	EN 455 Force at Break (N)	ASTM Tensile Strength (MPa)	ASTM Elongation (%)
Before Aging	Min. 6.0,	Min. 14	Min. 500
After Aging	Min. 6.0,	Min. 14	Min. 400

Packaging Dimensions

Inner	235 x 125 x 62 mm
Carton	325 x 265 x 250 mm

Powder Residue

Powder Free (mg/glove)	Max. 2
------------------------	--------

Instructions For Use

Description - ASAP X-Tra Thick T-Grip Nitrile Powder Free Disposable Gloves, Non-sterile, Single Use Only.

Intended Use - ASAP X-Tra Thick T-Grip nitrile glove is a disposable glove product worn to protect the hand of wearer against mechanical action whose effects are superficial, cleaning materials of weak action and easily reversible effects.

How To Don Gloves - Inspect the gloves to ensure there are no pinholes or tears. If gloves are ambidextrous, they can be worn on either hand. If not, align the glove's fingers and thumb with the proper hand before donning. Insert five fingers into the cuff and pull the cuff over the wrist. Check for a secure fit around the fingers and palm. The cuff should fit snugly around the wrist.

How To Doff Gloves - After use, users should visually check the glove and remove any contamination from the outer surface before removing the gloves from the hands. Grasp the outside edge of the glove near the wrist. Peel the glove away from the hand, turning it inside out. Hold it in the opposite glove hand. Slide an ungloved finger under the wrist of the remaining glove, be careful not to touch the outside of the glove. Peel the remaining glove off from the inside, creating a "bag" containing both gloves. Discard.

Disposal - Properly dispose of all used nitrile glove. Follow your institution's policies for use and disposal of these gloves.

Storage - Store in a dry place. Avoid excessive heat (30°C). Exposed product should be shielded from direct sunlight, intense artificial light, x-ray machines, and other source of ozone.

Shelf Life - Three years from the manufacturing date.

Warning - These gloves are for single and transient use only.

Caution - This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal or over 400mm - where the cuff is also tested) and relates only to the chemical tested.

It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on the temperature, abrasion, and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in the physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly.

For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. This product contains nitrile rubber, which may cause allergic reactions in individuals who are known or suspected to be allergic to nitrile rubber. If an allergic reaction occurs, stop using immediately and consult a physician. This product is not made of natural rubber latex.

EN ISO 374

Chemical Permeation (EN ISO 374-1:2016+A1:2018/Type B)	Level	Mean Degradation % (EN ISO 374-4:2019)
K 40% Sodium Hydroxide	5	10.0
P 30% Hydrogen Peroxide	2	36.5
T 37% Formaldehyde	5	17.5

Degradation levels indicate the change in Puncture Resistance of the glove after exposure to the challenge chemical.

EN ISO 374

EN 16523-1:2015+A1:2018 Classification of Permeation Performance Level

Measured Breakthrough Time (min)	>10	> 30	> 60	> 120	> 240	> 480
Permeation Performance Level	1	2	3	4	5	6

The penetration levels have been assessed under laboratory conditions and relates only to the tested specimen.

Resistance against Bacteria and Fungi - PASS
Resistance against Virus - PASS

ASAP INTERNATIONAL SDN BHD
No. 1, Jalan Sitar 33/6, Seksyen 33,
40400 Shah Alam, Selangor,
Malaysia.

T : +603 5191 0166
F : +603 5191 0702
E : info@whyasap.com
W : www.whyasap.com

ASAP INNOVATIONS LTD.
Unit 7, The Courtyard,
Fonthill Business Park,
Fonthill Road, Dublin,
D22 XA07, Ireland.

T : +353 1466 1660
E : info@whyasap.ie
W : www.whyasap.ie

ASAP INNOVATIONS (UK) LTD.
13, Diamond Court,
Opal Drive, Fox Milne,
Milton Keynes,
MK15 ODU, United Kingdom.

T : +44 (0) 1908 732700
E : info@whyasap.co.uk
W : www.whyasap.co.uk

