

# Label-Free Cell-based Assay for Drug Discovery using the Corning Epic®

## System and Corning HYPERFlask® Cell Culture Vessel

Lori Romeo, Jeff Scibek, Arron Xu, Kim Titus and Todd Upton, Corning Life Sciences, Corning Inc., Corning, NY 14830; Ralph Garippa, Theresa Truitt and Rachid Hamid, Hoffmann-La Roche, Nutley, NJ



### Introduction

#### Abstract

The Corning Epic® technology is an optical resonant waveguide-based, label-free detection system. It provides a label-free platform of high detection sensitivity for measuring both molecular interactions in a biochemical assay, as well as integrated cellular responses in a cell-based assay. Detection of live cell and time dependent cellular response in a pathway unbiased manner with the Epic® system may provide previously unattainable biological and pharmacological information for an integrated drug-stimulated cellular response. Moreover, the Epic system is capable of performing primary screening in both biochemical and cell based assays meeting a throughput requirement of approximately 80 384-well microplates/8 hr run. To support the demand of cells for HTS cell-based assays, the Corning HYPERFlask® vessel provides a solution with high capacity and high efficiency for cell culture with each HYPERFlask® vessel providing the number of cells equivalent to that from 10 standard T-175 flasks. The study reported here describes the successful use of the Epic system for a label-free screening to discover agonists targeting serotonin receptors in a cell-based GPCR assay. It is also demonstrated that the cells cultured in the conventional T-flask and those in the HYPERFlask® vessel show highly similar assay performance, thus further supporting the use of HYPERFlask® vessel for HTS cell-based assays.

#### Background

The Corning Epic® technology provides a label-free detection system that measures integrated cellular response in a cell-based assay or molecular interaction response in a biochemical assay. It enables high throughput screening for drug discovery in a label-free format with robust signal detection and sensitivity. Cell-based high throughput screens require a large supply of cells that often present challenges for conventional cell culture techniques. Conventional cell culture requires numerous T-flasks, significant incubator space and a high demand for skilled labor, resulting in substantial costs. The Corning HYPERFlask® vessel provides a solution to the need of high capacity cell culture in support of cell-based HTS screen with its novel design. We present here a comparative study of Epic® assay performance using cells cultured in a HYPERFlask® vessel and conventional T-150 flasks. The Epic® results demonstrated highly similar performance with cells cultured in the HYPERFlask® vessel and in conventional T-flasks. Results of this study showed a significant savings in materials, time and labor costs with equivalent assay performance.

### Materials & Methods

#### Cell Seeding

- HYPERFlask vessels and T-150 flasks were seeded with  $1 \times 10^4$  cells/cm<sup>2</sup> with engineered CHO cells
- Flasks were incubated at 37°C with 5% CO<sub>2</sub>
- Cells were grown for 72 hours at which time they were ~ 90% confluent


#### Harvesting:

- Cells were harvested using Trypsin-EDTA
- Cell Yields:
  - HYPERFlask® vessel:  $2.6 \times 10^4$  cells per vessel ( $1.5 \times 10^5$  cells/cm<sup>2</sup>)
  - T-150 Flask:  $0.22 \times 10^4$  cells per flask ( $1.5 \times 10^5$  cells/cm<sup>2</sup>)


#### Plate Seeding:

- Fibronectin-coated Epic® microplates were seeded with 7,500 cells/well using a Multi-Drop Combi


### Assay Flowchart




### Results




**Figure 1: Dose-dependent response of engineered CHO cells stimulated with serotonin.** Left panels: Epic® response profiles for engineered CHO cells stimulated with a dilution series of serotonin. Serotonin concentration range: 0.3nM – 5µM. Arrow indicates the point of compound addition. Right panel: Dose-response curves and EC<sub>50</sub> values from Epic response profiles.



**Figure 2: Representative cell-based screening data for engineered CHO cells cultured on a T-150 flask or HYPERFlask® vessel.** Compounds were screened at 5µM. Left panel: Scatter plot for cells from T-150 flasks. Right panel: Scatter plot for HYPERFlask® vessels. Red circle: Negative control (assay buffer, N=16) Green circle: Positive control (Serotonin, EC<sub>50</sub> concentration, N=16). The threshold was set at 40% of the response of the positive control. The same 5 compounds were identified as hits on both platforms. Z' > 0.7



**Figure 3: Agonist dose-response curves for a subset of the compounds that were identified as hits during the cell-based screen.** A dilution series of each compound was added to the cells and the response was measured at 20 minutes. Concentration range: 0.2nM – 12.5µM



**Figure 4: Inhibition curves for a subset of the compounds that were identified as agonist hits during the cell-based screen.** A dilution series of each compound was added to the cells and the cells were incubated for ~30 minutes. An EC<sub>50</sub> concentration of serotonin was added to the cells and the response was measured after an additional 20 minutes. Compound #5 inhibits the serotonin response most likely due to receptor desensitization. Compounds #2 and #6 do not inhibit the serotonin response. These compounds may be targeting other endogenously expressed receptors on the engineered CHO cells (not the serotonin receptor).

**Table 1: Summary of EC<sub>50</sub> and IC<sub>50</sub> values measured on the Epic® System**

Compound	EC <sub>50</sub> (nM)		IC <sub>50</sub> (nM)	
	T-150	HYPERFlask® Vessel	T-150	HYPERFlask® Vessel
1	5	10	1	5
2	539	578	3	No inhibition
3	1006	964	2	No inhibition
4	126	76	4	No inhibition
5	19	16	5	20
6	78	77	6	No inhibition
7	>10000	>10000	7	No inhibition
8	436	480	8	382
9	690	1008	9	No inhibition
10	776	1154	10	No inhibition
11	352	417	11	146
12	>10000	>10000	12	2491
13	233	309	13	No inhibition
14	54	57	14	30
15	375	570	15	No inhibition
16	585	530	16	No inhibition

- EC<sub>50</sub> and IC<sub>50</sub> values were calculated using GraphPad Prism.
- Similar EC<sub>50</sub> and IC<sub>50</sub> values were observed for the all compounds tested.
- Dose-response curves for compounds 2, 5, and 6 are shown in Figures 3 and 4.

### Summary and Conclusions:

- Cells cultured on the HYPERFlask® vessel and a T-150 flask exhibited similar morphology and Epic assay performance.
- Epic response profiles, EC<sub>50</sub> values and assay robustness (Z'>0.7) for cells from both vessels are in good agreement.
- The same compounds were identified as hits during the agonist cell-based screen with both sets of cells.
- The HYPERFlask® vessel can be used to replace multiple T-flasks for cell-based high throughput screening with no compromise to the results of the assay.
- Using the HYPERFlask® vessel will:
  - Increase cell yield
  - Generate less waste
  - Decrease labor costs
  - Require less incubator space

- Results demonstrated that Epic can be used as a robust label-free screening tool for high throughput drug discovery.
- The combination of the HYPERFlask® vessel and Epic provides a highly efficient and robust platform for high throughput cell-based screening in a label-free format.