

Determination of Apical to Basolateral Permeability for 21-Day Caco-2 Cell Monolayers for Falcon® 96 Multiwell Insert Systems

CORNING

Protocol

Purpose

This is a method to determine the apical to basolateral permeability of [¹⁴C]-mannitol in cell monolayers grown on Falcon 96 Multiwell Inserts using a manual (non-robotic) procedure.

Materials

- ▶ Falcon 96 Multiwell Insert System [1.0 μm, 0.0804 cm² membrane surface area] (Cat. Nos. 351130 and 351131 includes feeder tray, Cat. No. 353938 includes receiver plate)
- ▶ Receiver Plates:
 - Falcon 96 Square Well, Angled-Bottom Plates (Cat. No. 353925)
- ▶ Wash Buffer:
 - Hanks Balanced Salt Solution (HBSS) +10 mM HEPES, pH 7.4
- ▶ Receiver Solution:
 - HBSS +10 mM HEPES +0.5% DMSO
- ▶ Donor Solution:
 - 50 μM [¹⁴C or ³H]-Mannitol in HBSS +10 mM HEPES +0.5% DMSO

Example:

- ▶ 5 μL of [¹⁴C]-Mannitol stock (Amersham Cat. No. CFA238, 3.4 mM, 0.2 μCi/μL)
- ▶ 44.9 μL of 10 mM Mannitol (non-radioactive) in DMSO
- ▶ 5.1 μL of DMSO
- ▶ 9935 μL HBSS +10 mM HEPES

Procedure

Washing the Monolayers

1. Transfer the insert plate from the feeder tray or 96 well receiver plate to an empty 96 well receiver plate.
2. Wash the monolayers, one column at a time:
 - Using a yellow tip attached to a vacuum, aspirate the culture medium from the apical side of a column of monolayers.
 - Using a multichannel or repeating pipettor, add back 50 μL of wash buffer to the apical sides of the monolayers.
3. Place the insert plate into a feeder tray containing 30 mL of wash buffer, or a 96 well receiver plate containing 260 μL of wash buffer/well. Keep plate at 37°C.

Permeability Assay

1. Transfer the insert plate from the wash plate to an empty 96 well receiver plate.
2. Remove the apical wash buffer and add the donor solution, one column at a time:
 - Using a yellow tip attached to a vacuum, aspirate the wash buffer from the apical side of a column of monolayers.
 - Using a multichannel or repeating pipettor, deliver 50 μL of donor solution to the apical sides of the monolayers.

- Place the insert plate into a 96 well receiver plate containing 260 μL of receiver solution/well.
- Incubate the plate at 37°C, with orbital shaking at 50-100 rpm for 2 hours.

After 2 hours of incubation:

- Remove the insert plate from the receiver tray, place it into an empty 96 well receiver plate, and transfer 40 μL of the solution in the apical chamber (donor at t = 120 min.) to scintillation vials. Count.
- Transfer 100 μL of the solution in the basolateral chamber (receiver at t = 120 min.) to scintillation vials. Count.
- Count 50 μL of the original donor solution (t = 0 min.). This represents the amount of Mannitol added to the monolayers (donor). Determine pmol/dpm ratio.

Example:

52699 dpms counted in 50 μL of 50 μM (50 pmol/ μL) mannitol donor (50 μL) x (50 pmol/ μL) = 2500 pmol added to donor side 2500 pmol / 52699 dpms = 0.0474 pmol/dpm

- Calculate amount of pmol in receiver tray. Dpms in 100 μL of receiver chamber solution x volume factor (2.6) = total dpms in chamber x pmol / dpm ratio = pmol in receiver chamber
- Calculate Mannitol flux from donor side (apical) to receiver side (basolateral).

Flux = pmol in receiver chamber at 120 min. / 2500 pmol in donor chamber at t = 0-10.
Calculate Mannitol P_{app} at 120 min. (cm/sec).

$P_{app} = (\text{Flux of drug}) \times (\text{vol in donor chamber}) / (\text{sec of incubation}) \times (\text{surface area})$

$P_{app} = (\text{Flux}) \times (0.05 \text{ cm}^3 / [7200 \text{ sec} \times 0.0804 \text{ cm}^2])$

- Calculate Mannitol mass balance (donor and receiver chambers).

Mass Balance = $([\text{pmol in 120 min. receiver} + \text{pmol in 120 min. donor}] / 2500 \text{ pmol}) \times 100$

Corning acquired the Falcon® brand.

For additional Corning product, technical, or distributor information, please e-mail us at CLSTechServ@corning.com, visit our website www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call 978.442.2200.

For information on the acquisition, visit www.corning.com/discoverylabware.

**Corning Incorporated
Life Sciences**

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

**Worldwide
Support Offices**

ASIA/PACIFIC

Australia/New Zealand
t 0402-794-347

China
t 86 21 2215 2888
f 86 21 6215 2988

India
t 91 124 4604000
f 91 124 4604099

Japan
t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE
France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

**All Other European
Countries**
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA
Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589

CORNING | **FALCON®** **cellgro®** **AXYGEN®** **PYREX®** **GOSELIN™**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks.
Corning Incorporated, One Riverfront Plaza, Corning, NY 14831-0001