

# P150 Adaptation of Corning® HepatoCells for 3D Spheroid Culture and Hepatotoxicity Studies

CORNING

Feng Li<sup>1</sup>, Rongjun Zuo<sup>1</sup>, Mercyanne Andes<sup>1</sup>, Kirsten Cooper<sup>1</sup>, Sweta Parikh<sup>1</sup>, Sarah O'Donnell<sup>3</sup>, Ronald A. Faris<sup>2</sup>, Alice Gao<sup>4</sup>, Charles L. Crespi<sup>4</sup>, and Christopher J. Patten<sup>1</sup>


<sup>1</sup>Two Oak Park, Corning Incorporated, Life Sciences, Bedford, MA 01730 USA; <sup>2</sup>1 Science Center Drive, Corning Incorporated, Life Sciences, Corning, NY 14831 USA;

<sup>3</sup>6 Henshaw Street, Corning Incorporated, Life Sciences, Woburn, MA 01801 USA; <sup>4</sup>836 North Street, Bldg 300 Ste 3401, Corning Incorporated, Life Sciences, Tewksbury, MA 01876 USA

## Introduction


- In vitro* hepatic 3D cultures are known to better reflect the *in vivo* behaviors of liver cells; Industry is paying increasing attention to 3D *in vitro* systems.
- Corning® HepatoCells (Cat. No. 354881) are a renewable hepatic cell source derived from primary human hepatocytes with consistent phenotypes, drug metabolism enzyme/transporter activities.
- Adaptation of HepatoCells to spheroid culture (96-well or 384-well format) and tests with classic liver toxins demonstrate the applications of HepatoCells for ADME research and liver cytotoxicity studies.

## Correlation between Sizes of Spheroids and Initial Seeding Density



Corning HepatoCells were seeded at indicated cell densities/well in 384-well Corning Ultra-Low Attachment surface spheroid microplate (Cat. No. 4516). No medium changes were made for the first two days and spheroid cultures were then maintained with half medium change every other day. Bar = 100 µm.

## Higher CYP3A4 Activity in Corning HepatoCells Spheroid Microplate vs. 2D Culture (Normalized to ATP Levels)



Corning HepatoCells were seeded at 1.5K/well for spheroid culture in 96-well Corning Ultra-Low Attachment surface microplates (Cat. No. 4515). 100K/well HepatoCells were used to set up 2D culture in 96-well Corning BioCoat™ Collagen I microplate. CYP3A4 activity was measured with Promega Luciferin-IPA kit for day 7 spheroids (pool of 6 spheroids = 9K cells or 2D (100K cells/well) culture (chart, left). ATP levels (chart, middle) were measured for spheroids and 2D cultures and used for normalizing the CYP3A4 activity (chart, right). Error bars = Standard Deviation, n = 3.

## Inducible CYP3A4 Activity in Corning HepatoCells Spheroids

CYP3A4 Activity Fold Change*	Single Spheroid		5 spheroids	
	DMSO	RIF	DMSO	RIF
	167.1 (n=5)	535.7 (n=4)	983.9	3483.8
		3.2		3.5

CYP3A4 3-day induction was performed with 10 µM Rifampicin (RIF) for day 4, 5 and 6 for single or pooled 5 spheroids. 0.1% DMSO was included as negative control. CYP3A4 activity was measured on day 7 with Luciferin-IPA probe.

\* Corning HepatoCells 2D culture shows a CYP3A4 induction fold change of 19.6 (Poster 115 by Rongjun Zuo, et al).


\* CYP3A4 induction fold change in 3D culture (e.g., primary hepatocytes) is known to be lower than 2D culture due to increased basal enzyme activity.

## Conclusion


- Corning HepatoCells were successfully adapted to a high throughput 3D culture system using Corning Ultra-Low Attachment surface spheroid microplates.
- Tests with classic liver toxins such as Aflatoxin B1 and Acetaminophen, demonstrate that HepatoCells spheroids are a promising 3D model system for liver cytotoxicity studies.

## Uniform Corning HepatoCells Spheroid Formation on a Corning Ultra-Low Attachment Surface Spheroid Microplate


### Hematoxylin and Eosin Staining Shows Healthy Core Structure of a Corning HepatoCells Spheroid




### ATP Assay for Day 7 Spheroids



### Spheroid Size Measurement Over a 7-day Culture




### Close Correlation Between ATP Levels and Cell Numbers at Lower Seeding Density




Corning HepatoCells were seeded at indicated cell densities/well in a 384-well Corning Ultra-Low Attachment surface spheroid microplate. For H&E staining, day 7 spheroids from HepatoCells (1.5K/spheroid) were collected and fixed in 10% PBS buffered formalin for >24 hours prior to paraffin embedding and 10 µm sections were prepared and stained. Spheroid images were taken and sizes (diameter, µm) were measured during the culture period. ATP Levels (viability) were measured with a Promega CellTiter-Glo® 3D kit.

## Corning HepatoCells Spheroid Treatment with Classic Liver Toxins

### Day 4 Spheroid Culture vs. 2D Culture to Aflatoxin B1 (AFB1) Treatment



## EC<sub>50</sub> Values of Day 7 Corning HepatoCells Spheroid to Classic Liver Toxins and Comparisons to Literature



Toxin	EC <sub>50</sub> (µM) HepatoCells, 3D	Reported EC <sub>50</sub> (µM) (2D culture)	References
Aflatoxin B1	5.3	10 (HH, 2D) >100 (HepG2, 2D) 78.4 (HepatoCells, 2D)	Avior Y. et al, 2015 In-house data
Acetaminophen	14610	26200 (HH), 29755 (HepG2)	Jermolitz et al, 2008 Wang et al, 2002
Menthol	Non toxic	Non toxic	Avior Y. et al, 2015

### HH, human hepatocytes

Day 7 spheroid cultures (1.5K/spheroid) in 96-well format were treated for 24 hours with classic liver toxins and non-toxic control for 24 hours. ATP levels were measured at the end of incubation for cell viability. Non-linear regression and EC<sub>50</sub> values were calculated using Graphpad. Error bars = Standard Deviation, n = 3.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Use of Genetically Modified Microorganisms (GMM): Information for European Customers: These products are genetically modified as described in Corning Life Sciences technical literature. As a condition of sale, use of this product must be in accordance with all applicable local guidelines on the contained use of genetically modified microorganisms, including the Directive 2009/41/EC of the European Parliament and of the Council. © 2015 Corning Incorporated. All rights reserved.