



Short instructions Eppendorf Easypet® 3 English (EN)

These short instructions are a supplement to the operating manual and not a substitute. Please read the operating manual before using the Easypet 3 for the first time!

Battery status display




✓ ready to use



Indicates low battery
➔ connect the charger and continue working

Charging process




Charging battery


Light goes off
when fully charged

Short instructions
Eppendorf Easypet® 3
English (EN)

Speed control




Slow




Fast

Wall mounting device



Shelf stand (not available in the USA)



Your local distributor: www.eppendorf.com/contact
Eppendorf AG · 22331 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com