

Centrifuge 5702/5702 R/5702 RH, MiniSpin/MiniSpin plus

This document supplements the operation manual Centrifuge 5702, 5702 R, 5702 RH, MiniSpin, MiniSpin plus and does not replace it. Therefore, please also read the operating manual before starting up the rotors for the first time. You can find the current version on the Internet at www.eppendorf.com/worldwide.

Transport and storage

Transport

- ▶ Use the original packaging for transport.

	Air temperature	Relative humidity	Atmospheric pressure
General transport	-25 °C – 60 °C	10 % – 75 %	30 kPa – 106 kPa
Air freight	-20 °C – 55 °C	10 % – 75 %	30 kPa – 106 kPa

Storage

	Air temperature	Relative humidity	Atmospheric pressure
In transport packaging	-25 °C – 55 °C	10 % – 75 %	70 kPa – 106 kPa
Without transport packaging	-5 °C – 45 °C	10 % – 75 %	70 kPa – 106 kPa

Interference emission and noise immunity

Tab. 1: Centrifuge 5702

EMC: interference emission (radio interference)	EN 55011/B (230 V), EN 61000-3-2/A14 (230 V), EN 61000-3-3 (230 V), FCC 15 Part B (120 V)
EMC: noise immunity	EN 50082-1 (230 V), EN 61000-4-14 (230 V)

Tab. 2: Centrifuge 5702 R/RH

EMC: interference emission (radio interference)	EN 55011/B (230 V), EN 61000-3-2 (230 V), EN 61000-3-3 (230 V), FCC 15 Part B (120 V)
EMC: noise immunity	EN 61000-6-1 (230 V), EN 61000-4-14 (230 V)

Tab. 3: MiniSpin/MiniSpin plus

EMC: interference emission (radio interference)	EN 61326-1 Klasse B (230 V), EN 61000-3-2/A14 (230 V), EN 61000-3-3 (230 V), CFR FCC 15 subpart B class B (120 V)
EMC: noise immunity	EN 50082-1 (230 V)

eppendorf

eppendorf® is a registered trademark
01 16 900.017-00/112012

Your local distributor: www.eppendorf.com/worldwide

Eppendorf AG · 22331 Hamburg · Germany · Tel: +49 40 53801-0 · Fax: +49 40 538 01-556 · E-mail: eppendorf@eppendorf.com
Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA
Tel: +1 516 334 7500 · Toll free phone: +1 800-645-3050 · Fax: +1 516 334 7506 · E-mail: info@eppendorf.com

Application Support Europe: Tel: +49 1803 666 789 (Preis je nach Tarif im Ausland; 9 ct/min aus dem dt. Festnetz; Mobilfunkhöchstpreis 42 ct/min)
support@eppendorf.com

North America: Tel: +1 800 645 3050 · E-mail: techserv@eppendorf.com
Asia Pacific: Tel: +60 3 8023 6869 · E-mail: support_asiapacific@eppendorf.com