

# Centrifuge 5702/5702 R/5702 RH, MiniSpin/MiniSpin plus

This document supplements the operation manual Centrifuge 5702, 5702 R, 5702 RH, MiniSpin, MiniSpin plus and does not replace it. Therefore, please also read the operating manual before starting up the rotors for the first time. You can find the current version on the Internet at [www.eppendorf.com/worldwide](http://www.eppendorf.com/worldwide).

## Transport and storage

### Transport

- ▶ Use the original packaging for transport.

	Air temperature	Relative humidity	Atmospheric pressure
General transport	-25 °C – 60 °C	10 % – 75 %	30 kPa – 106 kPa
Air freight	-20 °C – 55 °C	10 % – 75 %	30 kPa – 106 kPa

### Storage

	Air temperature	Relative humidity	Atmospheric pressure
In transport packaging	-25 °C – 55 °C	10 % – 75 %	70 kPa – 106 kPa
Without transport packaging	-5 °C – 45 °C	10 % – 75 %	70 kPa – 106 kPa

## Interference emission and noise immunity

Tab. 1: Centrifuge 5702

EMC: interference emission (radio interference)	EN 55011/B (230 V), EN 61000-3-2/A14 (230 V), EN 61000-3-3 (230 V), FCC 15 Part B (120 V)
EMC: noise immunity	EN 50082-1 (230 V), EN 61000-4-14 (230 V)

Tab. 2: Centrifuge 5702 R/RH

EMC: interference emission (radio interference)	EN 55011/B (230 V), EN 61000-3-2 (230 V), EN 61000-3-3 (230 V), FCC 15 Part B (120 V)
EMC: noise immunity	EN 61000-6-1 (230 V), EN 61000-4-14 (230 V)

Tab. 3: MiniSpin/MiniSpin plus

EMC: interference emission (radio interference)	EN 61326-1 Klasse B (230 V), EN 61000-3-2/A14 (230 V), EN 61000-3-3 (230 V), CFR FCC 15 subpart B class B (120 V)
EMC: noise immunity	EN 50082-1 (230 V)

# eppendorf

eppendorf® is a registered trademark  
01 16 900.017-00/112012

Your local distributor: [www.eppendorf.com/worldwide](http://www.eppendorf.com/worldwide)

Eppendorf AG · 22331 Hamburg · Germany · Tel: +49 40 53801-0 · Fax: +49 40 538 01-556 · E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)  
Eppendorf North America, Inc. · 102 Motor Parkway · Hauppauge, N.Y. 11788-5178 · USA  
Tel: +1 516 334 7500 · Toll free phone: +1 800-645-3050 · Fax: +1 516 334 7506 · E-mail: [info@eppendorf.com](mailto:info@eppendorf.com)

**Application Support** Europe: Tel: +49 1803 666 789 (Preis je nach Tarif im Ausland; 9 ct/min aus dem dt. Festnetz; Mobilfunkhöchstpreis 42 ct/min)  
[support@eppendorf.com](mailto:support@eppendorf.com)

North America: Tel: +1 800 645 3050 · E-mail: [techserv@eppendorf.com](mailto:techserv@eppendorf.com)  
Asia Pacific: Tel: +60 3 8023 6869 · E-mail: [support\\_asiapacific@eppendorf.com](mailto:support_asiapacific@eppendorf.com)