

Corning® Repeater and Syringe Tips

CORNING

Corning Repeater

The Corning Repeater is suitable to dispense volatile, dense, viscous, or foaming liquids, and for all applications where long dosing series are required. It works with 8 syringe tip sizes to dispense doses from 1 μ L up to 5 mL. The lightweight construction, ergonomically shaped handle with a finger hook allows for single-handed volume setting, fast syringe tip filling, and convenient serial dispensing. The Repeater's precise stepping mechanism ensures repeatability of all dispensed doses, and the positive displacement allows for contamination-free operation.

Each device is tested for functionality and precision according to EN ISO 8655 and is supplied with an individual QC certificate.

Features

- ▶ Easy volume setting: Volume selection chart conveniently located on the device
- ▶ Fast and precise multiple dispensing: Up to 48 pipetting doses without refilling the tip syringe
- ▶ Wide range of polypropylene syringe tips: 8 different sizes, sterile and non-sterile versions available
- ▶ Positive displacement syringe tips: Dedicated to work with high density liquids or vapor pressure such as glycerol or ethanol
- ▶ Maintenance-free

Corning Repeater Starter Pack

In order to enable a smooth start, the Corning Repeater Starter Pack is conveniently packaged with a Corning Repeater, 2 syringe tip adapters (25 mL and 50 mL), and 24 syringe tips (0.05 mL, 2.5 mL, 5 mL, 12.5 mL, 25 mL, and 50 mL).



Corning Repeater is supplied with two syringe tip adapters (25 mL and 50 mL).



Corning Repeater Starter Pack

Ordering Information

Corning® Repeater and Syringe Tips

Cat. No.	Description	Qty/Cs
4217	Corning Repeater	1
6620	Syringe tips, 0.05 mL	100
6628	Syringe tips, 0.05 mL, sterile	100
6621	Syringe tips, 0.5 mL	100
6629	Syringe tips, 0.5 mL, sterile	100
6622	Syringe tips, 1.25 mL	100
6630	Syringe tips, 1.25 mL, sterile	100
6623	Syringe tips, 2.5 mL	100
6631	Syringe tips, 2.5 mL, sterile	100
6624	Syringe tips, 5 mL	100
6632	Syringe tips, 5 mL, sterile	100
6625	Syringe tips, 12.5 mL	100
6633	Syringe tips, 12.5 mL, sterile	100
6626	Syringe tips, 25 mL	25
6634	Syringe tips, 25 mL, sterile	25
6627	Syringe tips, 50 mL	25
6635	Syringe tips, 50 mL, sterile	25

Accessories

4218	Adapter for 25 mL syringe tips	1
4219	Adapter for 50 mL syringe tips	1

Corning Repeater Starter Pack

Cat. No.	Description	Qty/Cs
6636	Corning Repeater Starter Pack Includes: Corning Repeater (Cat. No. 4217); Adapter for 25 mL syringe tips (Cat. No. 4218) and 50 mL syringe tips (Cat. No. 4219); Syringe tips (5/pk) of 0.05 mL (Cat. No. 6620), 2.5 mL (Cat. No. 6623), 5 mL (Cat. No. 6624), 12.5 mL (Cat. No. 6625); Syringe tips (2/pk) of 25 mL (Cat. No. 6626) and 50 mL (Cat. No. 6627)	1 kit

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

ASIA/PACIFIC

Australia/New Zealand
t 61 427286832

China
t 86 21 3338 4338
f 86 21 3338 4300

India
t 91 124 4604000
f 91 124 4604099

Japan

t 81 3-3586 1996
f 81 3-3586 1291

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6572-9740
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2516-7500

EUROPE

CSEurope@corning.com

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 020 655 79 28
f 020 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117

All Other European Countries

t +31 (0) 206 59 60 51
f +31 (0) 206 59 76 73

LATIN AMERICA
grupoLA@corning.com

Brasil
t 55 (11) 3089-7400

Mexico
t (52-81) 8158-8400