

Why it is so important to "get the lead out.".....

Standardizing on lead-free products, such as 3M[™] Comply[™] Steam Indicator Tapes, makes good sense – by helping to simplify operations, avoid the risk of penalties for non-compliance, and enhancing your reputation for good environmental stewardship.

3M adhesives are engineered to fit our customers' needs and are used in such diverse products as airplanes and cell phones. The industry-standard product for exposure monitoring now has even more to offer.

- Does not require hazardous waste disposal
 - Redesigned packaging uses up to 30% less plastic
 - Convenient, individually wrapped, easy-to-open rolls
 - Sustainable, solvent-free manufacturing process





Best Practices

Always place an AAMI Class 1 Process Indicator* on the outside of every pack, wrapped tray, container or peel pouch, unless the internal chemical indicator is readily visible.

Unexposed



Unexposed



Exposed



The chemical indicator ink, printed in a diagonal striped pattern, changes color or turns darker when exposed to steam.

*As defined by ANSI/AAMI/ISO 11140-1:2005 (R)2010

Ordering Information

Indicator Tapes

These Comply Indicator Tapes (ANSI/AAMI/ISO Class 1 Process Indicators) show at-a-glance whether packs have been exposed to a steam sterilization process and function as a closure to securely seal the sterilization packs. They accept hand-writing and labeling for recordkeeping information.

Catalog Number	Product Name	Size	Rolls/ Intermediate	Rolls/Case
1322-12MM	3M™ Comply™ Lead Free Steam Indicator Tape (designed for use on all wraps)	0.47 in x 60 yd (12mm x 55m)	_	42
1322-18MM		0.70 in x 60 yd (18mm x 55m)	_	28
1322-24MM		0.94 in x 60 yd (24mm x 55m)	_	20
1322-48MM		1.89 in x 60 yd (48mm x 55m)	_	10
1355-18MM	3M™ Comply™ Lead Free Steam	0.70 in x 60 yd (18mm x 55m)	_	28
1355-24MM	Indicator Tape (designed for use on disposable wraps)	0.94 in x 60 yd (24mm x 55m)	-	20

Indicator Tape Dispenser

Catalog Number	Product Name	Description	Items/Box	Boxes/Case
C22	Heavy Duty Tape Dispenser	Holds 2 rolls	1	2
M52	Heavy Duty Tape Dispenser	With Tabber	1	1

Sustainable Solutions for a Healthier World

Environmental Stewardship

In partnership with communities around the globe, we continually strive to protect our planet's health by developing sustainable solutions, business practices and manufacturing processes. Protecting the environment is central to the mission of 3M.

Economic Success

By creating viable business opportunities, investing in local economies and supporting growth and development in emerging markets, we're always looking for innovative ways to enrich the vitality of local communities around the world.

Social Responsibility

We are committed to fostering vibrant communities by actively contributing our expertise, resources and financial support on both corporate and individual levels.



Infection Prevention Division

3M Health Care 2510 Conway Avenue St. Paul, MN 55144-1000 IISΔ

1-800-228-3957 www.3M.com/infectionprevention **Medical Division** 3M Canada Post Office Box 5757 London, Ontario N6A 4TI Canada

1 800 563-2921

Please recycle. Printed in USA. © 3M 2013. All rights reserved. Issued 10/13. 9438HB 70-2011-5436-9

3M and Comply are trademarks of 3M. Used under license by 3M subsidiaries and affiliates.

