

# STRIPSPIN™ 12D

## Digital Mini Centrifuge



- + Unique rotor accepts 12-position PCR strips
- + Digital control up to 5,500 rpm
- + “Mini” sized, less than 8 inches wide
- + Timer up to 99 mins.



Holds up to 4 x PCR strips or 48 x 0.2ml tubes

The StripSpin 12D from Benchmark is the only mini centrifuge designed for the spin down of 12-position PCR strips. The centrifuge features digital speed control up to 5,500 rpm, time control up to 99 min., and a safety interlock to stop the centrifuge if the lid is opened during operation.

With built in rcf conversion, speed can be set/displayed in rpm or rcf. The StripSpin 12D is the perfect solution for centrifugation of both 8 and 12-position PCR strips (or 48 x 0.2ml tubes) and includes a 2 year warranty.

### Technical Data

Speed:	5,500rpm
Max. g-force:	1,690 xg
Capacity:	4 x 12 position PCR strips 4 x 8 position PCR strips 48 x 0.2ml PCR tubes
Dimensions: (WxDxH)	7.6 x 9 x 4.7 in 19.4 x 23 x 12 cm
Weight:	4.4 lbs / 2 kg
Electrical:	100 to 240V, 50-60Hz
Warranty:	2 Years

### Also available...



myFUGE 12 Plus

Item: C2212



Includes COMBI-Rotor™ for tubes and strips

### Ordering Information

C2248\* *StripSpin* Mini Centrifuge, with rotor for up to 4 x PCR strips (8 or 12 position) or 48 x 0.2ml tubes

\*Includes US plug. For EU plug, add (-E) to the item number.

**Benchmark** Scientific

PO Box 709 Edison, NJ 08818  
Phone: 908-769-5555 Fax: 908-222-1864  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)