



## Labcold Flake Ice Makers



LITV-IQ50/IQ85C

**Labcold** Flake Ice Makers produce crushed ice ideal for filling packs and moulding to fit any item or object. The ice produced is both hygienic and clear so ideal for medical and laboratory applications.

Robustly constructed from stainless steel with well-insulated plastic collection bins and a lime scale reducing spray system.

**Labcold** Ice Makers require both a water and electricity supply and drainage 150mm below the level of the ice maker.



LITV-IQ135

| Part Number                 | LITV-IQ50                            | LITV-IQ85C                           | LITV-IQ135         | LITV-IQ200                                |
|-----------------------------|--------------------------------------|--------------------------------------|--------------------|---|
| Ice Storage                 | 20kg                                 | 20kg                                 | 40kg               | 80kg                                      |
| Daily Ice Production        | 66kg                                 | 90kg                                 | 156kg              | 234kg                                     |
| Style                       | Self Contained                       | Self Contained                       | Modular            | Modular                                   |
| Exterior Dimensions (HxWxD) | 855 (790 without legs) x 465 x 595mm | 855 (790 without legs) x 465 x 595mm | 1355 x 515 x 550mm | 575 x 515 x 550mm<br>Silo 870 x 555 x 685 |
| Weight                      | 58kg                                 | 58kg                                 | 75kg               | Ice maker 52kg<br>Silo 50kg               |
| Exterior Construction       | Stainless Steel                      | Stainless Steel                      | Stainless Steel    | Stainless Steel                           |
| Automatic Production        | ✓                                    | ✓                                    | ✓                  | ✓   |
| Feet                        | 4                                    | 4                                    | N/A                | N/A                                       |