





SpillTrays & Liners



Available in general purpose, biohazard and radiation safety formats, spilltrays provide a re-usable work area with the added benefit of safe containment of spillages and lab bench protection. All trays supplied with one free liner when purchased. Additional environmentally friendly APET liners provided in packs of 25.



spilltrays
provide a re-usable work area which provides both containment of spillages and lab bench protection

- SAVES ON RUNNING COSTS
- CONTAINS LIQUID SPILLAGES
- PROVIDES CLEARLY DEFINED WORK AREA
- EASILY CLEANED
- APET[†] - ENVIRONMENTALLY FRIENDLY WHEN INCINERATED
- SUITABLE FOR MOST HAZARDOUS SPILLAGES
- ANTI-STATIC AND NON-POROUS

WWW.CLEAVERSCIENTIFIC.COM

ORDERING INFORMATION

Tray Size	Radiation Hazard Tray, Yellow	BioHazard Tray, White	General Purpose Tray, White	APET liners, pk 25
46 x 26cm	CSR-TY1	CSR-TO1	CSR-TW1	CSR-TL1
54 x 34cm	CSR-TY2	CSR-TO2	CSR-TW2	CSR-TL2
57 x 54cm	CSR-TY3	CSR-TO3	CSR-TW3	CSR-TL3
68 x 54cm	CSR-TY4	CSR-TO4	CSR-TW4	CSR-TL4
70 x 46cm	CSR-TY5	CSR-TO5	CSR-TW5	CSR-TL5
113 x 54cm	CSR-TY6	CSR-TO6	CSR-TW6	CSR-TL6

[†] polyethylene terephthalate

CLEAVER SCIENTIFIC LTD

Unit 41, Somers Road Industrial Estate,
Rugby, CV22 7DH
United Kingdom

T_ +44 (0)1788 565300
E_ INFO@CLEAVERSCIENTIFIC.COM
W_ WWW.CLEAVERSCIENTIFIC.COM