

# TL23 Series Benchtop Turbidimeter

The new TL23 Series laboratory turbidimeters blend trusted technology of the 2100N/AN and improved features to simplify testing in the most demanding industrial and wastewater applications.

## What is the same?

The TL23 series uses the exact same optic block as the legacy 2100N/AN series turbidimeters.

No need for additional validation or evaluation

## What is new?

- Improved and intuitive design: large full colour display with guided procedures for calibration and verification
- A smart device for more reliable measurements: the instrument detects sample stability prior to taking a reading
- Easy to use. Easy to be right: with a USB port for data export, sample identification for traceability and self-diagnosis for quality assurance



## What is in the box:

TL2300 Turbidimeter, silicone oil, oiling cloth, USEPA filter assembly, 1-inch sample cells (30 mL) with caps (x6), Gelex secondary turbidity standardisation kit, StablCal calibration kit, power supply, power cord, dust cover.

For more information, please visit: [uk.hach.com/TL23](http://uk.hach.com/TL23)