

Cat. No.	10-090	10-092	10-103	15-090	16-405	90-090	90-091
Description	Liquid, 1x	Liquid, 1x	Liquid, 1x	Liquid 1X	Liquid, 1x	Powder	Powder
Units	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L	mg/L

<i>Other</i>							
Dextrose	3,151	3,151	3,151	3,151	3,151	3,151	3,151
Hypoxanthine, Sodium	2.385	2.385	2.385	2.385	2.385	2.385	2.385
DL-Thioctic (lipoic) acid	0.105	0.105	0.105	0.105	0.105	0.105	0.105
Methyl Lineoleate	0.044	0.044	0.044	0.044	0.044	0.044	0.044
Phenol Red, Sodium	8.1	8.1	8.1	8.1	--	--	8.1
Putrescine*2HCl	0.08	0.08	0.08	0.08	0.08	0.08	0.08
Sodium Pyruvate	110	110	110	110	110	110	110
Thymidine	0.365	0.365	0.365	0.365	0.365	0.365	0.365
HEPES	--	3574.8	--	--	--	--	--

Add							
NaHCO ₃ Powder (g/L)	--	--	--	--	--	2.438	2.438
7.5% Solution (mL/L)	--	--	--	--	--	32.50	32.50
L-glutamine Powder (mg/L)	--	--	--	--	--	365.10	365.10
200 mM Solution (mL/L)	--	--	--	--	--	12.50	12.50

"--" indicates not present

For more specific information on claims, visit the Certificates page at www.corning.com/lifesciences.

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Not for use in humans. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit www.corning.com/lifesciences or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

For a listing of trademarks, visit www.corning.com/clstrademarks. All other trademarks are the property of their respective owners.

CORNING

Corning Incorporated
Life Sciences

836 North St.
Building 300, Suite 3401
Tewksbury, MA 01876
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences