

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Aniline - ReagentPlus*, 99%

Product Number: 132934 CAS Number: 62-53-3

MDL: MFCD00007629

Formula: C6H7N Formula Weight: 93.13 g/mol



TEST Specification

Appearance (Color)

Material darkens in storage with no loss in purity.

Appearance (Form)
Infrared Spectrum
Purity (GC)

Specification: PRD.4.ZQ5.10000004245

Colorless to Dark Yellow and Colorless to Dark Brown-Yellow

Liquid

Conforms to Structure

> 98.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.