

Specification

DNA/RNA-ExitusPlus™ IF

A7409

Product Code:	A7409
Product Name:	DNA/RNA-ExitusPlus™ IF
Specifications:	<ul style="list-style-type: none"> : • removes DNA and RNA contaminations from surfaces : • for cleaning of PCR work stations and equipment : • for cleaning of electrophoresis equipment, pipetts, reaction tubes etc. : • delivered as spray bottles or as refill bottles (RF) : ♦ The special features : • indicator-free version of DNA-ExitusPlus™ : • All components of DNA-ExitusPlus™ IF are readily biologically degradable and not harmful or toxic for humans. : • Doesn't contain aggressive mineralic acids or alkaline substances : • No toxic fumes. Contains low concentration of alcohol. <p>Functionality (Strand breakage test): passes test</p> <p>Efficacy (ExitusPlus™ activity test): passes test</p>
WGK:	1
Storage:	RT protected from light
CS:	38220000

AppliChem GmbH