

Specification

Autoclave-ExitusPlus™

A7600

Product Code:	A7600
Product Name:	Autoclave-ExitusPlus™
Headline Comment:	<ul style="list-style-type: none"> • Autoclave-ExitusPlus™ is a registered trademark of AppliChem GmbH.
Specifications:	<p>: ♦ Additive for autoclaving</p> <ul style="list-style-type: none"> : • removes DNA and RNA contaminations from culture media waste : • for cleaning of culture vessels (e.g. Mini and Maxi prep) : • bottled <i>ready-to-use</i> : • delivered as powder mixture <p>: ♦ special, new features</p> <ul style="list-style-type: none"> : • nonenzymatic degradation of DNA and RNA by catalytic and cooperative effects of the components : • All components of Autoclave-ExitusPlus™ are readily biodegradable and not harmful or toxic for humans. : • Does not contain aggressive chemicals
WGK:	1
Storage:	RT
CS:	38220000

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340