

AT50/CWMD3

Orange medium duty clinical waste sacks 5kg



Features

For alternative treatment or incineration. On the roll.

Produced using co-extrusion technology, star-sealed and printed according to UN requirements (UN3291)

30 litres/ 5kg UN approved weight

Code	AT50/CWMD3
Colour	Orange
Sizes	360/560x640mm

AT50/CWMD3

Sectors:

- Medical, Vet, Dental, Pharmaceutical
- Waste Management and Local authority

Packing

50 bags per roll / 10rolls per case



For further technical information, please contact Customer Services on 033 3320 8550 or the Technical team on 033 3320 8555