

KERN 317-12

Single weight, cylindrical, polished stainless steel

KERN[®]



Construction

| | |
|-----------|--------------------------|
| Material: | polished stainless steel |
| Shape: | cylindrical |

Measuring system

| | |
|----------------|------|
| Nominal value: | 2 kg |
| OIML class: | E2 |
| Tolerance: | 3 mg |

Accessories

| | |
|--------------|-------------|
| Plastic box: | 317-120-400 |
| Wooden box: | 317-120-100 |

Service

| | |
|--------------------|---------|
| DAkkS Certificate: | 962-342 |
| Verification: | 952-342 |

Packaging & shipping

| | |
|-------------------------------|----------------|
| Delivery: | 24 h |
| Dimensions packaging (WxDxH): | 1 x 1 x 1 mm |
| Gross weight: | 2 kg |
| Net weight: | 2 kg |
| Shipping method: | Parcel service |