

KERN 317-02

Single weight, cylindrical, polished stainless steel

KERN[®]



Construction

Material:	polished stainless steel
Shape:	cylindrical

Measuring system

Nominal value:	2 g
OIML class:	E2
Tolerance:	0,040 mg

Accessories

Plastic box:	317-020-400
Wooden box:	317-020-100

Service

DAkkS Certificate:	962-332
Verification:	952-332

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	1 x 1 x 1 mm
Gross weight:	2 g
Net weight:	2 g
Shipping method:	Parcel service