

KERN 317-11

Single weight, cylindrical, polished stainless steel

KERN[®]



Construction

Material:	polished stainless steel
Shape:	cylindrical

Measuring system

Nominal value:	1 kg
OIML class:	E2
Tolerance:	1,600 mg

Accessories

Plastic box:	317-110-400
Wooden box:	317-110-100

Service

DAkkS Certificate:	962-341
Verification:	952-341

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	1 x 1 x 1 mm
Gross weight:	1 kg
Net weight:	1 kg
Shipping method:	Parcel service