

KERN 326-04

Single weight, compact shape, polished stainless steel



Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	1 x 1 x 1 mm
Gross weight:	10 g
Net weight:	10 g
Shipping method:	Parcel service

Construction

Material:	polished stainless steel
Shape:	compact shape

Measuring system

Nominal value:	10 g
OIML class:	F1
Tolerance:	0,200 mg

Accessories

Plastic box:	347-050-400
Wooden box:	317-040-100

Service

DAkkS Certificate:	962-434
Verification:	952-434