

IKA

designed for scientists



DB 3.1

/// Data Sheet

Used for test tube combination (6 mm / 12/13 mm / 25 mm)

Pore size: Ø 8.3 / Ø 13.9 / Ø 26.2 mm

No. of holes: 6 / 5 / 3

Hole depth: 48.4 / 48.4 / 48.4 mm

Dimensions (W x D x H): 95 x 76 x 51 mm



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	95 x 51 x 76
Weight [kg]	0.68
Permissible ambient temperature [°C]	5 - 40