

# IKA

designed for scientists



## DB 3.2

/// Data Sheet

Used for centrifuge tube combination (1.5 ml / 15 ml / 50 ml)

Pore size: Ø 11.1 / Ø 17.1 / Ø 29.0 mm

No. of holes: 4 / 3 / 2

Hole depth: 39.1 / 44.5 / 47.6 mm

Dimensions (W x D x H): 95 x 76 x 51 mm



designed for scientists

### Technical Data

Dimensions (W x H x D) [mm]	95 x 51 x 76
Weight [kg]	0.75
Permissible ambient temperature [°C]	5 - 40