

IKA

designed for scientists



DB 3.3

/// Data Sheet

Used for microtube combination (0.5 ml / 1.5 ml / 2.0 ml)

Pore size: Ø 7.9 / Ø 11.1 / Ø 11.5 mm

No. of holes: 6 / 10 / 5

Hole depth: 27.6 / 39.1 / 38.1 mm

Dimensions (W x D x H): 95 x 76 x 51 mm



designed for scientists

Technical Data

Dimensions (W x H x D) [mm]	95 x 51 x 76
Weight [kg]	0.817
Permissible ambient temperature [°C]	5 - 40