

ID-Nr. 30122-001
Model S90.129.120.WDAS



Safety storage cabinet S-CLASSIC-90
Model S90.129.120.WDAS, light grey (RAL 7035)

Safe and approved storage of flammable hazardous substances in working areas

Function / construction:

- **Robust construction and longevity:** triple hinged door, body made of powder-coated sheet steel, safety elements assembled outside the cabinet for increased protection against corrosion
- **User-friendly and safe (model .WDAS):** door open arrest system, convenient and safe removal or storage of containers
- **No unauthorised use:** doors lockable with profile cylinder (integration in an existing locking system possible)
- **Easy alignment:** adjusting aids to compensate for uneven floor
- **Ventilation:** integrated air ducts ready for connection (DN 75) to a technical exhaust system

Dimensions	
Width (external)	1196.00 mm
Depth (external)	616.00 mm
Height (external)	1298.00 mm
Width (internal)	1050.00 mm
Depth (internal)	524.00 mm
Height (internal)	1072.00 mm



The available internal dimensions can vary depending on the interior equipment package.

Weight / Load	
Weight without interior equipment	270 kg
Distributed load	450.00 kg/m ²
Maximum load	600 kg

Ventilation	
Extraction air	75 DN
Minimum exhaust air quantity 10 times	6 m ³ /h
Differential pressure (10 times air change)	5 Pa

Colouring / Surface finish	
Body colour code	RAL 7035
Body colour	light grey
Body surface finish	plain
Door / Drawer colour code	RAL 7035
Door / Drawer colour	light grey
Door / Drawer surface finish	plain

General Information	
Number of doors	2
Door opening angle	89.00°
Interior grid dimension	32 mm
Depth with open doors	1173 mm
Customs tariff number	94032080



Top view

