

Cole-Parmer®
stuart™

Cole-Parmer® Stuart™ SCF4 Microcentrifuge

- Personal size for quick spin downs
- Safely spin samples with dual locking lid
- Quickly change rotors with quick-release adapter
- Built-in recessed rotor bowl
- BioCote® antimicrobial protection





Cole-Parmer® SCF4 Microcentrifuge

The Cole-Parmer™ Stuart SCF4 microcentrifuge provides high speed centrifugation with enhanced safety features. The high speed 8500 rpm and 3884 RCF rating allows users to spin samples rapidly or be used in applications where greater g-force is required.

Prevent accidental opening with the innovative dual safety locking lid. The unit will only spin when the lid is closed and will automatically shut off after 60 minutes of continuous use. The rubber-coated base allows for enhanced stability and reduced vibration.

Technical Specification

| Specification | SCF4 |
|------------------------|---|
| Maximum Speed | 8500 rpm |
| Maximum RCF | 3884 x g |
| Speed Control | Fixed |
| Rotors Included | 8 x 1.5 to 2.0mL tubes, 4 x 0.2mL 8-tube strips |
| Noise level | 60 dB |
| Power | 100-240V 50 60Hz |
| Electrical Plug Type | Detachable US, UK, and EU plugs included |
| Dimensions (w x d x h) | 198 x 185 x 109 |

Ordering Information

| Description | Ordering Number | Series No. | Model No. | Legacy Sku. |
|---|-----------------|------------|-----------|-------------|
| Microcentrifuge | 17414-70 | CF-100 | CF-100 | SCF4 |
| Replacement 8 place rotor | 17414-72 | CF-100 | N/A | SCF4/2 |
| 0.5ml tube adaptors | 17414-73 | CF-100 | N/A | SCF4/1/05 |
| Replacement strip tube rotor, 4 x 0.2ml | 10097-80 | CF-100 | N/A | SCF2/2 |

coleparmer.com

USA:
+1.800.323.4340
+1.847.549.7600

Canada: +1.800.363.5900
China: +86.21.5109.9909
France: +33 (0) 1486 37800
Germany: +49 (0) 9377 92030

India: +1.800.266.1244
Italy: +39 (0) 1313 89513
UK: +44 (0) 1480 272279
All other countries: +1.847.549.7600

Cole-Parmer®
essentials

an antylia scientific company

Find out more at antylia.com