

Economical, High Grade Disinfection with Easy Application



Klercide™ 70|30 Alcohol Range Blended with DI

Sterile, spore free 70|30 alcohol blended with deionised water

A convenient range of either Isopropyl Alcohol (IPA) or Denatured Ethanol sprays which are comfortable and easy to use.

The fully adjustable trigger spray allows the liquid to be dispensed as either a jet or spray, enabling thorough wetting of the surface so that disinfection procedures are effective and minimise product wastage.

Klercide 70|30 Alcohol contains no propellants so does not require any special disposal procedures.

The range is:

- ▲ 0.2 micron filtered, filled and double bagged in a Grade C (ISO class 7) cleanroom
- ▲ Irradiated at no less than 25 kGy to ensure a spore free alcohol
- ▲ Fitted with the SDS trigger spray to give an in use shelf life of three months
- ▲ Available in a range of formats





Klercide™ 70|30 Alcohol Range

Blended with DI

Sterile, spore free 70|30 alcohol blended with deionised water

Why use the Klercide 70 | 30 Alcohol Range?

EFFECTIVE

Blended with deionised water, Klercide 70|30 IPA is 70% v/v Isopropyl Alcohol and Klercide 70|30 Denatured Ethanol is 70% v/v denatured ethanol. The alcohol blend is 0.2 micron filtered, filled and double bagged in a Grade C (ISO class 7) cleanroom, before irradiation using a validated process at no less than 25 kGy to ensure a sterile alcohol.

UNIQUE DELIVERY SYSTEM

The 500ml and 1L trigger spray use the SteriShield Delivery System which has been validated to enable an in-use shelf life of three months. So there is no need to discard unused product, therefore reducing waste.

The correct method of application ensures the optimum and thorough wetting of the surface with the product.

VALIDATION

Manufactured by an ISO 9001:2008 and ISO 14001:2015 accredited company, the Klercide Alcohol Range passes EN test 1276:1997 for bactericidal efficacy.

Certificates of Analysis and Conformity are supplied with every batch and the sterility test method is to an approved pharmacopial method.

The products have a two year shelf life. Full supporting documentation is available on these products.



Klercide 70|30 IPA Blended with DI
0.5L and 1L Sprays, 5L Capped



Klercide 70|30 Denatured Ethanol Blended with DI
0.5L and 1L Sprays, 5L Capped

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences

EUROPE HEADQUARTERS

Richtistr. 7
8304 Wallisellen
Switzerland

www.ecolablifesciences.com

POWERFUL.
EFFECTIVE.
ECONOMICAL.

APPLICATION AREA

For use on floors, walls and hard surfaces.

INSTRUCTIONS FOR USE

Spray: For optimum results hold bottle upright and spray from a distance of 10cm to 20cm. Always close the nozzle after use.

Wetting: For optimum results ensure complete wetting of the surfaces using the spray and wipe technique (max. 40ml/m²).

INGREDIENTS

Klercide 70|30 IPA 100ml contains:
70ml isopropyl in deionised water.

Klercide 70|30 Denatured Ethanol 100ml contains: 70ml denatured ethanol in deionised water.

PACKAGING INFORMATION

UNIT SIZE	UNIT OF SALE
500ml spray	6 bottles x 500ml
1.0 litre spray	6 bottles x 1.0L
5.0 litre capped	4 bottles x 5.0L

FOR FULL PRODUCT ORDER CODES, PLEASE SEE OUR PRODUCT CATALOGUE OR VISIT OUR WEBSITE

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

