



PRODUCT SPECIFICATION

DI Water 5 Litre RFU STERILE



PRODUCT INFORMATION

PRODUCT DESCRIPTION

Sterile and Non-sterile DI water in RFU format

PRODUCT FEATURES

- 5 Litre container complete with a screw on tamper evident air tight cap
- Convenient carrying size and weight
- Visible liquid level

TECHNICAL INFORMATION

MANUFACTURING

The product is filled and capped in ISO 14644-1 Class 5 (in use) unidirectional airflow cabinets.

Packaging is performed within an ISO 14644-1 Class 7 (in use) cleanroom.

PACKAGING

Primary packaging: 5 Litre container complete with a screw on air tight cap.
Secondary packaging: double bagged in polythene bags in ISO 14644.
2 double-bagged bottles in the box.

STERILISATION DETAILS

Gamma Irradiation (minimum 25kGy).

QUALITY ASSURANCE

Product supplied under BS EN ISO 9001 registered systems.

SHELF LIFE

2 years from date of manufacture.

STORAGE CONDITIONS

Store in a dry, well-ventilated area in ambient conditions.

HAZARD STATEMENT

Non-hazardous.

ORDERING/PACKING INFORMATION

Code	Description	Size	Pack Qty	Case Qty	Case Dimensions	Weight
ZSDI82005 - A/B/C/D/E/F/G/I	CleanGuard 6 DI Water 5 Litre RFU Sterile	5L	Each	2	335mm x 330mm x 265mm	15 Kg
ZNDI82005 - A/B/C/D/E/F/G/I	CleanGuard 6 DI Water 5 Litre RFU Non-sterile	5L	Each	2	335mm x 330mm x 265mm	15 Kg