




Sterilin™ 7mL Bijou Containers

- Ideal for small volume samples, such as pediatric samples and other small-volume clinical applications
- Aseptically manufactured
- Suitable for centrifugation at 7,200 x g
- Leak-tested in accordance with BS EN 14254 Annexe D
- Dims, Ø x h, mm: 22.5 x 50.2
- In vitro use only
- CE marked in accordance with the European Directive 98/79/EC

Supplied in packs of 700



Technical Data

Code	Description	Material Base/Cap	Sterility	Pack
CON7522	7mL bijou, no label	PS/PE	AS	700

00476