

GN01

## NITREX Extra Sensitive



Latex-Free

Ultimate Tactile Sensitivity

Strong Chemical Resistance

Passes Viral Penetration Tests (F1671)

Tested for use with Chemotherapy Drugs

NITREX Extra Sensitive examination gloves are manufactured from a high quality nitrile formulation resulting in a strong performing glove with a high level of tactile sensitivity.

The advanced thin film technology utilised in the manufacture process results in a glove, which provides outstanding barrier protection against viruses, bacteria, chemicals and has been tested for permeation performance against 13 chemotherapy drugs.

### Specification

Glove details	Specification
Material	Nitrile
Length	Min. 240mm
Protein Level	Contains no latex proteins
Surface	Micro-textured fingers
Shape	Ambidextrous
Colour	Blue
Sterilisation	Non sterile
Shelf Life	3 Years

### Physical Properties

Glove details	Properties
Force at Break (Newtons)	>6
Freedom from Holes	AQL 1.5
Length (mm)	Min. 240
Cuff Thickness (mm)	0.04
Palm Thickness (mm)	0.05
Finger Thickness (mm)	0.06

AQL (Acceptable Quality Level) refers to the maximum number of defective products that could be considered acceptable during the random sampling of an inspection, in this case freedom from holes in gloves. The lower the number, the fewer the holes and the higher the quality of gloves.

**Sizing**

Size	Palm width
Extra Small	75 ± 5mm
Small	85 ± 5mm
Medium	95 ± 5mm
Large	106 ± 5mm
Extra Large	116 ± 5mm

**Quantities**

Code	Box Qty	Case Qty
GN01E	200 Gloves	10 Boxes
GN01S	200 Gloves	10 Boxes
GN01M	200 Gloves	10 Boxes
GN01L	200 Gloves	10 Boxes
GN01X	200 Gloves	10 Boxes

**Re-order Codes**

Size	MPC	NHS Code
Extra Small	GN01E	FTE1744
Small	GN01S	FTE1745
Medium	GN01M	FTE1746
Large	GN01L	FTE1747
Extra Large	GN01X	FTE1748

**Quality standards:** Manufactured in a facility holding ISO 9001 and ISO 13485.

**Product standards:** The product is tested and manufactured to the following standards: EN455-1, EN455-2, EN455-3, EN455-4.

**Storage:** Store in a cool, dry place away from sources of heat or direct sunlight and shielded from ozone and UV light.

**Disposal:** Dispose of gloves as clinical waste. Collection case and transit case can be recycled as paper or disposed of as clinical waste. Please follow local regulation and guidance.

**Shelf life:** Three (3) years from date of manufacture.

**Country of origin:** Malaysia

**Notice:** This product is free from natural rubber latex and is suitable for use in latex-free environments. Chemical accelerators are used in the manufacture process, which can cause skin reactions in some users. Stop use of the product if skin reactions occur.



**ASTM  
F1671**

Find out more at [www.medicareproducts.com](http://www.medicareproducts.com)

Medicare Products Ltd, Unit B Dolphin Way, Purfleet, Essex RM19 1NZ, UK  
 e: [enquiries@medicareproducts.com](mailto:enquiries@medicareproducts.com) t: +44 (0)1708 863868