

Serological pipette 1ml ClearLine®
sterile in individual packaging


**Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: yellow

- Sterile (Gamma)
- DNase and RNase free
- Non-pyrogenic

- Volume 1.0 ml
- Graduations: 0.01ml
- Negative graduations 0.4ml

Total length (A)	271mm
Ø Tip (B)	1.9mm
Ø Tip (C)	4.5mm

**Conditioning**

- Individual paper/plastic packaging
- Box of 10 bags of 100 units

Units :

1000 Tips