


Serological pipette 5ml ClearLine® in bulk**Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: blue

- Non-sterile
- DNase and RNase free
- Non-pyrogenic

- Volume 5 ml
- Graduations: 0.1 ml
- Negative graduations 2.5 ml

Total length (A)	340 mm
Ø Tip (B)	3.2 mm
Ø Tip (C)	8.1 mm

**Conditioning**

- Bulk packaging
- Box of 20 bags of 25 units

Units :

500 Tips