


Serological pipette 25ml ClearLine® sterile in individual packaging**Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: red

- Sterile (Gamma)
- DNase and RNase free
- Non-pyrogenic

- Volume 25 ml
- Graduations: 0.2 ml
- Negative graduations 9 ml

Total length (A)	342 mm
∅ Tip (B)	4.7 mm
∅ Tip (C)	14.3 mm
∅ Tip (D)	8.1 mm

**Conditioning**

- Individual paper/plastic packaging
- Box of 4 bags of 50 units

Units :

200 Tips