


**Serological pipette 50ml ClearLine<sup>®</sup> sterile in individual packaging****Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: purple
  
- Sterile (Gamma)
- DNase and RNase free
- Non-pyrogenic
  
- Volume 50 ml
- Graduations: 0.5 ml
- Negative graduations 7 ml

Total length (A)	365 mm
∅ Tip (B)	5.6 mm
∅ Tip (C)	18 mm
∅ Tip (D)	8.2 mm

**Conditioning**

- Individual paper/plastic packaging
- Box of 4 bags of 25 units

**Units :**

100 Tips