


Serological pipette 50ml ClearLine® in bulk**Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: purple

- Non-sterile
- DNase and RNase free
- Non-pyrogenic

- Volume 50 ml
- Graduations: 0.5 ml
- Negative graduations 7 ml

Total length (A)	365 mm
Ø Tip (B)	5.6 mm
Ø Tip (C)	18 mm
Ø Tip (D)	8.2 mm

**Conditioning**

- Bulk packaging
- Box of 4 bags of 25 units

Units :

200 Tips