


**Serological pipette 10ml ClearLine® sterile in individual packaging****Features**

- Disposable serological pipettes
- Material: polystyrene
- Cotton pipette: yes
- Color code: orange
  
- Sterile (Gamma)
- DNase and RNase free
- Non-pyrogenic
  
- Volume 10 ml
- Graduations: 0.1 ml
- Negative graduations 2.5 ml

Total length (A)	344 mm
Ø Tip (B)	3.6 mm
Ø Tip (C)	10 mm
Ø Tip (D)	7.8 mm

**Conditioning**

- Individual paper/plastic packaging
- Box of 4 bags of 50 units

**Units :**

200 Tips